

QUICKVUE®

H. pylori gII® TEST

Test and Treat Today!

Now there's a tool that helps you determine the presence of *H. pylori* in a few minutes.

Before you write another prescription for an H-2 antagonist or proton pump inhibitor, take a few minutes to test your ulcer patients for *H. pylori*. *H. pylori* is

responsible for more than 80% of all gastric ulcers and is classified by the W.H.O. as a group I carcinogen.

Quidel has a rapid *H. pylori* test, the QuickVue® H. pylori gII® test, an in-office test that allows you to Test and Treat patients in the same office visit.



Rapid Diagnosis and Treatment

- Results in 5 minutes, or less.
- Patients can be tested and treated in the same office visit.

Easy to Use

- Simple, one-step procedure from only 2 hanging drops of whole blood; serum and plasma can also be used.

Excellent Performance

- 92% overall agreement with EIA.

Convenient

- Available in both 10 and 30 test kit sizes.
- Room temperature storage.
- 14-month shelf life.
- No costly instrumentation needed.
- Prompt diagnosis and treatment equals high patient satisfaction for your practice.

CLIA Classification


- CLIA Waived when using whole blood specimens.
- CLIA Moderately Complex when using serum or plasma.

QUIDEL®
CORPORATION


For more information on QuickVue H. pylori gII test kits, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517.

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • www.quidel.com

Procedure:



Results:



To order: Call your Quidel Account Manager at 1.800.874.1517. Outside the U.S. call: 858.552.1100.



QuickVue H. pylori gII – 10 Test Kit:
Catalog #0W009

QuickVue H. pylori gII – 30 Test Kit:
Catalog #0W010



For more information on QuickVue H. pylori gII test kits, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517.

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • www.quidel.com

Distributor Information: