

# QUICKVUE<sup>®</sup>

Influenza TEST


## Test and Treat Today!

Diagnosing influenza can be confusing. One third of all upper respiratory infections are caused by influenza. This leads to over 36,000 deaths, 114,000 hospitalizations and accounts for over \$12 billion spent in lost work and related medical expenses each year in the U.S.

Quidel's CLIA-waived QuickVue Influenza

test provides you with test results that allow you to Test and Treat in less than 10 minutes.

Prompt diagnosis and treatment of influenza is a solution to the overuse of antibiotics, allows for timely use of antivirals, and prevents healthcare costs related to unnecessary testing and treatment.



### Rapid Diagnosis and Treatment

- Results in 10 minutes or less (right in the exam room) for better patient management decisions.
- Most patients tested can be diagnosed and treated in just one office visit.

### Easy to Use

- 1 step (nasal wash/aspirate), 3 steps (nasal swab), 1 reagent, 30 seconds hands-on time.
- First CLIA-waived flu test.

### Accurate Detection

- QuickVue Influenza test detects both types A and B viral antigens.
- Nasal swab procedure:  
Sensitivity: 73% (95% C.I. 67%-78%)

Specificity: 96% (95% C.I. 93%-98%)

Accuracy: 87%

- Nasal wash/nasal aspirate procedure:  
Sensitivity: 81% (95% C.I. 72%-88%)  
Specificity: 99% (95% C.I. 93%-99%)  
Accuracy: 94%

### Convenient Packaging


- Each kit contains everything needed to perform the test, including both positive and negative controls.
- Room temperature storage with up to 24-month shelf life from date of manufacture.
- Transport swab set also available separately.

**QUIDEL**<sup>®</sup>  
CORPORATION


For more information on QuickVue Influenza test kits, visit our website or call your Quidel Account Manager at 1.800.874.1517 (U.S. only), or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [www.quidel.com](http://www.quidel.com)

## Nasal swab procedure:



## Nasal wash/nasal aspirate procedure:



Refer to package insert for complete instructions.



To order: Call your Quidel Account Manager or your local distributor

QuickVue Influenza – 25 Test Kit: Catalog #00317  
 QuickVue Swab Transport Tube Set – Catalog 20171



For more information on QuickVue Influenza test kits, visit our website or call your Quidel Account Manager at 1.800.874.1517 (U.S. only), or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [www.quidel.com](http://www.quidel.com)

Distributor Information: