



DR Panel Accessory Items

Drop deterrence for your DR Panel

800.752.4918 | www.reinaimaging.com

Reina
IMAGING

X-Ray Cassette and Repair Co. Inc.

Protect-A-Grid® DRP Premium Grid Encasement

The Protect-A-Grid DR panel premium grid encasement allows for safe and effective gridded examinations within the direct radiographic modality. The PAG-DPR is a drop-on style that snaps on and off your existing DR panel quickly and effectively. Allow for superlative image detail and contrast with Reina Imaging's PAG-DRP premium grid encasement.


Product advantages:

- Sturdy heavy duty aluminum frame
- Familiar Drop-On style
- User friendly

Models available to accommodate all DR panel makes and sizes

Technical Data:

- Dimensions:
- Length – 19 ½ in.
 - Width – 19 in.
 - Height – 1 in.
 - Weight – 4.85 lbs.



Trans-Former™ DRP Holder

Reminiscent to your old wall bucky; this mobile or wall-mounted DR panel holder is made to order for every OEM DR panel on the market. Have the instant ability to fully and effectively utilize your DR panel in every exam room for any scenario. Digital imaging allows for wider exposure latitude than conventional radiography, now the DRP Holder from Reina Imaging gives the department the same width of versatility throughout the entire x-ray department.

Product Advantages:

- Mobile or wall-mounted versions
- Panel can be positioned portrait or landscape
- DR Panel always perpendicular to ground
- Holder can extend from floor to ceiling
- Pushbutton counterbalanced pulley system
- Non-PBL
- No electrical outlet needed
- Removable grid

Technical Data:

- Dimensions:
- Length – 23 in.
 - Width – 20 in.
 - Height – 54 in.
 - Weight – 4.85 lbs.
 - Highest position – 73 in.
 - Lowest position – 36 in.

DRP Weight Distribution Cap

All medical personnel are well aware of two pressing issues w/in the radiographic industry: Morbid obesity and high costs of DR technology. The DRP Weight Distribution Cap protects your imaging investment by distributing the weight of the patient away from the sensitive DR panel. Engineered for every OEM's DR panel, the DRP-WD has a quick and easy snap-on application creating a simple user-friendly system.



Product advantages:

- Snap-On application
- Lightweight
- Durable
- Can withstand 750lbs.

Technical Data:

- Dimensions:
- Length – 19.50 in.
 - Width – 19 in.
 - Height – 1 in.
 - Weight – 4.25 lbs.

DRP Table Top Holder

Taking shoot-through exams was difficult enough with conventional radiography. Now with the advent of direct radiography the image receptor has become larger, more cumbersome and increasingly heavy. The DRP Table Top Holder from Reina Imaging allows the radiographer to position the DR panel with ease and precision. Custom fabricated for all DR panels the DRP-TTH can also be angled caudal and cephalic depending on the examination.

Product Advantages:

- Sturdy construction
- Accommodates all DR panel sizes
- Holder rotates 360 degrees
- Extends 27 in.
- User friendly





Reinventing Radiography...One Image at a Time™