

FILM LOADING BINS

ALL STEEL CONSTRUCTION / WELDED SEAMS AND CORNERS

DOUBLE STEEL LIGHT TRAP

FULL WIDTH STAINLESS STEEL DOOR HINGE

SIZES AND TYPES TO MEET MOST REQUIREMENTS



These loading bins eliminate the problem of fogging film due to improper storage. A deep double steel light trap around the door gives maximum protection and assures film safety. Bins have a stainless steel full width hinge to insure trouble free action. All corners and seams are welded for strength, guaranteeing no light leakage that can fog film.

The interior of the hinged bin is divided into five or six compartments, each capable of holding up to 14"x17" films. Each bin comes with inserts to accommodate smaller film sizes and can be arranged to suit any film assortment.

With the optional built in light limit switch **127R** installed in a film loading bin and wired to the room lighting, the lights will automatically be turned off whenever the loading bin is opened.

The **125X Film Loading Bin** was designed for areas where space is limited, measuring only 18.375" wide.

The **129 Film Loading Bin** is lined with .031" lead.

Choose the Film Loading Bin most appropriate to your needs and design a darkroom by selecting any combination of cabinets from either the 15" Compact Darkroom Bench Cabinets (see page G-9) or the Heavy Duty 21.5" Deluxe Darkroom Bench Cabinets (see page G-7). See chart below for the proper base to use with the 21.5" bench cabinet.

Bases add 3.5" to the height of a bin, tops add 1", and top and base add 4.5" to the total height of a bin.

Finish: All Bins and Bases: Standard finish - light gray baked enamel. Formica tops are white.

* **Optional Items:** Room Light Control **127R**

FILM LOADING BINS

NO. OF COMPART.	EXTERIOR DIM			APPROX. SHIP WT	LOADING BIN	CATALOG NUMBERS				
	WIDTH	HGT	DEPTH			BIN W/15" DEEP BASE	BIN W/15" DEEP TOP	BIN W/TOP & BASE	BIN W/BASE TO MATCH TO 21.5" CAB	
5	22.5"	33"	15"	90 lbs	125	125B	126	126B	125 & 68R	
6	22.5"	33"	15"	90 lbs	125K	125KB	125KT	125KTB	125K & 68R	
5	18.375"	33"	12.5"	75 lbs	125X	125XB	125XT	125XTB	125X & 68RX	
.031" LEAD LINED										
5	22.5"	33"	15"	130 lbs	129	129B	129T	129BT	129 & 68R	

TABLE TOP FILM LOADING BINS



ALL STEEL CONSTRUCTION / WELDED SEAMS AND CORNERS

SMALL AND COMPACT

IDEAL FOR MAMMOGRAPHY Model# 128 holds 24x30cm and smaller

Model# 128A holds 14"x17" and smaller

FULL WIDTH BLACK ENAMEL DOOR HINGE

BUILT TO LAST

These small, space saving counter top Film Loading Bins are ideal for any darkroom. The level top can be used as an alternate work surface to minimize loss of counter space. The **128** is designed specifically to hold 24x30cm and 18x24cm (or smaller) films making it a must for the mammography darkroom. The **128A** holds 14"x17" (or smaller) film and will complement any darkroom where space is limited. A deep double steel light trap around the door gives maximum protection and assures film safety. All corners and seams are welded for strength guaranteeing no light leakage that can fog film. The door has a full width hinge to insure trouble free action for the life of the cabinet.

The interior of the **128** is divided into two compartments by a welded shelf and is finished in black baked enamel. Each compartment holds up to two full boxes of film. The interior of the **128A** has five compartments (each capable of holding one box of 14"x17" film).

An optional built in light limit switch is available for either the 128 or the 128A. With the **127A** Light Limit Switch installed in the Film Loading Bin and wired to the room lighting, the lights will automatically be turned off whenever the bin is opened.

TABLE TOP FILM LOADING BINS

NO. OF COMPART.	EXTERIOR DIM			APPROX. SHIP WT	CATALOG NO.
	WIDTH	HGT	DEPTH		
2	13"	10.75"	15"	25 lbs	128
5	17.5"	13.25"	20.75"	40 lbs	128A

Finish: Black baked enamel.

* **Optional Items:** Room Light Control **127A**