

FILM IDENTIFICATION PRINTERS AND CARDS

AUTOMATIC or MANUAL MODELS
 TABLE TOP or RECESSED MODELS
 CLEAR, SHARP IDENTIFICATION

AUTOMATIC I.D. PRINTERS

The automatic timing circuit and capacitor discharge system built into the Automatic I.D. Printer consistently provides sharp, clear patient data on every X-Ray film. The printer is constructed of a heavy duty all steel body, a Stainless Steel door with a foam rubber contact pad and a Stainless Steel face plate with guide rails for perfect positioning. Stainless Steel parts are polished to a satin finish, steel parts are finished in light gray baked enamel. Recessed model has a baked enamel flange for recessing into a counter top.

AUTOMATIC I.D. PRINTERS OPERATE ON 118 volts 60HZ AC. Other voltages available on special order.

AUTOMATIC I.D. PRINTER

DESCRIPTION	EXTERIOR DIM			APPROX. SHIP WT	CATALOG NO.
	WIDTH	LENGTH	DEPTH		
Table Top	8"	8"	2.5"	5 lbs	144
Recessed	10"	10"	2.5"	6 lbs	144X



MANUAL I.D. PRINTERS

Exposure in these models is not automatically timed. The duration of the exposure is determined by the amount of time the printer door is held in the down position. The printer is constructed of a heavy duty all steel body, a Stainless Steel door with a foam rubber contact pad and a Stainless Steel face plate with guide rails for perfect positioning. Stainless Steel parts are polished to a satin finish, steel parts are finished in light gray baked enamel. Recessed model has a baked enamel flange for recessing into a counter top.

MANUAL I.D. PRINTER

DESCRIPTION	EXTERIOR DIM			APPROX. SHIP WT	CATALOG NO.
	WIDTH	LENGTH	DEPTH		
Table Top	8"	8"	2.5"	5 lbs	143
Recessed	10"	10"	2.5"	6 lbs	143X

Custom personalized imprint

J. Smith M.D.	
100 Park Ave. New York, N.Y. 10010	
Patient	
No.	Date
Exam.	

143C
 Identification Card

Name	_____
No.	Date
Name	_____
Address	_____
Case No.	_____
Submitted for:	Radiographs Fluoroscopy
Suspected Pathology:	_____
Reason to be examined:	_____
Drawings	may not be removed
Referred by	M.D.
Radiologist	M.D.

IDENTIFICATION CARDS

Available for use with the above printers are special 3" x 5" translucent cards which will assure clear, sharp transfer of information onto the X-Ray film. Cards are available with a standard imprint format or a custom personalized two line imprint.

Identification cards are available in multiples of **500 only**. No partial carton orders.

Custom personalized imprints are available with a **minimum order of 1000** (for over 1000 order must be in multiples of 500). For personalized imprint a sample card or written order is required.

143C - Identification Cards, 500 per CTN