

Multi-specialty Exam Table

Standard Ergonomic Features

- 13" x 13" Manual imaging drop section for cardiac exams
- 13" x 13" Wedge Section for patient support
- Dropping footboard
- Retractable steel stirrups
- Trendelenburg & reverse (+/- 15°)
- Electric Fowler position
- Electric height adjustment 25"- 37"
- Width = 30"; Length = 74"
 - Also available in a 27" width
- 5" Total lock casters
- Solid steel frame
- Powder coat finish
- Infinite position electric hand control
- Storage shelf
- Lift/Load Capacity 600 lbs
- UL Approved 2601-1 & CSA

Available Options:

- Retractable safety side rails
- Detachable arm extensions for upper extremity vascular studies
- IV pole
- Power controlled imaging drop section
- Power controlled footboard
- Ultra-low column height adjustable from 20"-36"
- 5" Central locking casters



Risk factors addressed by this exam table:

Twisting and bending of the trunk
Unnatural and static postures
Reach and arm abduction
Wrist flexion and extension

Item # MSET-DLR