

Image-Vet® *SuperVista*

USB Intraoral Camera System

- Superior image quality
- Take dental diagnostics to the next level with this easy-to-use imaging tool
- The convenience of an intraoral camera that integrates with a digital dental radiography system
- Seamlessly compatible with the EVA-Vet® Digital Dental Imaging System and ViewAll® imaging software
- USB 2.0 connectivity - No external power supplies or proprietary capture cards



*The perfect companion for the
EVA-Vet Digital Dental Sensor!*

Distributed by CAN-med Healthcare
East 1 866 277 3577 / West 1 800 661 5671
www.canmedhealthcare.com

www.ImageWorksCorporation.com


ImageWorks
Generations of Imaging

ImageWorks Dental Imaging Solutions For Every Veterinary Practice

The ImageWorks SuperVista intraoral camera is a simple to use handheld USB camera. This wand-style camera uses photography for diagnostic evaluation and client presentation. Simply position the camera close to a region of interest and use the single-click shutter button to capture as many images as you would like. After you snap a picture it saves in the software and the camera turns live again, ready for the next picture.

Using ViewAll with the SuperVista means that you can save dental photographs alongside dental radiographs in a patient's file. This expansion of the typical ViewAll usage provides more information to referring doctors and to clients, and allows an office to quickly explain through visual records.

The SuperVista has been widely regarded as a critical aid in documentation and explanation.

Once people start taking pictures they never stop!

SuperVista Specifications:

Working distance:	5 ~ 35mm
Light output:	≥ 500 Lux
Resolution:	640 X 480 ~ 1280 X 1024
Field of view:	90° ± 5%
Power:	5 Vdc (via USB 2.0)
Video output:	USB 2.0 (full motion digital)
Imager:	1/4 inch CMOS
Illumination:	(6) LED's
Temp & humidity:	5° ~ 68° C / ≤ 80%
Handpiece weight:	40 grams
Operating system:	Windows 7/Vista/XP
Drivers:	Windows 7: 64 bit & 32 bit Windows Vista & XP: 32 bit

P/N: 9992470700

Each SuperVista Camera System Includes:

- (1) USB Camera
- (1) USB Cable
- (1) Wall/Table mounting holster
- (2) Mounting screws
- (1) Driver CD
- (1 Pkg) Protective sheaths
- (1) User Manual



Camera



Holster



No Film or Chemicals!

Veterinary Dental Imaging Solutions from ImageWorks

EVA-Vet® Digital Dental Sensors • Image-Vet® DC Dental X-Ray • ViewAll® Software

Image-Vet® SuperVista Intraoral Camera • Image-Vet® SAINT 3D Tomography System • Image-Vet® 4G Dental X-Ray



ImageWorks
250 Clearbrook Road, Suite 240
Elmsford, NY 10523 USA

1.800.592.6666
P 914.592.6100
F 914.592.6148

www.ImageWorksCorporation.com

Distributed by CAN-med Healthcare
East 1 866 277 3577 / West 1 800 661 5671
www.canmedhealthcare.com

Products are continuously under review in the spirit of technical advancement. Actual specifications are subject to improvement or modification without notice. Copyright © 2011 ImageWorks. All rights reserved. Reproduction in whole or any part of these contents without written permission is prohibited. Printed in U.S.A.

SMO-165-00