

CM-GP2

Components

- Generators – various configurations
- Elevating Table
- FMT
- Wallstand



CM-GP2

S-HL Elevating Table

- 750# Load Rating
- Vertical Travel – 19.5" – 35.5"
- Bucky Travel – 20"
- Tabletop Longitudinal – 43.3"
- Tabletop Lateral – 9.5"
- Collision Protection



CM-GP2

S-FMT

- Vertical Travel – 59"
- Lateral Travel – 11.8"
- Column Rotation
- 86.6" Floor Track



CM-GP2

Wall Stand

- Floor to Wall Mount
- Vertical Travel 16" to 75"
- Electromagnetic Locks
- Reciprocating Bucky
- 103L, 10:1 ratio Grid
- Manual Tray

