

H/F Generators

Features –

- Wide selection –
 - 32Kw to 80Kw
 - Battery Powered
 - R/F Generators
 - Touch Screen upgrades
- 2 Point or 3 Point



H/F Generators

Features -

- 2560 APR Techniques Possible
 - 8 Regions, 8 Anatomical
 - 8 Projections, 20 Thicknesses
- Technique adjustment by cm.
- Actual mAs
- Last Exposure/Actual Time



Hercules H/F Generator

- 150 kHz
- 50, 65 or 80kW
- Touch Screen Console
- Remote Diagnostics
- Hi Speed Starter
- 14.1" x 19.5" x 16.5"



Hercules H/F Cont'

- Auto Tube Warm up & Self-diagnostics at Start up
- Multiple Languages
- Multiple Users
- 480v. 3 phase line required

