

Advantage Plus

by **TXR**

Distributed By:


OPTI-RAD Radiographic Suite




Advantage Plus

by **TXR**

OPTI-RAD Radiographic Suite



DIMENSIONS



TXR • 5481 Skyland Blvd., E • Cottondale, AL 35453

Phone: 205-556-3803 • Fax: 205-556-3824 • 800-897-9729 • www.txr.com • sales@txr.com

COMPLIES WITH FDA AND PRODUCT SAFETY REQUIREMENTS. • SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE

0510 ORRS

- High Quality
- Hospital, Clinic or Private Practice
- DR, CCD, CR or Film Solutions

Advantage Plus

by **TXR**

OPTI-RAD Radiographic Suite

OPTI-RAD

The OPTI-RAD is an affordable general purpose radiographic system that can be configured with a choice of table and tube stand to accomplish a wide range of table, wall stand, and wheel chair examinations.

The OPTI-RAD was designed to deliver high performance imaging for a demanding environment. The solid design ensures outstanding reliability for the busy hospital and private practice settings.



AFFORDABILITY, DURABILITY AND FLEXIBILITY



LATERAL EXAMS

The tube stand of the OPTI-RAD can be easily moved and rotated so that the positioning for lateral exposures is very fast and accurate.



Advantage Plus

by **TXR**

OPTI-RAD Radiographic Suite

FAR REACHING DOWN TRAVEL



The counter-balanced tube stand provides smooth vertical movement. For exposures using the wall stand, the tube column can be moved past the table on either the right or left-hand side. On either side, the vertical range of movement of the tube assembly can be lowered to a source-to-floor distance of 15.75 inches.

STANDARD FEATURES

FLOOR MOUNT TUBE STAND:

- 11.8" Transverse travel
- 180° Column rotation
- Vertical travel 15.7" to 74.75"
- Electromagnetic locks
- 86.6" Floor track provides 78.7" longitudinal tube travel

FIXED HEIGHT 4 WAY FLOAT TOP TABLE:

- 78.75" x 28.5" Table top
- 440 Lb. patient load
- 23.5" Longitudinal and 6.7" lateral table top movement
- 20" Bucky travel
- 63.4" Radiographic coverage
- Super speed bucky with 10:1 ratio, 103 line grid
- Electromagnetic locks

WALL STAND:

- Floor-to-Wall mounted
- 15.75" to 75" Vertical travel
- Electromagnetic locks
- Super speed bucky with 10:1 ratio, 103 line grid

EXTENSIVE COVERAGE

- Wheelchair Patients
- Chest exams
- Leg exams
- Weight bearing exams
- Table exams



VERSATILITY



OPTIONAL ELEVATING TABLE

- 86.5" x 31.5" (220 cm x 80 cm) Table top
- 43.3" (110 cm) Longitudinal and 9.5" (24 cm) lateral table top movement
- 20.25" (52 cm) Vertical travel



OPTIONS:

- 118" Floor track for floor mount tubestand
- Digital readout of SID and angulation on floor mount tube stand
- 750 lb. patient load elevating table
- Various grids
- Lateral cassette holder
- Hand grips
- Ion chambers
- DR panels