

Universal Radiographic Systems

Three System Configurations

URS – TS



URS – T LP



URS – T LP Plus



Four Capture Options

- ***Film, CR, CCD, Flat Panel***

URS –TS



- *Straight Arm*
- *Most Compact*
- *Ideal For Chest, Extremity, & Ambulatory Care*
- *Motorized SID And Vertical Positioning*
- *Manual Angulation Of Receptor*
- *Degree Of Angle Displayed At Handle Bar*
- *Recumbent Views Possible With Mobile Table*
- *No Lateral Cross Table Ability*
- *CCD not an option*
- *Cost Effective*

URS -T LP



- **C-arm Design For Lateral Cross-Table Exams**
- **Fully Motorized Movements**
 - **SID**
 - **Vertical**
 - **Rotation**
 - **Image Receptor**
- **Dual Speed Motorized Movements**
- **Touch Pad Controls At Both Receptor & Tube head**
- **Digital Display of SID and Rotation**
- **Multiple Collision Avoidance Sensors**
- **CCD an Option**

URS –T LP Plus

Most Advanced Single Detector System Available



Advanced Features Layered Upon Same Features of The LP

- **IR Remote Control**
- **Auto Collimation from APR**
- **Advanced Feature Large LCD Command Center**
 - **Full Generator & APR Integration**
 - **Positioning Guide**
 - **Selection of up to Ten Programmed Positions**
 - **Worklist Query with Digital**

URS -T LP Plus



- Technique is Displayed at URS Control Display
- Automatic Positioning of Structure for Selected View
- Automatic Position Requires “Dead Man” Switch
- Technique and Position may be Fine Tuned

Mobile Table Options



**Top – Laminate or Carbon Fiber
441 lb Weight Capacity (both
models)**

78.7" x 25.6" x 17.6"

Locking Castors