



## Digital Support Cart

24" Wide Base with Laminate Top for CR readers, digitizers, viewers, monitors, keyboards, and other support equipment.

**NEW!**



- Use anywhere you need a computer, compact equipment or mobile equipment.
- Overall size (not including the top): 24"wide x 20"deep x 38 1/4"tall,
- Height of open shelf is 9".
- Height of shelf w/dividers is 20".
- Height of three open shelves are 9"/10"/9".
- Width of tall divided spaces are 6"/6"/12".
- Laminate top includes a hole for adding the DSM-MCP monitor mount.

### TOP "C"

25"wide x 25"deep

- Compatible suggestions:
- Vidar Diagnostic ProAdvantage
  - Kodak DryView 5800/5850
  - Fuji DryPix 2000
  - Konica/Minolta DryPro 832
  - Codonics® Horizon

Cat.# DSC-2412C  
\$695.00



Cat.# DSC-2421C  
\$695.00



Cat.# DSC-2430C  
\$695.00



### TOP "D"

30"wide x 30"deep

- Compatible suggestions:
- Kodak DryView 5800/5850
  - Array 2905
  - AllPro ScanX 14
  - Agfa CR30 / Quantum CR30
  - Kodak CR120/140/260/500
  - Sony Filmstation

Cat.# DSC-2412D  
\$745.00



Cat.# DSC-2421D  
\$745.00



Cat.# DSC-2430D  
\$745.00



### TOP "E"

36"wide x 21"deep

- Compatible suggestions:
- iCRCO 6125I/2000SL
  - Agfa DryStar 5300
  - AllPro ScanX 14
  - iCRCO ICR3600

Cat.# DSC-2412E  
\$785.00



Cat.# DSC-2421E  
\$785.00



Cat.# DSC-2430E  
\$785.00



# 10 Imaging Support Carts (digital support carts)



**Digital Support Cart**  
 18" Wide Base with Laminate Top for **NEW!**  
 CR readers, digitizers, viewers, monitors, keyboards, and other support equipment.

- Use anywhere you need a computer, compact equipment or mobile equipment.
- Overall size (not including the top): 18"wide x 20"deep x 38<sup>1</sup>/<sub>4</sub>"tall,
- Height of open shelf is 9".
- Height of shelf w/dividers is 20".
- Height of three open shelves are 9"/10"/9".
- Width of tall divided spaces are 3"/3"/12".
- Laminate top includes a hole for adding the DSM-MCP monitor mount.

- DIGITAL SUPPORT CARTS**
- Two sizes of carts:
    - 18"wide x 20"deep x 38<sup>1</sup>/<sub>4</sub>"tall,
    - 24"wide x 20"deep x 38<sup>1</sup>/<sub>4</sub>"tall.
  - Three interior configurations:
    - open shelf over large shelf w/dividers,
    - large shelf w/dividers over open shelf,
    - three open shelves.
  - Four sizes of laminate tops (choose the best for your application).
  - All carts feature three cable holes in the back at each shelf level.



Unsure of your application? Call us for recommendations or custom design configurations.

**TOP "A"**  
 20"wide x 22"deep

- Compatible suggestions:
- Vidar Sierra Plus
  - Radlink LaserPro 16
  - LogicRad
  - AllPro Imaging SCAN12
  - Protec ProScan35



Cat.# DSC-1812A  
 \$575.00



Cat.# DSC-1821A  
 \$575.00



Cat.# DSC-1830A  
 \$575.00

**TOP "B"**  
 25"wide x 25"deep

- Compatible suggestions:
- Vidar Diagnostic ProAdvantage
  - Kodak DryView 5800/5850
  - Fuji DryPix 2000
  - Konica/Minolta DryPro 832
  - Codonics® Horizon



Cat.# DSC-1812B  
 \$595.00



Cat.# DSC-1821B  
 \$595.00



Cat.# DSC-1830B  
 \$595.00

# 12 Imaging Support Carts (digital support cart options)

## Retrofit Door For 18" Wide Digital Support Cart



- Secure computers, plates and equipment behind the locking door.
- Two hinges easily mount to existing cart frame.
- Includes cam style lock.

Cat.# DSD-18  
\$235.00  
(cart must be ordered separately)

Cat.# DSD-24  
\$265.00  
(cart must be ordered separately)

## Retrofit Door For 24" Wide Digital Support Cart



**NEW!**



## Monitor Support For 18" or 24" Wide Digital Support Carts

- Monitor post mounts into a hole in the cart's laminate top.
- Height adjusts with a screw crank so that stability is maintained during rolling vibrations.
- Includes standard VESA mount.
- Leaves room for a keyboard on top of the cart.




(equipment not included)

Cat.# DSM-MCP  
\$325.00  
(cart must be ordered separately)



## Protective Carrying Cases

For carrying CR-plates, wireless DR-plates, laptops, sensitive equipment and parts.

- Maximum protection and storage of sensitive plates and small equipment.
- Interdepartment transport and other mobile applications.
- Water proof • dust proof • air proof • crush proof • corrosion proof.
- "empty case" does not include foam inside; "with foam insert" includes a customizable foam interior.
- Injection molded ultra high impact copolymer resin with a lifetime warranty. 
- Many other sizes and styles available.

**NEW!**  
*Virtually Indestructible!*



### laptop case

Cat.# PXC-P30 \$235.00 (empty case)  
Cat.# PXC-P31 \$255.00 (with foam insert)

- Easy-open, key-lockable latches.
- Fits notebook computers up to 17".
- Automatic pressure equalization valve.
- Meets airline carry-on regulations.
- Interior: 18.2" x 12.1"x 5.2".
- Exterior: 20" x 14.7"x 5.8".
- Ships @ 8 lbs



### medium case

Cat.# PXC-M40 \$215.00 (empty case)  
Cat.# PXC-M41 \$240.00 (with foam insert)

- Stainless steel hardware and padlock protectors.
- "O" ring & automatic pressure equalization valve.
- Layered soft-grip handle.
- Interior: 17.8" x 12.7"x 6.75".
- Exterior: 19.1" x 15.5"x 7.6".
- Ships @ 11 lbs.



### top loading case

Cat.# PXC-T70 \$299.00 (empty case)  
Cat.# PXC-T71 \$315.00 (with foam insert)

- Two edge casters and a retractable handle.
- Easy-open, key-lockable latches.
- Layered soft-grip handle.
- Interior: 17.1" x 7.5"x 16".
- Exterior: 19.7" x 12"x 18".
- Ships @ 8 lbs.



### dr/cr-plate case

Cat.# PXC-D91 \$425.00 (with foam insert)

- Perfect for securely carrying expensive plates from place to place.
- Edge casters and fold-down handles.
- Interior: 21.6" x 16.6"x 10.5".
- Exterior: 24.7" x 19.6"x 11.9".
- Ships @ 23 lbs



### large mobile case

Cat.# PXC-L80 \$425.00 (empty case)  
Cat.# PXC-L81 \$465.00 (with foam insert)

- Four wheels and fold-down handles are internal and will reduce the usable interior area.
- Stainless steel hardware and padlock protectors.
- Layered soft-grip handle.
- Interior: 28.5" x 17.4"x 10.5".
- Exterior: 32.5" x 20.5"x 11.3".
- Ships @ 31 lbs.

# 14 Imaging Support Carts (cr-plate protective transporters)

- Each compartment has special impact-resistant padding.
- All impact surfaces have protective rubber edging.
- Space saving carts have sturdy 4" plate casters.
- All steel construction with durable powder coat finish.

**Protect your investment!**

All "CRC" CR-Plate Units have special, **impact resistant padding** on shelves and vertical surfaces to **prevent damage** to expensive CR-Plates.

## Tall CR-Plate Cart



Cat.# CRC-38 \$794.00

- Eight 6" wide compartments @ 16½" vertical clearance.
- Nominal size: 24" x 18" x 42½"h.
- Shipping weight: 130 lbs.



All CRC Units ship:



Medical Grade Casters:  
two locking swivel  
two stationary

## Mini CR-Plate Cart



Cat.# CRC-24 \$794.00

- Four 6" wide compartments @ 16½" vert. clearance.
- Nominal size: 24" x 18" x 36"h.
- Shipping weight: 100 lbs.

**ALMOND is our new default paint color!**



## Top Loading CR-Cart

Cat.# CRC-46  
\$794.00

Cat.# CRC-45  
\$794.00

Two compartments:  
8"w x 24"l x 13"vert.

Four compartments:  
6"w x 16"l x 13"vert.

- Nominal size: 24" x 16" x 36"h.
- Shipping weight: 90 lbs.



- All steel construction with durable powder coat finish.
- Space saving carts have sturdy 4" plate casters.
- All units have rubber mat inside each compartment.
- Carts have padded bar-style handles and scuff resistant tops.

## Standard Cassette & Film Carts

All "CRM" Cassette & Filing Carts have corrugated mat on the shelves to prolong the life of film cassettes and filing jackets.



### Tall Cassette Cart



Cat.# CRM-38 \$680.00

- Eight 6" wide compartments @ 16<sup>3</sup>/<sub>4</sub>" vertical clearance.
- Nominal size: 24" x 18" x 42<sup>1</sup>/<sub>2</sub>"h.
- Shipping weight: 120 lbs.

### Mini Cassette Cart



Cat.# CRM-24 \$680.00

- Four 6" wide compartments @ 16<sup>3</sup>/<sub>4</sub>" vert. clearance.
- Nominal size: 24" x 18" x 36"h.
- Shipping weight: 95 lbs.



Medical Grade Casters:  
two locking swivel  
two stationary

## CR-Plate & Cassette Covers

- Disposable; prevents cross-contamination.
- Protects technologists, physicians, and patients.
- All varieties are non-sterile.
- 2 mil polyethylene is transparent to x-ray.
- Keeps cassettes clean and dry from exposure to bodily fluids, solutions, plaster, etc.
- Easy to use sandwich bag style folding flap.



Compact Dispenser Box

cat.#	for cass./plate size	quantity	price
CCS-08	8" x 10"	1 box of 100	\$59.00
CCS-10	8" x 10"	1 case of 500	\$268.00
CCM-10	10" x 12"	1 box of 100	\$59.00
CCM-12	10" x 12"	1 case of 500	\$268.00
CCL-14	14" x 17"	1 box of 100	\$59.00
CCL-17	14" x 17"	1 case of 500	\$268.00

14" x 17" also available with ZIPLOCK™ CLOSURE!

- Disposable bag securely seals in CR-plates and cassettes.
- Provides Maximum Protection from cross contamination.



No Image Compromise!



cat.#	for cass./plate size	quantity	price
CCZ-14	14" x 17"	1 box of 100	\$59.00
CCZ-17	14" x 17"	1 case of 500	\$268.00

# 16 Imaging Support Carts (top loading for cr/cassette/film)

- Reinforced wear points.
- All steel construction with durable powder coat finish.
- Not recommended for routine elevator use (see page 40).

- Sturdy 4" plate casters: two stationary, two locking swivel.
- Padded bar-style handle and rubber corner bumpers.
- Rubber mat inside each compartment.

Medical Grade Casters:  
two locking swivel  
two stationary



## Large Activity Cart

Cat.#41AC  
\$680.00

- Nominal size: 24" x 24" x 36"h.
- Four compartments:  
each is 6" x 24" x 13"vert.
- Shipping weight: 127 lbs.



## Small Activity Cart

Cat.#46AC  
\$680.00

- Nominal size: 24" x 16" x 36"h.
- Two compartments:  
each is 8"w x 24"l x 13"vert.
- Shipping weight: 85 lbs.



Convenient  
POP-UP Handle!

## Short Activity Cart

Cat.#47AC  
\$680.00

- Ideal for use while seated.
- Nominal size: 24" x 16" x 27"h.
- Two compartments:  
each is 8"w x 24"l x 13"vert.
- Shipping weight: 82 lbs.

(shown with handle UP)



All Toploaders ship:



- Slim-Lines are designed for cassettes, folders, caddies, film.
- Extra heavy 14 gauge steel, all welded construction.
- Large 5" plate casters (two are locking).

Slim-Line Carts ship:



- Slanted shelves with impact resistant stainless steel edges.
- Rubber mat in each compartment & protected wear points.
- Rubber corner bumpers and dual padded handles.

## two-tier open top

Cat.#SL-1  
\$794.00

- Open top.
- Two tier.
- Two dividers per shelf on 6 1/2" centers.
- Lower shelf @ 15 3/4" vertical clearance.
- Nominal size: 20"w x 17"d x 41"h.
- Shipping weight: 127 lbs.



## Slim Line Carts

## two-tier top shelf

Cat.#SL-2  
\$850.00

- Shelf top.
- Two tier.
- Lower shelf @ 15 3/4" vertical clearance, upper @ 15".
- Two dividers per shelf on 6 1/2" centers.
- Nominal size: 20"w x 17"d x 41"h.
- Shipping weight: 130 lbs.



Padded Handles

## three-tier open top

Cat.#SL-3  
\$1080.00

- Open top.
- Three tier.
- Lower shelves @ 15 3/4" vertical clearance.
- Four dividers per shelf on 4" centers.
- Nominal size: 20"w x 17"d x 49"h.
- Shipping weight: 160 lbs.



# 18 Imaging Support Carts (side loading for film and cassettes)

- All steel construction with durable powder coat finish.
- Sturdy 4" plate casters: two stationary, two locking swivel.
- Padded bar-style handle and scuff resistant top.
- Rubber mat inside each compartment.



Cat.#412C \$749.00

## Cassette / Film Cart (see page 35 for CR-Plate carts)

- Nominal size (not incl. handle): 36" x 18" x 40"h.
- Top shelf: six 6" wide compartments @ 15<sup>3</sup>/<sub>8</sub>" shelf clearance.
- Bottom shelf: four 9" wide compartments @ 15<sup>3</sup>/<sub>8</sub>" shelf clearance.
- Lower shelf is modified to accept 14" x 17" (35 x 43cm) cass./plates.
- Shipping weight: 135 lbs.

Made for  
TALLER FILES.  
Shelf height  
clearance  
15<sup>3</sup>/<sub>8</sub>".



Cat.#415MV \$908.00

## Organization Cart

- Nominal size (not incl. handle): 36" x 18" x 40"h.
- Sixteen 4<sup>1</sup>/<sub>2</sub>" wide compartments @ 15<sup>3</sup>/<sub>8</sub>" shelf clearance.
- No single large compartment.
- Shipping weight: 167 lbs.



NEW!  
Medical Grade Casters:  
two locking swivel  
two stationary

412C  
with locking  
DOORS!



Cat# DR412C \$1049.00

415MV  
with locking  
DOORS!



Cat# DR415MV \$1249.00

- Locks offer attractive security.
- Each door locks independently.
- Flush-mounted paddle handles.
- Heavy duty construction.

Comply with HIPAA  
Privacy Regulations

- All Transporters will withstand elevators and rugged terrain.
- All steel construction with durable powder coat finish.
- 5" industrial casters: two stationary, two swivel.
- Easy movement of heavy film and cassette loads.
- Padded bar-style handle adds 3" to cart length.
- Assembled and ready for years of use.

All Transporters ship:  




Cat.# HD-838 \$1195.00

- 102 linear filing inches.
- Nominal cart dimensions: 36" x 18" x 57"h.
- Each compartment: 6" x 18" x 15<sup>3</sup>/<sub>8</sub>" vert. clearance.
- Shipping weight: 150 lbs.

## Heavy-Duty Transporters




Cat.# HD-460 \$794.00

- 38 linear filing inches.
- Nominal cart dimensions: 24" x 18" x 35"h.
- 2 upper bins @ 9" x 24" x 13"vert.
- 4 lower compartments: 6" x 18" x 14<sup>3</sup>/<sub>4</sub>" vert. clearance.
- Shipping weight: 110 lbs.

## Standard Utility Cart

Cat.# UC-70 \$257.00




- All-purpose utility cart. Ideal for sponges, linens, accessory items, etc.
- 150 pound capacity.
- Unassembled.
- 25" x 19" x 30"h. 
- Shipping weight: 36 lbs.

## Stainless Steel Utility Cart

Cat.# BCS-73 \$349.00



- All-purpose utility cart. Ideal for sponges, linens, accessory items, etc.
  - 300 pound capacity.
  - Pre-assembled.
  - 24" x 15" x 32"h.
  - Shipping weight: 36 lbs. 
- UPS or FedEx may add significant charges to ship this product. Call for a freight quote.