


Prices are not current in this File

Safelights

- Safelights allow a technologist to "see" what they're doing under darkroom conditions without destroying valuable film stores.
- Increases darkroom efficiency. 



safe-spot
Cat.# SSP-15 \$179.00

- Adjustable swivel & pivot base.
- Safe for high speed film.
- Filter is GBX-2 compatible
- Pull chain switch.
- 15-watt incandescent lamp.
- Shipping weight: 3 lbs.

Replacement Filter (Dual Round)
Cat.# SFR-45 *GBX-2 compatible*
\$88.00/ea.
(Interior bulb can be replaced with 15-watt incandescent bulb.)



fluorescent safelight: single lamp model
Cat.# SLS-15 \$225.00

- Ceiling or wall mount.
- 15-watt, 18" fluorescent lamp.
- Blue/Green film type filter will not dry or crack with age.
- Pull chain switch and 6 foot power cord.
- Shipping weight: 8 lbs.

Replacement Filter (Fluorescent)
Cat.# SFT-40 *GBX-2 compatible*
\$117.00/ea.
(White bulb can be replaced with any 15-watt cool white or warm white fluorescent tube.)



mammo-safelight
Cat.# LM-09 \$270.00

- Swivel base for reflective light.
- Specially designed for mammo film.
- Filter is mammography compatible.
- In-line switch in 6 foot power cord.
- Special intensity fluorescent lamp.
- Shipping weight: 5 lbs.

Replacement Filter (Mammography)
Cat.# SFM-47 *mammo compatible*
\$129.00/ea.
Replacement Lamp (Mammography)
Cat.# SLB-48
\$44.00/ea.

"IN USE" Door Signs



Illuminated Door Signs

- Lighted sign focuses attention on exam room status.
- Enhances Radiation Safety by suggesting the risk of accidental entry into an active exam room.
- NOT a plug-in sign; must be hardwired @ 120VAC.
- Each sign measures 11 1/2" x 2 1/2" x 7 1/4" h.
- Shipping weight of each sign: 4 lbs.

English

cat.#	description	price
SLU-XR	X-RAY IN USE	\$169.00
SLU-DK	DARKROOM IN USE	\$169.00
SLU-CT	CT IN USE	\$169.00
SLU-MR	MRI IN USE	\$169.00

Spanish

cat.#	description	price
SLU-XS	X-RAY IN USE	\$169.00
SLU-DS	DARKROOM IN USE	\$169.00



44 General Accessories (general information signs)

Silk Screened Signs

- All signs are constructed of 1/8" thick Trovicel™.
- Double-faced tape is attached to the back of each sign for mounting.
- Large Signs are 8" x 10".
- Small Signs are 4" x 10".
- The large MRI "Danger" Signs are 12" x 18".



Cat.# NFD-22B
\$19.00




Cat.# NFD-22R
\$19.00

NEW!

English


NEW!




Cat.# NCP-24B \$13.00



Cat.# NCP-24R \$13.00



Cat.# PS-04EB \$19.00




Cat.# PS-04ER \$19.00




Cat.#CXR-13EB \$13.00




Cat.#CXR-13ER \$13.00



Cat.# PS-08B \$19.00



Cat.# PS-08R \$19.00



Cat.# XRR-19B \$13.00




Cat.# XRR-19R \$13.00



Cat.# CDK-11B \$13.00



Cat.# CDK-11R \$13.00



Cat.# VS-10EB \$19.00



Cat.# VS-10ER \$19.00


Spanish




Cat.# CXR-13SB \$13.00



Cat.# CXR-13SR \$13.00



Cat.# PS-04SB \$19.00




Cat.# PS-04SR \$19.00




Cat.# CDR-12B \$13.00



Cat.# CDR-12R \$13.00



Cat.# PS-08B \$19.00



Cat.# PS-08R \$19.00

- Signs give critical information to patients & technologists.
- Many states adopt certain signs as law.
- Check your local laws or program requirements.

- Enhances safety with a front line written warning.
- Part of a good quality assurance program.
- Custom signs available for quote (best price with multiples).



Caution



Cat.# CRA-15 \$19.00



Cat.# CRA-17 \$13.00



Cat.# CBH-14 \$13.00




Cat.# CRE-18 \$13.00

MRI



Cat.# CMF-01 \$19.00



Cat.# CMH-06 \$27.00



Cat.# MBN-13E \$50.00



Cat.# MBN-14S \$50.00



Cat.# MBN-15E \$30.00



Cat.#MBN-16S \$30.00

Pricing shown is outdated and not valid

- Large Signs are 8" x 10".
- Small Signs are 4" x 10".
- The large MRI "Danger" signs are 12" x 18".

Custom Engraved Signs

- Engraved 1/16" thick Rowmark™.
- Stock colors: White letters in Red plate, White letters in Black plate, Red letters in Yellow plate.
- Double-faced tape attached to the back for mounting.
- Letters and numbers only; no symbols.
- Other unusual sizes and styles available for quote.

Example: 5" x 9" = 45 sq. in.

Cat.# SI-48 @ \$32.00



Example Sign

size	cat.#	price
48 sq. inches or less	SI-48	\$32.00
over 48 sq. in. and up to 144 sq. in.	SI-144	\$38.00