

52 General Accessories (stainless steel MRI furniture)

MRI Furniture

- Each piece is non-magnetic and safe to use in the MRI environment.
- Some weights do not reflect OverSize UPS/FedEx Rate.



= MR Conditional (no hazards up to 3-Tesla);



= MR Safe (no known hazards);



= MR Unsafe (danger/known hazard).

MRI Office Chair

Cat.# MBC-50CC
\$1395.00



New !

- Offers comfort and versatility for MRI console, MRI room, or any other area within the MRI department.
- Features a unique non-magnetic gas lift mechanism.
- Seat height range: 19"-to-24".
- Back height 27".
- Five-point base with ez-glide casters.
- Weight: 49 lbs.;
- ships @ 120 lbs, UPS/FedEx.



MRI Comfort Stool

Cat.# MBC-40CC
\$795.00



New !

- The most comfortable task stool available, even for extended periods of time.
- Features a unique non-magnetic gas lift mechanism.
- Seat height range: 19"-to-24".
- Padded seat is 16" x 4" thick.
- Five-point base with ez-glide casters.
- Weight: 22 lbs.;
- ships @ 37 lbs, UPS/FedEx.



MRI Step Stool w/ Handrail

- Stainless steel construction.
- **500#** weight capacity.
- Step: 16"w x 12"d x 8 $\frac{1}{2}$ "h.
- Handrail height: 41".
- "MRI Green" mat; rubber tips.
- Ships @ 13 lbs (+oversize) UPS/FedEx.



MRI 2-Step Stool Wide w/ Handrail

Cat.# MBD-68
\$999.00

- Stainless steel.
- **400#** weight capacity.
- Top Step: 36"w x 12"d x 16"h, Bottom Step: 36"w x 12"d x 8"h.
- Handrail height: 41".
- "MRI Green" mat; rubber tips.
- Ships @ 48 lbs, motor freight.



MRI Two-Step Stool Narrow


- Stainless steel construction.
- **500#** weight capacity.
- Top Step: 16"w x 12"d x 16 $\frac{1}{4}$ "h, Bottom Step: 16"w x 12"d x 8 $\frac{1}{2}$ "h.
- Handrail height: 48 $\frac{3}{4}$ ".
- "MRI Green" mat; rubber tips.
- Shipping weight: stool 38 lbs.; handrail 9 lbs.



MRI Platform Stool

Cat.# MBS-60
\$570.00

- Stainless steel construction.
- **550#** weight capacity.
- Step: 30"w x 12"d x 8"h.
- "MRI Green" mat; rubber tips.
- Ships @ 16 lbs, UPS/FedEx.



cat.#	description	price
MBS-65	stool w/handrail	\$325.00
MBS-66	step stool only	\$260.00

cat.#	description	price
MBD-63	stool w/handrail	\$810.00
MBD-64	step stool only	\$750.00

Pricing is not current

MRI Chairs & Stools



- 4" padded seat.
- Screw-post height adjustment.
- Seat upholstered in "MRI Green".
- Two styles of five-leg bases: 2" swivel casters, or stationary rubber tips.
- Ships via UPS/FedEx.



cat.#	height range	features	base	ship wt.	price
CHAIRS					
MBC-30CC	17"-23"	chair only	casters	22	\$710.00
MBC-31CC	17"-23"	chair+back	casters	24	\$845.00
MBC-30CT	15"-21"	chair only	tips	22	\$625.00
MBC-31CT	15"-21"	chair+back	tips	24	\$760.00
STOOLS					
MBC-30SC	23"-29"	stool only	casters	22	\$755.00
MBC-31SC	23"-29"	stool+back	casters	24	\$880.00
MBC-30ST	21"-27"	stool only	tips	23	\$660.00
MBC-31ST	21"-27"	stool+back	tips	25	\$790.00

MRI Wheel Chairs



Standard Wheel Chairs

Bariatric Wheel Chair

- 24" diameter wheels, 8" front casters.
- Seat upholstered in "MRI Green".
- Smooth ride, rust-proof hand rims.
- Retractable wheel locks.
- Quick release seat belt.
- Sling style seat (bariatric: solid padded seat).
- Ships via UPS/FedEx (bariatric: motor freight).



cat.#	seat width overall	weight cap.	ship wt.	price
MBW-18	18"	26"	250#	53 \$2225.00
MBW-20	20"	28"	250#	55 \$2975.00
MBW-22	22"	30"	350#	62 \$3315.00
MBW-26B	26"	35"	550# bari.	67 & 13 \$4680.00
MBW-10	I.V.Pole & O ₂ E-tank mount		7	\$290.00



MRI Transport Gurney

Cat.# MBG-6SS
\$5500.00

- 600# weight capacity.
- Stainless steel construction.
- Features flip-down side rails.
- Backrest adjusts from 0° to 75°.
- Non-marring bumper pads.
- 70"l x 30"w x 39"h.
- Four locking casters.
- Ships @ 300 lbs, motor freight.



MRI PVC Transport Gurney

350# CAPACITY **600# CAPACITY**
Cat.# MBG-3PV Cat.# MBG-6PV
\$1890.00 \$2140.00


- Economic PVC construction.
- Two weight capacities.
- Features flip-down side rails.
- Three position backrest.
- 79"l x 28"w x 33 1/2".
- Four locking casters, two stationary.
- Either ships @ 150 lbs, motor freight.



Transfer Boards

- Facilitates safe and efficient transfer of patients.
- Semi-flexible, slick material makes positioning under the patient easy.
- Durable, easy to clean, white high-density polyethylene.
- Ideal for general imaging, MRI, and CT.
- Transfer boards are radiolucent can be used in the imaging field.
- Conveniently located 1 1/2" x 4 1/2" hand-holds allow optimal movement.

New!



cat.#	description	price
TBS-2372	standard	185.00
TBN-1872	narrow	165.00
TBB-3272	bariatric	235.00
TBH-2332	half	110.00



MRI IV Pole

- Extends from 48" to 92 1/2".
- 24" diameter base.
- 2" swivel casters.
- Ships via UPS/FedEx:
Pole @ 3 lbs,
Base @ 7 lbs.



cat.#	description	price
MBI-70	2-hook iv pole	\$379.00
MBI-75	4-hook iv pole	\$395.00



MRI Conditional Stethoscope

Cat.# MBP-ST
\$330.00

- Equipped with three sound collecting bells suitable for auscultation of infants, children and adults.
- Optimal sensitivity in a wide range of frequencies.
- CE Certified and latex-free.
- Includes multi-pocket bag.
- Ships @ 4 lbs, UPS/FedEx.



MRI Conditional Sphygmomanometer

Cat.# MBP-BP
\$450.00

- Strong with triple mounted shockproof movement.
- Impact resistant.
- Adult cuff (infant cuff available).
- Latex-free/mercury-free.
- Ships @ 2 lbs, UPS/FedEx.



MRI Personal Products

Earplugs: 29db

Cat.# MBE-29
\$170.00



- Economical protection.
- Latex-free.
- Use w/ headset or alone.
- Bulk packaged as 1000 pcs.



Eye Shades

Cat.# MBY-MK
\$210.00



- Shields patients from visual distractions.
- Convenient and disposable.
- 100 individually wrapped shades per bag.



Prism Glasses

Cat.# MBY-PG
\$270.00



- Allows patients to see room outside of magnet during scan.
- Used alone or with sound system to keep patient relaxed.



54 General Accessories (facility MRI accessories)

Metal Detector

Cat.# MBM-MD
\$370.00

- Save patients and staff from embarrassing and hazardous situations.
- Features large scan surface with sensitivity reduction switch.
- Ships @ 15 lbs, UPS/FedEx.



Screening Magnet

Cat.# MBM-SM
\$225.00

- For screening devices and equipment for magnetic properties.
- Equivalent strength of a 0.7-T MRI.
- Lightweight, yet thorough and compact.
- Ships @ 2 lbs, UPS/FedEx.



MRI Illuminated Sign

Cat.# SLU-MR
\$169.00



- Focuses attention on room status.
- Suggests risk of accidental entry.
- Must be permanently mounted and hardwired.
- Measures 11 1/2" x 2 1/2" x 7 1/4".
- Ships @ 4 lbs, UPS/FedEx.

MRI Carpet Floor Mat

Cat.# MBN-12C
\$350.00


- 3ft x 5ft bold warning.
- Rubberized backing prevents slippage.
- Ships @ 9 lbs, UPS/FedEx.



MRI Warning Signs

- Durable plastic silk-screened signs.
- Double-faced tape is on the back for mounting.
- Additional signs are available on pages 44 & 45.

cat.#	sign #	description	size	price
CMR-02	A	yellow, 2 symbols	8"x10"	\$19.00
CMF-01	B	yellow, 1 symbol	8"x10"	\$19.00
CMM-06	C	yellow, 6 symbols	8"x10"	\$27.00
MBN-13E	D	restricted access (Eng)	12"x18"	\$50.00
MBN-14S	E	restricted access (Sp)	12"x18"	\$50.00
MBN-15E	F	mag. always on (Eng)	4"x10"	\$30.00
MBN-16S	G	mag. always on (Sp)	4"x10"	\$30.00




MRI Zone Signs

- Complies with ACR safe practice guidelines for identifying visible safeguards and boundaries to help protect patients, staff and anyone else near or within the MRI department.
- Durable plastic silk-screened signs.
- Double-faced tape is on the back for mounting.

cat.#	sign #	description	size	price
MBZ-FD	H	Zone General Diagram	8"x12"	\$87.00
MBZ-01	I	Zone I Sign	8"x12"	\$49.00
MBZ-02	J	Zone II Sign	8"x12"	\$49.00
MBZ-03	K	Zone III Sign	8"x12"	\$49.00
MBZ-04	L	Zone IV Sign	8"x12"	\$49.00

New !



Diagrams the four MRI ZONES in a generic format. Explains "who" is allowed and restricted "where" and how signs are posted.

pricing is not current

44 General Accessories (general information signs)

Silk Screened Signs

- All signs are constructed of 1/8" thick Trovicel™.
- Double-faced tape is attached to the back of each sign for mounting.
- Large Signs are 8" x 10".
- Small Signs are 4" x 10".
- The large MRI "Danger" Signs are 12" x 18".



Cat.# NFD-22B
\$19.00




Cat.# NFD-22R
\$19.00

NEW!

English


NEW!




Cat.# NCP-24B \$13.00



Cat.# NCP-24R \$13.00



Cat.# PS-04EB \$19.00




Cat.# PS-04ER \$19.00




Cat.#CXR-13EB \$13.00




Cat.#CXR-13ER \$13.00



Cat.# PS-08B \$19.00




Cat.# PS-08R \$19.00



Cat.# XRR-19B \$13.00



Cat.# XRR-19R \$13.00



Cat.# VS-10EB \$19.00



Cat.# VS-10ER \$19.00



Cat.# CDK-11B \$13.00



Cat.# CDK-11R \$13.00


Spanish




Cat.# CXR-13SB \$13.00



Cat.# CXR-13SR \$13.00



Cat.# PS-04SB \$19.00




Cat.# PS-04SR \$19.00




Cat.# CDR-12B \$13.00



Cat.# CDR-12R \$13.00



Cat.# PS-08B \$19.00



Cat.# PS-08R \$19.00

pricing is not current