

Mammography Marker Sets



Screening ACR Mammo Set

(standard set shown)

ITEM	Cat.#	PRICE
• Standard Set	ACR-525	\$160.00
• Velcro™ Set	ACR-526	\$175.00
• Suction Cup Set	ACR-527	\$190.00
• Magnetic Set*	ACR-528	\$185.00



Basic ACR Mammo Set

(suction cup set shown)


ITEM	Cat.#	PRICE
• Standard Set	ACR-510	\$295.00
• Velcro™ Set	ACR-511	\$300.00
• Suction Cup Set	ACR-512	\$300.00
• Magnetic Set*	ACR-513	\$300.00

ACR STANDARDIZED ABBREVIATIONS

Sets are designed respective of the ACR abbreviations.

PREFIX (7/16" letters)	R (Right) L (Left) M (Magnification)
BASE (1/4" letters)	MLO (Mediolateral Oblique) CC (Craniocaudal) ML (Mediolateral) LM (Lateromedial) XCCL (Exaggerated Craniocaudal) CV (Cleavage) AT (Axillary Tail) TAN (Tangential) ID (Implant Displaced) SIO (Superolateral to Inferomedial Oblique) LMO (Lateromedial Oblique) FB (From Below)
SUFFIX (1/4" letters)	RL (Roll Lateral) RM (Roll Medial) RS (Roll Superior) RI (Roll Inferior)
DEGREES (1/4" letters)	(no standard)
INITIALS (1/4" letters)	(no standard)

*Magnetic Sets include markers with magnetic backing that stick to a white metal wallboard which is attached with Velcro™ to the wall or equipment. Can be mounted horizontally or vertically.



*Magnetic Sets



ALL NEW
MAGNETIC
MAMMO SETS!


NEW COLOR!
Mammography Pink!

Mammography Markers custom or replacement

Any marker can be created with tech's initials!

		STANDARD		VELCRO™		SUCTION CUP		MAGNETIC*	
description	example	cat.#	price	cat.#	price	cat.#	price	cat.#	price
single letter	R	CPT-01	\$13.00	CPT-01V	\$14.00	CPT-01S	\$15.00	CPT-01M	\$14.00
2-3 letters	R ML	CPT-23	\$14.00	CPT-23V	\$15.00	CPT-23S	\$16.00	CPT-23M	\$15.00
4-5 letters	L 50°	CPT-45	\$17.00	CPT-45V	\$18.00	CPT-45S	\$19.00	CPT-45M	\$18.00
6-7 letters	R CC MBT	CPT-67	\$20.00	CPT-67V	\$21.00	CPT-67S	\$22.00	CPT-67M	\$21.00
8+ letters	L XCCL MBT	CPT-80	\$24.00	CPT-80V	\$25.00	CPT-80S	\$26.00	CPT-80M	\$25.00

CPT-67 with initials



Cat.#CPT-45



Cat.#CPT-23



- Same style as those on the ACR sets shown above.
- Aluminum-backed, high density letters.

- R & L are 7/16" - all other letters are 1/4".
- Comes with a hole for hanging (except magnetic version).

Economy Nipple-Artifact Markers

Cat.# LNA-20
\$50.00



- Encapsulated 2mm LEAD ball.
- Used to mark nipple, lesion, scar, artifact, etc.
- Neutral color, peel-off, self-adhesive backing.
- Non-sterile.
- 100 markers mounted on ten convenient strips.

pricing is not current