


## Transfer Boards **NEW!**


- Facilitates safe and efficient transfer of patients.
- Semi-flexible, slick material makes positioning under the patient easy.
- Ideal for general imaging, MRI, and CT.
- Transfer boards are radiolucent can be used in the imaging field.
- Conveniently located 1 1/2" x 4 1/2" hand-holds allow optimal movement.

Your Choice  
of  
Four Sizes!




### standard board

Cat.# TBS-2372  
23" w x 72" l x 3/16" thick  
1-4 boards: \$185.00 ea.  
5-9 boards: \$160.00 ea.  
10+ boards: \$135.00 ea.



### narrow board

Cat.# TBN-1872  
18" w x 72" l x 3/16" thick  
1-4 boards: \$165.00 ea.  
5-9 boards: \$150.00 ea.  
10+ boards: \$130.00 ea.



### bariatric board

Cat.# TBS-3272  
32" w x 72" l x 1/4" thick  
\$235.00 ea.



### half board

Cat.# TBH-2332  
23" w x 32" l x 3/16" thick  
\$90.00 ea.



Cat.# TBR-2P  
\$53.00

### WALL STORAGE RACK

- Two 6" pegs are angled to keep board from slipping.
- Mounting holes are on vertical centers.
- Steel with white enamel.
- Ships @ 4 lb

- Regular Boards are made of durable HDPE (high density polyethylene).
- Anti-static boards are made of UHMW (ultra high molecular weight polyethylene).
- Color: regular boards = **white** / anti-static boards = **black**.
- Due to size, transfer boards ship at OverSize rates.

### ANTI-STATIC TRANSFER BOARDS

cat.#	description	price
TBS-2372AS	standard	255.00
TBN-1872AS	narrow	245.00
TBB-3272AS	bariatric	575.00
TBH-2332AS	half	110.00

Pricing is not current