

Step Stool with Handrail

Cat.# BSH-65
\$99.00

- All steel 350 lb. capacity.
- Skid resistant black mat.
- 9" high stool with 11" x 14" top; 34" taal with handrail.
- Handrail is 7/8" dia. chrome-plated tubular steel.
- Easy assembly.
- Shipping weight: stool 8 lbs.; handrail 4 lbs.



Step Stool without Handrail



Cat.# BSS-60
\$67.00

- Similar to stool with BSH-65.
- NO HOLES for handrail.
- Shipping weight: 8 lbs.

Bariatric Step Stools

- Bariatric 600 lb. capacity.
- Steps have a skid resistant black mat surface.
- Height:
Single Step: 8³/₄" high step; 34¹/₂" with handrail.
Double Step: 8" & 16" high steps; 43" with handrail.
- Handrail is 1" dia. chrome-plated 16g steel.
- Easy assembly.
- Shipping weight:
Single Step: stool 8 lbs.; handrail 6 lbs.
Double Step: stool 20 lbs.; handrail 5 lbs.



Cat.# BBS-61
\$165.00




Cat.# BBD-63
\$336.00



Cat.# BBS-62
\$124.00



Cat.# BBD-64
\$269.00



IV Stands

- 1" and 7/8" telescoping, chromed steel tubes.
- Adjustable height range: 47" to 84"; Twist Lock.
- Coated, bar steel base with scuff guard tips resists corrosives.
- Two or four rams-horn style IV hooks.
- Weight: 17 lbs.

cat.#	description	price
BIV-70	2-hook iv pole	\$99.00
BIV-75	4-hook iv pole	\$186.00

Waste Receptacle

- 32-quart waste cans with foot operated lid.
- Use Biohazard Bags as safety liners.
- Easy to clean plastic shell and lid.
- Available in White and Red (includes biohazard label).
- 15¹/₂w x 14"d x 19" h.
- Weight: 8 lbs.

cat.#	description	price
BPT-6W	White Waste Can	\$95.00
BPT-6R	Red Waste Can	\$95.00



Biohazard Bags

- High density, 12 microns.
- 10 gallon (24" x 24").
- 250 bags per box.

Cat.# NTB-30
\$58.00/bx



48 General Accessories (stools, instrument stand & lamps)



Screw Adjustable Stool

Cat.# BAS-30
\$199.00



- Solid steel, machined thread screw adjusts height by rotation of the stool cushion.
- Adjustable height range: 18" to 24½".
- 14½" diameter seat cushion; 4" thick (black).
- 4-leg, chrome plated base of 14 gauge steel.
- Tubular steel footring for comfort and stability.
- Weight: 20 lbs.; Ships @ 30 lbs.

Pneumatic Stool

Cat.# BPS-20
\$385.00



- Pneumatic air cylinder easily adjusts height with convenient under-cushion ring.
- Adjustable height range: 17½" to 22½".
- 16" diameter seat cushion; 4" thick (black).
- 5-leg, die cast aluminum base.
- Base is 22" diameter for excellent stability.
- Weight: 21 lbs.; Ships @ 30 lbs.

Cat.# BRF-45
\$17.00



- Set of four (4) rubber tips for BAS-30 stool ONLY.
- Also fits standard and bariatric step stools.



Instrument Stand (Mayo Stand)


Cat.# BMS-80
\$220.00

- Legs and base are 1" chrome plated steel.
- Adjustable height range: 34" to 53".
- Solid steel, easy-slide extension posts.
- Stainless steel tray is 12" x 19" x ¾" deep.
- Weight: 12 lbs.; Ships @ 30 lbs.

Standard Exam Lamp

- Flexible gooseneck arm with 360° adjustable steel shade.
- Adjustable height range: 50" to 74".
- Chrome-plated, weighted steel base.
- Uses one 10 to 40 watt standard bulb (not included).
- Weight: 15 lbs.

Exam Lamp
Cat.#BEL-90
\$155.00



Halogen Lamp
Cat.#BEH-91
\$299.00

Halogen Exam Lamp

- Flexible gooseneck arm with 540° adjustable head.
- Beige coated, weighted base.
- Uses one high intensity 35 watt halogen lamp.
- 1500 Footcandles @ 24"; 3000° Kelvin.
- Weight: 16 lbs.

