

40 General Accessories (mobile cr-plate/cassette holders)

Standard Tilt & Rotate (mobile cr-plate & cassette)

Standard Tilt & Rotate
(complete mobile unit)
Cat.# MCH-33B
\$879.00


- Standard arm length provides maximum stability.
- Fully articulating holder.





Keep it secure
with Posi-lock
height adjustment!

ALL MOBILE UNITS INCLUDE:

- Vertical range from bottom of cassette/plate: floor-to-54".
- Arm slides vertically on 60" chromed steel post.
- Holder locks onto post with spring-loaded pin.
- 18" square steel shelf base with protective mat.
- Heavy base prevents tipping.
- E-Z-Roll swivel casters (two are locking).
- Neutral Almond baked enamel finish.
- Ships in two cartons: 1@15#, 1@70#.




Standard and Extension models include
the tilt & rotate, fully articulating holder.



Extension Arm Tilt & Rotate (not for dr)

Extension Arm Tilt & Rotate
(complete mobile unit)
Cat.# MCH-37B
\$943.00

- Extends plate/cassette 11" to 23" from post.
- Fully articulating holder.



Component Selection


Once you have a unit that includes the mobile base/post, you may choose to add additional holders to your system for a greater range of utility.

Cat.#	Description	Price
MCA-B31B	Mobile Base (no post)	\$355.00
MCA-P31B	Posi-lock Post	\$345.00
MCA-S33B	Standard Holder	\$240.00
MCA-E37B	Extension Arm Holder	\$295.00
MCA-T39B	Tall Cassette Holder	\$280.00
MCA-D40B	DR-Plate Holder	\$240.00

DR-Plate Holder (non-tilt or rotate)

DR-Plate Lateral
(complete mobile unit)
Cat.# MCH-40B
\$882.00

- Holds up to 14" x 17" DR-Plates.
- Fixed position only (non-articulating).



Tall Cassette Holder (non-tilt or rotate)

Tall Cassette Lateral
(complete mobile unit)
Cat.# MCH-39B
\$920.00

- Holds 14" x 36" or 51" cassettes.
- Fixed position only (non-articulating).



Tabletop Lateral Cassette/CR-Plate Holder



Improved!
Wider Channels
are Adjustable from
Base to 24"



24" Variable Height
Plate/Cassette Holder
Cat.# 35
\$216.00

- Eliminates hand held plate/cassette hazards.
- Enhances Technologist's Safety.
- Securely holds plate/cassette in a stable, upright position.
- Both 6" wide channels adjust to fit cassette.
- Holds up to 14" x 17" (35 x 43cm) cassettes.
- Accepts Grid Cassettes and CR-plates (channel is 1 1/4"d).
- Base plate slides under patient or pad.
- Upper channel includes a rotating R & L marker.
- Both Upper & Lower Channels are lined with rubber to keep all types of plates/cassettes in a stable position.
- Shipping weight: 7 lbs.

