



Alcohol Prep Pads & Swabsticks

- Topical antiseptic/anti-infective for preparation of skin prior to injections
- Saturated with 70% isopropyl alcohol
- Sterile alcohol pads gamma sterilized to meet AAMA/ANSI/ISO 11137 standards
- Soft, absorbent lint-free non-woven pad; cotton-tipped swabsticks
- Medium prep pad: 1-1/8 x 2-3/8 inches
- Large prep pad: 1-3/4 x 3-1/4 inches
- Swabsticks: 4 inches