

# NEW 5900 Digital Pediatric Scale

provides ease of use, comfort and accuracy




## ALL STEEL

## MODULAR DESIGN

The modular design provides a comfortable work surface height of 35 inches and plenty of handy storage space.



5900



## ACCURATE MEASUREMENT

Comes standard with stainless steel measuring rule with paddles by SECA® and digital load cell technology by SR Scales®. Both provide accurate measurement and ease of use.

## LASTING QUALITY LOW MAINTENANCE

The UMF 5900 table is designed to provide years of trouble-free use and comes with a 5-year warranty.

### 5900 Digital Pediatric Scale • Features

#### Standard Features

- Seamless upholstered top
- Padded front rail for infant safety
- Paper roll holder for 18" paper
- Digital scale operated by 6 D cell batteries (not included) with automatic shut-off feature
- SECA® stainless steel measuring rule with paddles
- Weight capacity of 50 lbs.
- Weighs in both lbs./oz. and kg.
- Tare weight capability

#### Dimensions

- Overall: 48" W x 44.5" H x 25" D
- Work surface: 35" H
- Drawers: 16.5" W x 17" D x 4" H
- Door opening: 18.875" W x 16" H



*UMF Quality, Design and Innovation*

**EXPECT MORE**

**United Metal Fabricators, Inc.**  
1316 Eisenhower Blvd. • Johnstown, PA 15904  
800-638-5322 • 814-266-1870 fax  
[www.umf-exam.com](http://www.umf-exam.com)