



UMF Quality, Design and Innovation

EXPECT MORE



Modular Cabinetry from UMF
Flexible • Functional

MODULAR CABINETRY

Durable • Functional • All-steel construction

UMF modular cabinetry allows you to be creative in setting up an exam room. Begin the process by selecting from the many different base cabinet configurations offered. UMF will be happy to assist you in designing your exam room and can custom manufacture cabinetry for you if our standard sizes and/or configurations do not match your needs.


Our base cabinets are 36 inches high with 19 inch deep countertops that have a 4 inch high backsplash. One-piece continuous countertops can be fabricated up to 12 feet in length. All cabinetry can be fitted with optional locks. Changing drawer and door fronts is quick and easy on UMF modular cabinets ~ a feature available only from UMF.



6002
30" Base Cabinet
with four drawers



6004
24" Base Cabinet
with four drawers



6006
18" Base Cabinet
with four drawers

Cabinet body, door and drawer fronts are 20 gauge.

Cabinet shelves are 20 gauge and adjustable in half-inch increments.

All drawers have 16 inch full extension drawer slides.

All cabinetry can be fitted with optional locks.

-501 single door

-502 double door

-503 sliding door

-505 single drawer



6010
36" Base Cabinet
with two doors and two
adjustable shelves



6014
24" Base Cabinet
with two drawers
and two doors



6017
18" Base Cabinet
with two drawers
and one door

6012 30" Base Cabinet with two drawers and two doors is also available

BASE CABINETS WITH SINKS

Available in the most useful sizes, configurations

You can easily add one or more sinks to any room. UMF base cabinets with sinks come in the most useful sizes and configurations. Sinks are available in 24, 30 or 36-inch base cabinets or our popular 30 x 30-inch corner cabinet. All come standard with a stainless steel sink, stopper basket and chrome plated gooseneck faucet with wrist blade handles.



6024
24" Base Cabinet
with sink



6048
Corner Base Cabinet
with sink

6022 • 30" Base Cabinet with Sink

- Two doors
- Stainless steel sink & faucet

6023 • 36" Base Cabinet with Sink

- Two doors
- Stainless steel sink & faucet

6024 • 24" Base Cabinet with Sink

- Two doors
- Stainless steel sink & faucet

6042 • Corner Base Cabinet

- One bi-fold door hinged right
- One adjustable shelf
- No sink
- Top size 30" x 30"

6048 • Corner Base Cabinet with Sink

- One bi-fold door hinged right
- Stainless steel sink & faucet
- Top size 30" x 30"

Standard UMF sink size:

12" W x 10" L x 5" D

30.48 cm W x 25.40 cm H x 12.80 cm D

MODULAR WALL CABINETS

Easy to install • All-steel construction

UMF wall cabinets are designed to complement the modular base cabinets and are 20 inches high and 11 inches deep. The cabinets are manufactured with a welded on “J” hook and are provided with wall hanger strips with holes at 1/2-inch intervals for ease of installation and leveling of cabinets. All wall cabinets have one adjustable shelf and concealed hinges that can be easily adjusted.



Modular wall cabinets are available in the following sizes:

6070 • 36” Wall Cabinet

6072 • 30” Wall Cabinet

6074 • 24” Wall Cabinet

6076 • 18” Wall Cabinet

6082 • 30” x 30” Corner Wall Cabinet



“J” hook and wall hanger strips with holes at 1/2 inch intervals make installation easy.

MODULAR DESK UNITS

Modular desk units are available in the following sizes:

6052 • 30” Desk Unit

6054 • 24” Desk Unit

- Predrilled to mount at 31 inches high
- Mounting brackets included



6052

TREATMENT CABINETS

Popular sizes

Designer colors

All-steel construction

Selection, style, value

UMF Treatment Cabinets are available in the most popular sizes and in an array of designer colors. The value is enhanced with all-steel construction. All UMF Treatment Cabinets can be purchased with optional locks on doors and drawers. Surfaces are Kydex in white leather color.



6104



6105



6110

6104 • Treatment Cabinet

- One removable non-adjustable metal shelf
- Leveling glides
- Kydex work surface

Options:

- 501 Single door lock
- 505 Single drawer lock

Overall dimensions:

20" W x 34" H x 16.25" D

50.8 cm W x 86.36 cm H x 41.28 cm D

6105 • Treatment Cabinet

- Lock on upper door
- Two adjustable glass shelves in upper section
- One removable non-adjustable metal shelf in bottom cabinet
- Leveling glides
- Kydex work surface

Options:

- 501 Single door lock
- 505 Single drawer lock

Overall dimensions:

20" W x 65" H x 16.25" D

50.8 cm W x 165.1 cm H x 41.28 cm D

6110 • Treatment Cabinet

- One removable non-adjustable metal shelf
- Leveling glides
- Kydex work surface

Options:

- 502 Double door lock
- 505 Single drawer lock

Overall dimensions:

25.25" W x 34" H x 16.25" D

64.14 cm W x 86.36 cm H x 41.28 cm D

TREATMENT CABINETS



6112

6112 • Treatment Cabinet

- Four storage drawers
- Two large doors
- Leveling glides
- Optional locks

Overall dimensions:

32" W x 34" H x 16.25" D

81.28 cm W x 86.36 cm H x 41.28 cm D



6113

6113 • Treatment Cabinet

- Lock on upper door
- One adjustable glass shelf in upper section
- Large convenient storage areas in lower cabinet
- Leveling glides
- Optional locks

Overall dimensions:

32" W x 65" H x 16.25" D

81.28 cm W x 165.1 cm H x 41.28 cm D



6116

6116 • Treatment Cabinet

- One removable non-adjustable metal shelf
- Leveling glides
- Optional locks

Overall dimensions:

32" W x 34" H x 16.25" D

81.28 cm W x 86.36 cm H x 41.28 cm D

TREATMENT CABINETS



6117

6117 • Treatment Cabinet

- Lock on upper door is standard
- One adjustable glass shelf in upper section
- Bottom section has two drawers, two doors with one removable non-adjustable metal shelf
- Leveling glides

Options:

- 502 Double door lock
- 505 Single drawer lock

Overall dimensions:

32" W x 65" H x 16.25" D

82.28 cm W x 165.1 cm H x 41.28 cm D



6120

6120 • Treatment Cabinet

- One removable non-adjustable metal shelf
- Leveling glides

Options:

- 502 Double door lock
- 505 Single drawer lock

Overall dimensions:

27" W x 34" H x 16.25" D

68.58 cm W x 86.36 cm H x 41.28 cm D

TREATMENT CABINETS & MOBILE TREATMENT CABINETS



6130

6130 • Treatment Cabinet

- Lock on upper door is standard
- One adjustable glass shelf in upper section
- Bottom section has one drawer, two doors with one removable non-adjustable shelf
- Leveling glides

Options:

- 502 Double door lock
- 505 Single drawer lock

Overall dimensions:

25.25" W x 65" H x 16.25" D
64.14 cm W x 165.1 cm H x 41.28 cm D



6204

6204 & 6214 • Mobile Treatment Cabinet

- 2-inch twin wheel carpet casters
- Kydex top with rail
- Push handle
- Four drawers (6204)
- Two drawers, two doors (6214)

Overall dimensions:

25" W x 34.75" H x 18" D
63.5 cm W x 88.27 cm H x 45.72 cm D



6620

6620 • EKG Cart

- Steel base with five twin wheel carpet casters
- Guard rail on three sides of Kydex work surface
- Two large drawers
- Stainless steel handle can be used to store leads

Overall dimensions:

20" W x 33.5" H x 16.25" D
50.8 cm W x 85.09 cm H x 41.28 cm D

STORAGE AND SUPPLY CABINETS

Designer colors • All-steel construction



7142 • Storage & Supply Cabinet

- Three adjustable glass shelves in top section
- Lock is standard in upper section
- 110 volt light in upper section automatically lights when door is opened
- One removable, non-adjustable metal shelf in bottom section
- Leveling glides

Options:

- 502 Double door lock
- 505 Single drawer lock

Overall dimensions:

31.5" W x 67" H x 15.5" D

80.1 cm W x 170.18 cm H x 39.37 cm D

Drawer size:

11.25" W x 12.5" L x 4.5" D

28.58 cm W x 31.75 cm H x 11.43 cm D



7165 • Storage & Supply Cabinet

- Three adjustable glass shelves in top section
- One drawer in bottom section
- One removable, non-adjustable metal shelf in bottom section
- Leveling glides
- Lock is standard in upper section

Options:

- 501 Single door lock
- 505 Single drawer lock

Overall dimensions:

19.5" W x 63.5" H x 15.5" D

49.53 cm W x 161.29 cm H x 39.37 cm D

Drawer size:

13" W x 12.5" L x 4.5" D

33.02 cm W x 31.75 cm H x 11.43 cm D

NARCOTICS STORAGE CABINETS

Available in designer (and custom) colors

7780 • Double Door

- 20 gauge double walled outer door / single inner door
- Separate locks on each door, keyed differently
- Three adjustable shelves
- Can be wall mounted or mounted inside a cabinet

Overall dimensions: 18" W x 24" H x 12" D

45.72 cm W x 60.96 cm H x 30.48 cm D

7781 • Double Door

- Separate locks on each door, keyed differently
- Two adjustable shelves
- Can be wall mounted

Overall dimensions: 12" W x 15" H x 8" D

30.40 cm W x 38.10 cm H x 20.32 cm D

7782 • Single Door Single Lock

- One shelf
- Extra wide end allows door to open fully
- Can be wall mounted or mounted inside a cabinet

Overall dimensions: 14" W x 12" H x 8" D

35.56 cm W x 30.48 cm H x 20.32 cm D

7784 • Single Door Single Lock

- No shelf
- Can be wall mounted

Overall dimensions: 12" W x 7.75" H x 7" D

30.48 cm W x 19.69 cm H x 17.78 cm D

7785 • Single Door Double Lock

- Five adjustable shelves
- Wall mounted
- Two locks keyed differently

Overall dimensions: 18" W x 20.25" H x 13.5" D

45.72 cm W x 51.44 cm H x 34.29 cm D

7787 • Single Door Double Lock

- Double walled door
- Two locks, keyed differently
- Wall mount or cabinet mount
- One adjustable shelf

Overall dimensions: 14" W x 12" H x 8" D

35.56 cm W x 30.48 cm H x 20.32 cm D

7788 / 89 • Single Door Double Lock

- No shelf

7788 Overall dimensions: 14" W x 8" H x 10" D

35.56 cm W x 20.32 cm H x 25.40 cm D

7789 Overall dimensions: 12" W x 7.75" H x 7" D

30.48 cm W x 19.69 cm H x 17.78 cm D



7780



7782



7785



7788

Option: -530 **Warning Light**. Available on all Narcotics Storage Cabinets. Factory installed only.

STORAGE AND SUPPLY CABINETS



7811 • Storage & Supply Cabinet

- Floor cabinet offers large amount of storage space
- Leveling glides
- Five removable metal shelves, adjustable on .5" centers

Option: -503 Sliding door lock

Overall dimensions: 47" W x 84" H x 16" D

119.38 cm W x 213.36 cm H x 40.64 cm D

Door opening: 18" W x 76" H

45.72 cm W x 193.04 cm H



7816 • Storage & Supply Cabinet

- Floor cabinet offers large amount of storage space
- Glass doors
- Leveling glides
- Five removable metal shelves, adjustable on .5 inch centers

Option: -503 Sliding door lock

Overall dimensions: 47" W x 84" H x 16" D

119.38 cm W x 213.36 cm H x 40.64 cm D

Door opening: 18" W x 76" H

45.72 cm W x 193.04 cm H



8123 • Utility Table

- One drawer with double walled front
- 2" casters
- Painted steel top and guard rail

Overall dimensions: 20" W x 32" H x 16" D

50.80 cm W x 81.28 cm H x 40.64 cm D



UMF Cabinetry redefines durability, style and functionality.

Designed with health care professionals in mind, the all-steel construction provides a virtual lifetime of trouble free use. UMF's designer colors allow you to coordinate your casework with UMF's entire line of exam tables, seating and accessories.

UMF Quality, Design and Innovation

EXPECT MORE

www.umf-exam.com