

Offering patient comfort and the convenience of an electro-mechanical actuator

- UMF Ultra Comfort Seamless Top
- 3 Year Limited Warranty
- Auto-Position Armrests
- Easy Access

Phlebotomy/ENT chairs from UMF help put your patients at ease while improving efficiency of operation. The UMF Ultra Comfort Seamless Top provides security and comfort for patients and also allows for easy cleaning and helps prevent cross-contamination.



**8678**  
Shown with optional fold-away table and articulating headrest

The stable base, auto-position armrests, and an electro-mechanical actuator make patient positioning easy and convenient. Both models also recline to a nearly horizontal position.

**Standard Features:**

- UMF Ultra Comfort Seamless Top
- 13 Upholstery color options (see page 12)
- Auto-position armrests
- Electro-mechanical actuator system
- Reclines to a nearly horizontal position
- Two function hand control
- 375 pound (170.1 kg) weight capacity

**Table Specifications:**

- Height (at seat):  
8677: 24-36" (61-91.4cm)  
8678: 19-35" (48.3-88.9cm)
- Width:  
31.5" (80cm)
- Seat Width:  
18" (45.7cm)
- Length (fully extended):  
67" (170.2cm)
- Arm Length:  
23" (58.4cm)
- Electrical:  
115 VAC, 60 HZ

**Options:**

- (see page 11 for popular option details)
- 032 Articulating headrest
  - 033 Factory installed fold-away table
  - 067 Simultaneous hand and foot control
  - 068 Foot control
  - 143 Easy-Glide locking casters
  - Fire retardant foam
  - 220 VAC, 50/60 HZ