



ULTRARAY
MEDICAL PRODUCTS INC.

RADIATION PROTECTION & ACCESSORIES

[HOME PAGE](#)

TEL: (905) 338-6857
FAX: (866) 338-8990
TOLL FREE: (877) 338-6857

 **PROUD TO BE CANADIAN MADE** 

Since its inception in 1983 UltraRay® has continued to be an industry leader in the field of radiation protection. We manufacture Canadian made products using the latest technology while maintaining the highest safety standards. Our products are manufactured using the strictest quality control procedures in accordance with Health Canada Medical Devices Regulations. We are a Class 1 Certified manufacturer and distributor operating under our Medical Establishment License 2312.

Our distributor network of 50 sales representatives across Canada provide on-site services such as samples, measurement and needs analysis consultation to ensure that clients receive the products best suited to their needs. Client input is utilized in the ongoing improvement & development of our products.

All our products carry a 2 year Warranty and 100% satisfaction guarantee.

Copyright © 2005-2010 UltraRay Medical Products Inc., All Rights Reserved

FILM/CASSETTE BLOCKERS



Cassette Blockers are made in a variety of sizes and available in .5mm & 1mm lead equivalency. Custom sizes are available upon request.

Coated with washable tan vinyl.

BL-1	Half Film - 5" x 8"	BL-6	Half Film - 5.5" x 14"
BL-2	Half Film - 4" x 10"	BL-7	Half Film - 7" x 15"
BL-3	Half Film - 6" x 10"	BL-8	Half Film - 7" x 17"
BL-4	Half Film - 5" x 12"	BL-9	Half Film - 8.5" x 14"
BL-5	Half Film - 7" x 11"	BL-10	Custom Made Sizes