



[HOME PAGE](#)

RADIATION PROTECTION & ACCESSORIES

TEL: (905) 338-6857
FAX: (866) 338-8990
TOLL FREE: (877) 338-6857

 **PROUD TO BE CANADIAN MADE** 

Since its inception in 1983 UltraRay® has continued to be an industry leader in the field of radiation protection. We manufacture Canadian made products using the latest technology while maintaining the highest safety standards. Our products are manufactured using the strictest quality control procedures in accordance with Health Canada Medical Devices Regulations. We are a Class 1 Certified manufacturer and distributor operating under our Medical Establishment License 2312.

Our distributor network of 50 sales representatives across Canada provide on-site services such as samples, measurement and needs analysis consultation to ensure that clients receive the products best suited to their needs. Client input is utilized in the ongoing improvement & development of our products.

All our products carry a 2 year Warranty and 100% satisfaction guarantee.

Copyright © 2005-2010 UltraRay Medical Products Inc., All Rights Reserved

PATIENT TRANSFER BOARDS



The Slider Board is manufactured from 1/4" PVC plastic and is used to slide a patient from a stretcher onto a bed or vice versa. Handles are conveniently located for access from all sides and a hole is drilled at the top for hanging when not in use.

SL-1 SLIDER BOARD (2'wide x 6.5' long)

SL-C CUSTOM SLIDER BOARD - ***YOU SPECIFY SIZE***