



ULTRARAY
MEDICAL PRODUCTS INC.

RADIATION PROTECTION & ACCESSORIES

[HOME PAGE](#)

TEL: (905) 338-6857
FAX: (866) 338-8990
TOLL FREE: (877) 338-6857

 **PROUD TO BE CANADIAN MADE** 

Since its inception in 1983 UltraRay® has continued to be an industry leader in the field of radiation protection. We manufacture Canadian made products using the latest technology while maintaining the highest safety standards. Our products are manufactured using the strictest quality control procedures in accordance with Health Canada Medical Devices Regulations. We are a Class 1 Certified manufacturer and distributor operating under our Medical Establishment License 2312.

Our distributor network of 50 sales representatives across Canada provide on-site services such as samples, measurement and needs analysis consultation to ensure that clients receive the products best suited to their needs. Client input is utilized in the ongoing improvement & development of our products.

All our products carry a 2 year Warranty and 100% satisfaction guarantee.

Copyright © 2005-2010 UltraRay Medical Products Inc., All Rights Reserved

APRON STORAGE

Click on item for specifications



[AR-1 Wall Mounted Rack](#)

[One Apron & Pair of Gloves](#)



[AR-5 MOBILE APRON RACK](#)

[5-8 Aprons](#)



[AR-6 WALL MOUNTED RACK](#)

[5 or 10 Aprons](#)



[AR-8 MOBILE SPIDER RACK](#)

[8 Aprons - Porta Rack](#)



[AR-10 MOBILE APRON RACK](#)

[10-15 Aprons](#)



[ARH APRON RACK HANGERS](#)

[Closed or Open Loop](#)



HAR-V HALF APRON RACK

Holds 4 VELCRO® Patient Aprons



HAR-H DEMI RING APRON RACK

Holds 4 Hoop Patient Aprons



AR-1 Apron Rack

Designed for individual room use. Unit is wall mounted and holds a pair of gloves and up to 2 aprons.

Features:

- Hold up to 2 aprons and a pair of gloves
- Chromed Steel
- Easy Wall Mount Installation



AR-5 Mobile Apron Rack

Great portability allowing easy transportation of multiple aprons and shields to any location in the facility.

Features:

- Holds 5-7 Aprons
- Heavy Duty Chromed Steel Tubing
- Steering Enabled Casters
- 24" x 22" Base
- 46" Height
- HANGERS NOT INCLUDED

AR-6R - Right Wall Mounted Apron Rack

AR-6L - Left Wall Mounted Apron Rack

Arms fold flush against the wall and are available in right and left folding models. Dimensions folded – 28"W x 4 1/8"D, arms 90° to wall- 26"D, overall height – 9". 15 lb. Arm length- 27 1/2".

Right and left can be combined to create one large rack for 10 aprons (5 fold right, 5 fold left).

Features:

- Holds 5 Aprons
- Heavy Duty Chromed Steel Tubing





AR-8 Spider Style Porta Rack

Great portability allowing easy transportation of multiple aprons and shields to any location in the facility. Uses 24" of Floor Space. Ideal for O.R. and Cath Labs

Features:

- Holds 8 Aprons
- Heavy Duty Powder Coating
- 5 Leg Safety Design
- 24" Base

Copywrite © 2005-2010 UltraRay Medical Products Inc., All Rights Reserved



AR-10 Mobile Apron Rack

Great portability allowing easy transportation of multiple aprons and shields to any location in the facility.

Features:

- Holds 10-13 Aprons
- Heavy Duty Chromed Steel Tubing
- Steering Enabled Casters
- 36" x 22" Base
- 54" Height
- HANGERS NOT INCLUDED

ARH Apron Hanger Closed Loop

Designed for use with the AR-5, AR-10. Hanger is made of solid steel chrome plated construction.

Also available with open loop as shown.



ARHO