

# 104 Vital Signs Devices

MULTIPARAMETER SPOT CHECK DEVICES AND MONITORS

## Spot Vital Signs® LXi


Take vital signs in about half the time.

- > SureBP™ technology measures blood pressure in about 15 seconds as the cuff is inflating—resulting in time savings for you and improved comfort for your patients; pulse rate and MAP included
- > Choice of SureTemp® Plus or Braun ThermoScan® PRO 4000 ear thermometer
- > Choice of Masimo® or Nellcor® pulse oximetry
- > Ability to manually enter weight, height, respiration rate, and pain level; calculates Body Mass Index
- > Connects to select Healthometer® or Tanita® digital scales
- > 50-reading memory allows retrieval of previous data for review, to print or to download later
- > Comes connectivity-ready, allowing you to connect to the electronic medical record system you have now or one you'll get in the future
- > Error messages give you descriptions and real solutions—not just error codes
- > Bigger, brighter LCD display provides greater visibility from farther distances



**45NT0-E1**  
Spot Vital Signs LXi with SureBP, Nellcor SpO<sub>2</sub>, and SureTemp Plus Thermometer

- 450T0-E1** Spot Vital Signs LXi with SureBP and SureTemp Plus Thermometer
- 450E0-E1** Spot Vital Signs LXi with SureBP and Braun ThermoScan PRO 4000 Ear Thermometer
- 45MT0-E1** Spot Vital Signs LXi with SureBP, Masimo SpO<sub>2</sub>, and SureTemp Plus Thermometer
- 45NT0-E1** Spot Vital Signs LXi with SureBP, Nellcor SpO<sub>2</sub>, and SureTemp Plus Thermometer
- 45ME0-E1** Spot Vital Signs LXi with SureBP, Masimo SpO<sub>2</sub>, and Braun ThermoScan PRO 4000 Ear Thermometer
- 45NE0-E1** Spot Vital Signs LXi with SureBP, Nellcor SpO<sub>2</sub>, and Braun ThermoScan PRO 4000 Ear Thermometer



Spot Vital Signs® LXi Blood Pressure Accessories



**45NE0-E1**  
Spot Vital Signs LXi with SureBP, Nellcor SpO<sub>2</sub> and Braun ThermoScan PRO 4000 Thermometer

Two-Piece Cuffs and Bladders

- 4500-01** Child Cuff and Bladder, two tubes [size 9]
- 4500-02** Adult Cuff and Bladder, two tubes [size 11]
- 4500-03** Large Adult Cuff and Bladder, two tubes [size 12]
- 4500-04** Thigh Cuff and Bladder, two tubes [size 13]
- 4500-200** Cuff Kit (child, adult and thigh cuffs, bladders, two tubes)

Durable One-Piece Cuffs

- 5082-82-4MQ** Infant Two-Tube Cuffs [size 7], (5/box)
- 5082-83-4MQ** Small Child Two-Tube Cuffs [size 8], (5/box)
- 5082-84-4MQ** Child Two-Tube Cuffs [size 9], (5/box)
- 5082-85-4MQ** Small Adult Two-Tube Cuffs [size 10], (5/box)
- 5082-86-4MQ** Adult Two-Tube Cuffs [size 11], (5/box)
- 5082-87-4MQ** Large Adult Two-Tube Cuffs [size 12], (5/box)
- 5082-88-4MQ** Thigh Two-Tube Cuffs [size 13], (5/box)
- 4500-201** Adult Cuff Kit (child, small adult, large adult two-tube cuffs and bladders)
- 4500-202** Pediatric Cuff Kit (small child, child, small adult two-tube cuffs and bladders)

Disposable One-Piece Cuffs

- 5082-92-4MQ** Infant Two-Tube Cuffs [size 7], (5/box)
- 5082-93-4MQ** Small Child Two-Tube Cuffs [size 8], (5/box)
- 5082-94-4MQ** Child Two-Tube Cuffs [size 9], (5/box)
- 5082-95-4MQ** Small Adult Two-Tube Cuffs [size 10], (5/box)
- 5082-96-4MQ** Adult Two-Tube Cuffs [size 11], (5/box)
- 5082-97-4MQ** Large Adult Two-Tube Cuffs [size 12], (5/box)
- 5082-98-4MQ** Thigh Two-Tube Cuffs [size 13], (5/box)

Miscellaneous Blood Pressure Accessories

- 4500-30** Blood Pressure Hose (5 ft, 1.5 m)
- 5200-08** Calibration "T" Connector

SURETEMP PLUS THERMOMETRY ACCESSORIES



**05031-101**  
Disposable  
Probe Covers

BRAUN THERMOSCAN PRO 4000 THERMOMETRY ACCESSORIES



**53020-0000**  
Rechargeable  
Battery Pack



**05075-800**  
Disposable  
Probe Covers



## Spot Vital Signs® LXi Thermometry Accessories

### SureTemp® Plus Thermometry Accessories

- 02895-000** Oral Temperature Probe and Well Kit (9 ft, 2.7 m)
- 02895-100** Rectal Temperature Probe and Well Kit (9 ft, 2.7 m)
- 02894-0000** Oral Probe Well (blue)
- 02894-1000** Rectal Probe Well (red)
- 05031-101** Disposable Probe Covers (1,000 covers, 25/box)
- 01802-110** 9600 Plus Calibration Tester
- 06138-000** SureTemp Plus Calibration Key

### Braun ThermoScan® PRO 4000 Thermometry Accessories

- 04000-200** Braun ThermoScan PRO 4000 Ear Thermometer (without rechargeable batteries)
- 53020-0000** Braun ThermoScan PRO 4000 Rechargeable Battery Pack
- 4500-53** Braun ThermoScan PRO 4000 Unlocking Key
- 05075-800** Braun ThermoScan PRO 4000 Disposable Probe Covers (800 covers, 20/box)
- 01802-110** 9600 Plus Calibration Tester



Spot Vital Signs® LXi Pulse Oximetry Accessories

Masimo Pulse Oximetry Accessories

<b>LNCS-DCI</b>	Adult Reusable Sensor, 3 ft
<b>LNCS-DCIP</b>	Pediatric Reusable Sensor, 3 ft
<b>LNCS-ADTX</b>	Adult Adhesive Sensors, 18 inches (20/box)
<b>LNCS-PDTX</b>	Pediatric Adhesive Sensors, 18 inches (20/box)
<b>LNCS-INF-L</b>	Infant Adhesive Sensors, 3 ft (20/box)
<b>LNCS Neo-L</b>	Neonatal Adhesive Sensors, weight < 3 or > 40 kg, 3 ft (20/box)
<b>LNCS NeoPt-L</b>	Neonatal Adhesive Sensors, weight < 1 kg, 3 ft (20/box)
<b>LNC-4</b>	LNC Patient Cable, 4 ft (14 Pin Connector)
<b>LNC-10</b>	LNC Patient Cable, 10 ft (14 Pin Connector)
<b>LNOP-ADT</b>	Disposable Adhesive Finger Sensor—Adult (case of 20)
<b>LNOP-PDT</b>	Disposable Adhesive Finger Sensor—Pediatric (case of 20)
<b>LNOP-NEO</b>	Disposable Adhesive Finger Sensor—Neonate (case of 20)
<b>LNOP-NEOPT</b>	Disposable Finger Sensor—Neonate, Sensitive Skin (case of 20)
<b>LNOP-DCI</b>	Reusable Finger Sensor—Adult
<b>LNOP-DCIP</b>	Reusable Finger Sensor—Pediatric
<b>PC-04</b>	4 ft Cable with Sensor Connector
<b>PC-08</b>	8 ft Cable with Sensor Connector
<b>PC-12</b>	12 ft Cable with Sensor Connector

Nellcor® Pulse Oximetry Accessories

<b>DS-100A</b>	Durasensor Adult Oxygen Transducer
<b>DOC-10</b>	Extension Cable (10 ft)
<b>D-YS</b>	Dura-Y Oxygen Transducer (1 sensor, 40 wraps)
<b>D-YSE</b>	Ear Clip (use with Dura-Y sensor)
<b>D-YSPD</b>	PediCheck Pediatric Spot Check (use with Dura-Y sensor)
<b>MAX-AI</b>	OxiMax Adult Sensor (single use, case of 24)
<b>MAX-PI</b>	OxiMax Pediatric Sensor (single use, case of 24)
<b>MAX-II</b>	OxiMax Infant Sensor (single use, case of 24)
<b>MAX-RI</b>	OxiMax Adult Nasal Sensor (single use, case of 24)
<b>OXCILIQ-PI</b>	Pediatric Oxygen Transducer (for use with OC-3 Cable)
<b>OC-3</b>	OxiCliq Sensor Cable
<b>OXI-A/N</b>	Oxiband Adult/Neonatal Transducer (1 sensor, 50 wraps)
<b>OXI-P/I</b>	Oxiband Pediatric/Infant Transducer (1 sensor, 50 wraps)
<b>SRC-MAX</b>	Portable Oximetry Tester



**45NT0-E1**  
Spot Vital Signs LXi  
with SureBP, Nellcor  
SpO<sub>2</sub>, and SureTemp  
Plus Thermometer

**4500-60**  
Mobile Stand with Basket

### Spot Vital Signs® LXi Miscellaneous Accessories

<b>4500-60</b>	Mobile Stand with Basket
<b>4500-62</b>	Wall Mount with Basket
<b>4500-101A</b>	AC Power Transformer
<b>4500-400</b>	Line Cord (USA/Canada/Japan version)
<b>4500-910</b>	Linear Barcode Scanner Kit (Includes: Scanner, Mounting Bracket, Hardware and Instructions)
<b>4500-915</b>	2D Barcode Scanner Kit (Includes: Scanner, Mounting Bracket, Hardware and Instructions)
<b>4500-925</b>	Connectivity Cable, USB for Spot LXi
<b>4500-500</b>	External Printer Kit (Includes: Printer, Cable, Mounting Bracket, Hardware, Paper Roll and Instructions)
<b>4500-505</b>	Rechargeable Replacement Battery for Printer
<b>4500-510</b>	Spot LXi Printer Paper (10 pack)
<b>4500-520</b>	Label Paper and Cleaning Kit (10 rolls/1 cleaning kit)
<b>4500-84</b>	Lead Acid Battery
<b>4500-100</b>	Carrying Case



**45NT0-E1**  
Spot Vital Signs LXi  
with SureBP, Nellcor  
SpO<sub>2</sub>, and SureTemp Plus  
Thermometer

**4500-500**  
External Printer Kit



**45ME0-E1**  
Spot Vital Signs LXi with  
SureBP, Masimo SpO<sub>2</sub> and  
Braun ThermoScan PRO  
4000 Thermometer

**4500-915**  
Barcode Scanner Kit



**45NE0-E1**  
Spot Vital Signs LXi with  
SureBP™, Nellcor SpO<sub>2</sub>  
and Braun ThermoScan  
PRO 4000 Thermometer

**4500-62**  
Wall Mount with Basket



**42MTB-E1**  
Spot Vital Signs with NIBP, Masimo Pulse Oximetry and SureTemp Thermometer

## Spot Vital Signs®

Vital signs in seconds with hospital-grade technology that's easy to use, accurate and affordable.

- > Measures noninvasive blood pressure (NIBP), pulse rate, temperature and SpO<sub>2</sub> in about 30 seconds
- > Optional—SureTemp® thermometry, Masimo® or Nellcor® pulse oximetry
- > IR Connectivity to electronic medical records
- > Lightweight and portable with a variety of mounting options
- > All configurations include NIBP, pulse rate and MAP

- 4200B-E1** Spot Vital Signs with NIBP only (pulse rate and MAP)
- 420TB-E1** Spot Vital Signs with NIBP and SureTemp Thermometer
- 42M0B-E1** Spot Vital Signs with NIBP and Masimo Pulse Oximetry
- 42N0B-E1** Spot Vital Signs with NIBP and Nellcor Pulse Oximetry
- 42MTB-E1** Spot Vital Signs with NIBP, Masimo Pulse Oximetry and SureTemp Thermometer
- 42NTB-E1** Spot Vital Signs with NIBP, Nellcor Pulse Oximetry and SureTemp Thermometer
- 4200B-E1-M** Spot Vital Signs with NIBP only (pulse rate and MAP) and Mobile Stand
- 420TB-E1-M** Spot Vital Signs with NIBP and SureTemp Thermometer with Mobile Stand
- 42N0B-E1-M** Spot Vital Signs with NIBP and Nellcor Pulse Oximetry with Mobile Stand
- 42NTB-E1-M** Spot Vital Signs with NIBP, Nellcor Pulse Oximetry and SureTemp Thermometer with Mobile Stand
- 42MTB-E1-M** Spot Vital Signs with NIBP, Masimo Pulse Oximetry and SureTemp Thermometer with Mobile Stand
- 42M0B-E1-M** Spot Vital Signs with NIBP and Masimo Pulse Oximetry with Mobile Stand

## Spot Vital Signs® Blood Pressure Accessories

### Durable One-Piece Cuffs

- 5082-203-3** Small Child One-Tube Cuff [size 8]
- 5082-204-3** Child One-Tube Cuff [size 9]
- 5082-205-3** Small Adult One-Tube Cuff [size 10]
- 5082-206-3** Adult One-Tube Cuff [size 11]
- 5082-207-3** Large Adult One-Tube Cuff [size 12]
- 5082-208-3** Thigh One-Tube Cuff [size 13]

### Two-Piece Cuffs and Bladder

- 5200-01** Adult Cuff and Bladder [size 11]
- 5200-02** Large Adult Cuff and Bladder [size 12]
- 5200-03** Child-Print Cuff and Bladder [size 9]
- 5200-10** Thigh Cuff and Bladder [size 13]

*(Spot Vital Signs Blood Pressure Accessories continued on next page)*

## Spot Vital Signs® Blood Pressure Accessories (continued)

### Disposable One-Piece Cuffs

**5082-93-3** Small Child One-Tube Cuffs [size 8], (box of 5)

**5082-94-3** Child One-Tube Cuffs [size 9], (box of 5)

**5082-95-3** Small Adult One-Tube Cuffs [size 10], (box of 5)

**5082-96-3** Adult One-Tube Cuffs [size 11], (box of 5)

**5082-97-3** Large Adult One-Tube Cuffs [size 12], (box of 5)

**5082-98-3** Thigh One-Tube Cuffs [size 13], (box of 5)

### Replacement Cuffs & Bladders

**5200-06** Small Child Bladder (latex-free, includes connector)

**5200-04** Adult Bladder (latex-free, includes connector)

**5200-05** Large Adult Bladder (latex-free, includes connector)

**5200-11** Thigh Bladder (latex-free, includes connector)

**5082-18** Child-Print Cuff [size 9]

**5082-01** Adult Cuff [size 11]

**5082-16** Large Adult Cuff [size 12]

**5082-64** Thigh Cuff [size 13]


### Miscellaneous Blood Pressure Accessories

**5200-12** Straight Pressure Hose (8 ft, 2.4 m)


**5200-19** Straight Pressure Hose (5 ft, 1.5 m)

**5200-08** Calibration "T" Connector

MASIMO PULSE OXIMETRY ACCESSORIES



NELLCOR PULSE OXIMETRY ACCESSORIES



Spot Vital Signs® Pulse Oximetry Accessories

Masimo Pulse Oximetry Accessories

- LNCS-DCI** Adult Reusable Sensor, 3 ft
- LNCS-DCIP** Pediatric Reusable Sensor, 3 ft
- LNCS-ADTX** Adult Adhesive Sensors, 18 inches (20/box)
- LNCS-PDTX** Pediatric Adhesive Sensors, 18 inches (20/box)
- LNCS-INF-L** Infant Adhesive Sensors, 3 ft (20/box)
- LNC-4-WA** 4 ft Patient Cable with 9 Pin Connector for LNCS
- LNC-10-WA** 10 ft Patient Cable with 9 Pin Connector for LNCS
- LNOP-DCI** Reusable Finger Sensor—Adult
- LNOP-DCIP** Reusable Finger Sensor—Pediatric
- LNOP-ADT** Adhesive Finger Sensor—Adult (20/case)
- LNOP-PDT** Adhesive Finger Sensor—Pediatric (20/case)
- LNOP-INFL** Adhesive Finger Sensor—Infant (20/case)
- PC-04-WA** 4 ft Cable with DB-9 Connector for LNOP
- PC-08-WA** 8 ft Cable with DB-9 Connector for LNOP

Nellcor Pulse Oximetry Accessories

- DS-100A** Durasensor Adult Oxygen Transducer
- DEC-4** Extension Cable (4 ft)
- DEC-8** Nellcor - 8 ft Differential Extension Cable, Pulse Oximeter Sensor (use only with Nellcor sensors)
- D-Y** Nellcor Dura-Y Oxygen Transducer (1 sensor, 40 wraps)
- D-YSE** Nellcor Ear Clip (use with Dura-Y sensor)
- D-YSPD** PediCheck Pediatric Spot Check (use with Dura-Y sensor)
- MAX-AI** OxiMax Adult Sensor (single use, case of 24)
- MAX-PI** OxiMax Pediatric Sensor (single use, case of 24)
- MAX-II** OxiMax Infant Sensor (single use, case of 24)
- MAX-RI** OxiMax Adult Nasal Sensor (single use, case of 24)
- OXICLIQ-PI** Pediatric Oxygen Transducer (for use with OC-3 Cable)
- OC-3** OxiCliq Sensor Cable
- OXI-A/N** Oxiband Adult/Neonatal Transducer (1 sensor, 50 wraps)
- OXI-P/I** Oxiband Pediatric/Infant Transducer (1 sensor, 50 wraps)
- SRC-MAX** Portable Oximetry Tester

Thermometry Accessories (Spot Vital Signs only)

- 02678-100** Oral Probe (9ft, 2.7 m)
- 02679-100** Rectal Probe (9ft, 2.7 m)
- 05031-101** Disposable Probe Covers (1,000 covers, 25/box)
- 06137-000** Temperature Calibration Key
- 01802-110** 9600 Plus Calibration Tester



**4200-62**  
Wall Mount  
with Basket

### Spot Vital Signs® Mounting Options

- 4200-60** Mobile Stand with Basket
- 4200-62** Wall Mount with Basket
- 4200-64** I.V. Pole Mount with Basket
- 4200-70** Anti-Theft Kit

### Transformers and Line Cords

- 5200-101A** AC Power Transformer
- 76400** Line Cord

### Miscellaneous Accessories

- 4200-84** Lead Acid Battery
- 4200-100** Carrying Case
- 4200-170** Connectivity Accessory Kit
- 4200-350** Pediatric Accessory Kit (Includes: Small Child, Child and Small Adult Cuffs and Pediatric Masimo SpO<sub>2</sub> Sensor)
- 4200-250** Pediatric Accessory Kit (Includes: Small Child, Child and Small Adult Cuffs and Pediatric Nellcor SpO<sub>2</sub> Sensor)
- 4200-00B** One-year extended warranty, Model 4200B
- 4200-0TB** One-year extended warranty, Model 420TB
- 4200-M0B** One-year extended warranty, Model 42MOB
- 4200-MTB** One-year extended warranty, Model 42MTB
- 4200-N0B** One-year extended warranty, Model 42NOB
- 4200-NTB** One-year extended warranty, Model 42NTB