

97200-MS: Classic Smart Set

• Coaxial Ophthalmoscope

- Halogen HPX™ illumination provides 30% brighter and whiter light than traditional ophthalmoscope lamps, for truer tissue color rendition
- Compact size makes access to all controls easy and convenient, which makes learning faster
- Triple-sealed scope keeps dirt and dust away from optics, for years of maintenance-free viewing
- Includes 28 built-in focusing lenses

• MacroView Otoscope/Throat Illuminator

- New optics system provides nearly twice the field of view and 30% higher magnification than traditional diagnostic otoscopes
- Focusing wheel allows you to fine tune the view
- Fiber optics transmit cool, white light to the distal tip of the otoscope in a 360 degree bundle for accurate, shadow-free exams and easier cleaning
- Halogen HPX™ illumination provides 30% brighter and whiter light than regular halogen, for truer tissue color rendition

• Lithium Ion Rechargeable Power Handle

- Weighs half as much as a Ni-Cad handle, delivering ergonomic comfort for any sized hand
- Low-battery indicator light provides ample warning that battery will need to be recharged soon
- Rugged Cycloy thermoplastic resin construction provides durability and lightweight balance
- AC charging module has retractable prongs and is easily detached from the handle