

## **97800-MS: PanOptic Smart Set**

### **. PanOptic™ Ophthalmoscope**

- Provides a 5X larger view of the fundus than ever before achieved with a standard ophthalmoscope in an undilated eye
- Enables a 25° field-of-view vs. the standard 5° field-of-view
- Increases magnification by 26% over a standard scope
- Greater working distance between practitioner and patient improves comfort for all

### **. MacroView™ Otoscope/Throat Illuminator**

- New optics system provides nearly twice the field of view and 30% higher magnification than traditional diagnostic otoscopes
- Focusing wheel allows you to fine tune the view
- Fiber optics transmit cool, white light to the distal tip of the otoscope in a 360 degree bundle for accurate, shadow-free exams and easier cleaning
- Halogen HPX™ illumination provides 30% brighter and whiter light than regular halogen, for truer tissue color rendition

### **. Lithium Ion Rechargeable Power Handle**

- Weighs half as much as a Ni-Cad handle, delivering ergonomic comfort for any sized hand
- Low-battery indicator light provides ample warning that battery will need to be recharged soon
- Rugged Cycloy thermoplastic resin construction provides durability and lightweight balance
- AC charging module has retractable prongs and is easily detached from the handle