

Welch Allyn Green Series™ Exam Light IV



Highly Flexible Light Pipe—
Twice as long as conventional pipes
and does not droop or wander

**Lightweight Focusing
Sleeve—**Provides a variable
homogenous spot

1 LED—3X the typical light
output; pinpoint accuracy;
true tissue color rendition



Touchless On/Off

**Disposable Sheath
Available—**Part #52640



Easy Power
Connector

Whiter, Brighter, and Greener

The GS Exam Light IV features:

- A single LED that provides whiter, brighter light than halogen, so you see truer tissue color rendition during an exam
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- A compact, unobtrusive head, ensuring pinpoint accuracy
- Minimal heat dissipation
- Touchless on/off
- 5-year limited warranty

Easy Assembly



WelchAllyn
Advancing Frontline Care™



Green Series Medical Exam Lights

Welch Allyn Green Series™ Exam Light IV



ORDERING INFORMATION

- 48810 GS Exam Light IV with Table/Wall Mount
- 48950 Mobile Stand for GS Exam Light IV
- 52640 Disposable Sheath for GS Exam Light IV (25/box, 5 boxes/case)

SPECIFICATIONS

# of LEDs	1 White LED
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	20,000 lux (1,858 footcandles)
LED Life	50,000 hours
Working Distance	16 in (40 cm)
Spot Size at Working Distance	4 in (10 cm) to 10 in (25 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 0.30 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 5.30 lbs (2.40 kg) Total Mobile Stand: 12.40 lbs (5.63 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Focusing Sleeve: Length: 5.16 in (13.1 cm) Diameter: 1.58 in (4 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 47.24 in (120.0 cm) Total Mobile Stand: 60.24 in (153.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 49.61 in (126 cm)
Mounting Options	Table/Wall Mount Mobile Stand
Control Box	Touchless On/Off Control
Accessory Options	Disposable Sheath #52640
Warranty	5-year limited

IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- OB/GYN
- OUTPATIENT FACILITIES
- FAMILY PRACTICE
- LABOR AND DELIVERY SETTINGS
- DERMATOLOGY
- URGENT CARE CENTERS
- EMERGENCY ROOMS

MOUNTING OPTIONS AVAILABLE:



Table Mount

Wall Mount

Mobile Stand

To place an order, contact our Customer Service Department at 1.800.535.6663

