

# Welch Allyn Green Series™ 600 Minor Procedure Light



Disposable Sheath  
Available—Part #52630



Easy Power  
Connector

Easy Assembly



3 LEDs



Touchless On/Off  
with Dimming Control

## Whiter, Brighter, and Greener

The GS 600 Minor Procedure Light features:

- Three LEDs that provide a whiter, brighter light than halogen
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- Wider head, designed specifically for demanding viewing conditions
- Uniform spot and superior brightness for minor procedures and exams in any healthcare environment
- Minimal heat dissipation
- Touchless controls with dimming features
- 5-year limited warranty

**WelchAllyn**  
Advancing Frontline Care™



Green Series Medical Exam Lights

# Welch Allyn Green Series™ 600 Minor Procedure Light

## ORDERING INFORMATION

- 44600 GS 600 Minor Procedure Light with Mobile Stand
- 44610 GS 600 Minor Procedure Light with Table/Wall Mount
- 52630 Disposable Sheath for GS 600/GS 900 (50/box, 5 boxes/case)

## SPECIFICATIONS

# of LEDs	3 White LEDs
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	50,000 lux (4,645 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	6 in (15 cm) to 7 in (18 cm)
Electrical	<b>Input:</b> 100-240 V~ 50 Hz-60 Hz, 0.65 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	<b>Weight:</b> Total Table/Wall Mount: 6.45 lbs (2.93 kg) Total Mobile Stand: 17.35 lbs (7.87 kg) <b>Control Box:</b> Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) <b>Luminaire:</b> Height: 13.58 in (34.5 cm) Width: 7.56 in (19.2 cm) Depth: 3.39 in (8.6 cm) <b>Mobile Stand:</b> Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) <b>Total Wall Mount:</b> 33.27 in (84.5 cm) <b>Total Mobile Stand:</b> 58.66 in (149.0 cm) <b>Table/Wall and Mobile Stand to Luminaire Center:</b> 31.1 in (79 cm)
Mounting Options	Table/Wall Mount Mobile Stand
Control Box	Touchless On/Off Control with Dimming
Accessory Options	Disposable Sheath #52630
Warranty	5-year limited



## IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- PRIMARY CARE CLINICS
- OB/GYN
- OUTPATIENT FACILITIES
- DERMATOLOGY
- HOSPITAL—ER AND ICU
- URGENT CARE CENTERS

## MOUNTING OPTIONS AVAILABLE:



Table/Wall Mount



Mobile Stand

To place an order, contact our Customer Service Department at 1.800.535.6663

