

White Cassette Carts

All steel construction with powder coat finish, units have rubber mat inside each compartment and padded bar-style handles and scuff resistant tops.



Tall Cassette Cart

Eight 6" wide compartments with 16 3/4" vertical clearance, 24" L x 18" W x 42 1/2" H. Weight 120 lbs. #80-038



Stationary version of the tall Cassette Cart. (no casters). #80-037



Mini Cassette Cart

Four 6" wide compartments with 16 3/4" vertical clearance, 24" L x 18" W x 36" H. Weight 95 lbs. #80-024

CR-Plate Carts

Carts are available in "Tall", "Stationary" and "Mini". They look exactly like the Cassette carts shown above except for the addition of the 3/8" padding.

Tall CR-Plate Cart #80-138 Stationary CR-Plate Cart #80-137
Mini CR-Plate Cart #80-124



Mammography Cassette Rack

All steel construction with white powder coat finish. Mounts to wall or for countertop use. 9 3/4" L x 9" H x 7 1/4" W. 8 lbs. # MCH-20



Multiple Plate/Cassette Rack

Made of steel with a white powder coat finish, the rack keeps CR-plates or cassettes clean and organized in four compartments. 24" L x 9" W x 39" H, 43 lbs. #80-220



File Hanging System

Available in letter or legal size and your choice of translucent or opaque plexi front.

Description	W	H	D	Trans.	Opaque
Letter Size	12"	9"	2 1/2"	21611C	21611O
Legal Size	14"	12"	3"	21612C	21612O



Wolf X-ray Utility File

Our Utility File is available with an opaque or translucent plexi front and in two sizes. Wall mounted via 2 keyholes, the file will hold your largest cassettes, films or jackets.

File with Translucent Bronze Plexi Front

Description	Width	Height	Depth	Stock #
Single	22 1/4"	22 3/8"	5 1/4"	21609
Double	22 1/4"	43"	5 1/4"	21610

File with Opaque Plexi Front

Single	22 1/4"	22 3/8"	5 1/4"	21609O
Double	22 1/4"	43"	5 1/4"	21610O

Top Loading CR-Cart



Like the "Tall and "Mini" CR Carts our Top Loading carts feature 3/8" high impact padding on surfaces that could damage CR-Plates. Available with 2 compartments (Length) or 4 compartments (width)



Two compartments
8"W x 24"L x 13" deep
#80-046

Four Compartments
6"W x 16"L x 13" deep
#80-045

Overall cart size 24"L x 16"W x 26" H. 90 lbs.

Activity Carts

All steel construction with white powder coat finish with rubber mat inside each compartment. Our Activity Carts have sturdy 4" plate casters: two stationary, two swivel, padded bar-style handle and rubber corner bumpers.



Large Cart

Four compartments, each is 6" x 24" x 13" deep. Overall dimensions 24"L x 24"W x 36"H. 127 lbs. #80-414C



Small Cart

Two compartments:
8"Wx24"Lx13"deep
← #80-464C

Four compartments:
6"Wx16"Lx13"deep
#80-454C →



Our small activity cart is available with the interior configured into two or four compartments. Like all our carts, a rubber mat at the bottom of the cart insures that whatever you put in it is cushioned and secure. Overall dimensions: 24" L x 16" W x 36" H. 85 lbs.

Short Cart



Ideal for use while seated, this cart comes with 2 or 4 interior files.
2 files: 8"W x 24"L x 13" D. #80-474C
4 files: 6"W x 16"L x 13" D. #80-484C
The left hand picture shows the handle up, the right, handle down.



Overall size: 24"L x 16"W x 27" H. 82 lbs.



Wolf Heavy Duty Cassette Cart

Designed for rugged use, the heavy duty cassette cart is constructed of heavy gauge steel, with a tough washable white finish. The cart rolls about easily on 4 sturdy rubber casters that lock in position. This model has 3 levels divided into 5 compartments that hold any size cassette from 14 X 17" on down (maximum 125 cassettes). 50" H x 17" D x 20" W. 100 lbs. # 25515



Wolf Cassette Carts

Constructed of anodized aluminum with gray panels, the Cassette Carts carry up to 64 cassettes. Dividers keep cassettes separated, and 1 inch tread strip prevents slippage. Choose from bi-level, or single level (not shown) .

Bi-Level (64 cassettes)
50" H x 32¹/₂" L x 17" W. 64 lbs # 25514

Single-Level (32 cassettes)
42" H x 32" L x 17" W
Approx. shipping weight. 50 lbs.
25513



Film Filing Carts

Handy for hospital rounds. It's constructed of heavy satin anodized aluminum extrusions and laminated gray composition board. Four casters provide total mobility, prevent tipping. Heavy duty wire dividers facilitate filing. Weight 20 lbs.

Description	Dimensions			Number
	L	W	H	
Stand-Up (10 Dividers)	32 ¹ / ₂ "	21 ¹ / ₂ "	40 ¹ / ₂ "	25505
Desk Height (10 Dividers)	29"	21 ¹ / ₂ "	28 ¹ / ₄ "	
Desk Height	24 ³ / ₄ "	17 ¹ / ₄ "	28 ¹ / ₄ "	25503



Wolf MRI Cart and Stool

Constructed entirely of PVC tubing and plastic, no metal parts to interfere with the magnetic field. The 3 level cart holds all your MRI accessories, rolls easily on casters. 35" H x 29" D x 16" W. 40 lbs. # 22105

The stool features rubber "grip" tips on the legs and top surface to prevent slipping. 8" High x 12" D x 15¹/₂" W. 10 lbs. # 22106



Stools

Our Foot Stool is made of steel with a white finish and a rubber tread top. 17¹/₂" W x 13" D x 8" H. 8 lbs. # 22102

The **Foot Stool with Handle** is the same as above, with a high handle. 17¹/₂" W x 13" D x 39" H. 15 lbs. # 22103

Kick Stool rolls wherever you want it. Rollers retract when you stand on it. 16" D x 16" H. # 22104



File Stak

File-Stak takes only 19¹/₂" x 32" of valuable floor space, stacks as high as you want. Made of sturdy 20-gauge steel. Individual units are secured to each other by bolts and come with 3 adjustable dividers. Paneled doors slide easily at the touch of a finger for ready access. Overall finish is in attractive, easy to clean white baked on enamel.

Cabinet 32¹/₄" long, 19¹/₂" wide
16¹/₂" high 50 lbs. # 25125
Base 34¹/₄" long, 19¹/₂" wide,
6" high 12¹/₂ lbs. # 25126
Dividers # 25128