



### Standard Tilt and Rotate Cassette Holder

- Straight Vertical
- Incline Backwards
- Incline Forward
- Rotate or Oblique Tilt
- Face Up
- Face Down

The industry standard, this mobile unit allows all the positions bullet pointed above. Vertical range from bottom of cassette/plate: floor-to-54". The arm slides vertically on a 60" chrome steel post. Available in neutral almond baked enamel. Base is 18" square. **#80-133**



### Extension Arm Tilt and Rotate Cassette Holder

Essentially the same unit as the "Standard", the Extension Arm features an "extending arm." Cannot be used for DR. The arm slides vertically on a 60" chrome steel post. Available in neutral almond baked enamel. Base is 18" square. **#80-237**



### DR-Plate Holder

Our DR-Plate holder does not tilt or rotate and is for lateral positions only (non articulating). It holds DR-plates up to 14"x17". The arm slides vertically on a 60" chrome steel post. Available in neutral almond baked enamel. Base is 18" square. **#80-140**



### Tabletop Lateral Cassette CR-Plate Holder

This 24" variable height cassette/plate holder accommodates up to a 14"x17" cassette. The two 6" wide channels are fully adjustable and are lined with rubber to firmly grip any cassette/plate. 7 lbs. **#80-035**


### Disposable Cassette & CR Plate Covers



Our disposable Cassette/CR Plate covers prevent cross-contamination and protect physicians and patients. They dispense easily from the box and have a sandwich bag style folding flap, can also be used to cover positioning sponges.

Cassette Size	Cat. # 100/box	Cat. # 500/box
8 X 10"	80-008	80-010
10 X 12"	80-012	80-112
14 X 14"	80-014	80-017

### Ortho Disc



The Ortho Disc allows you to establish a reference point of known size to calibrate or scale objects on an x-ray. Perfect for sizing hip and knee joints and judging angles in x-rays of the hip and knee joints. The kit comes with the large height adjustable disc which has an imbedded radiopaque disc that indicates a known location, a small disc and a baton. **Stock Number 55400**