

Film Bins



Floor Standing

This light-proof cabinet is designed for film storage of sizes 5" x 7" to 14" x 17". It's hinged at the bottom, providing easy access without bending over. The non-moveable dividers provide an easy filing system. This film bin's height is ideal for placement under a counter top. An optional switch is available for automatic room light cut-off when the door is opened. Holds 400 pieces of film. (4 boxes)

Specifications: 32½" H x 21½" W x 12" D
Approx. shipping weight: 60 lbs. #25118
With automatic cut-off switch #25119



Table Top

These bins can also be mounted on the side of a desk via 2 keyholes. Wolf manufactures 2 sizes: the "standard" (shown) holds up to 150 14 x 17" films; Specifications: 19" L x 15½" W x 11½" H #25120.

The "mammo" (not shown) holds up to 150 mammography films. Specifications: 13¼" L x 12½" W x

10½" H. #25121. On both feature the door closes automatically after films are removed.

Econo-Top

One for films up to 10" x 12" it holds up to 150 10" x 12" films. 13" Wide, 18" deep, 3" high. #25130.



One for films up to 14 x 17", it holds up to 150 14" x 17" films.

Dimensions 17" wide, 24"

deep, 3" high. #25131.

Safelights

Swivel



Perfect for direct spot lighting this wall mounted unit swivels to any position. Easy to install via 2 keyholes, the Swivel Safelight includes a GBX type filter, light bulb, and a 6 foot cord with a three prong plug. Specification: 7½" H x 6" W, Approx. shipping weight 5 lbs. #19198.



Angular



The Angular is ideal where a direct light is wanted. Easy to install via 2 keyholes, the Angular includes a GBX type filter, light bulb, and a 6 foot cord with a three prong plug. Specifications: 10" H x 7½" W x 5½" D. Approximate shipping weight 5 lbs. #19103.

The Swivel and Angular Safelights shown above are also available with an automatic light Sensor which automatically turns the safelight on when the room light is shut off.

Note: the Safelight can still be shut off with the on/off switch. Swivel #19298, Angular #19203.

Combination



This unique two-compartment unit features a safelight at the bottom and an incandescent room light at the top. The safelight is filtered by our 8' x 10' GBX type filter and comes with a 6 foot cord and a 3 prong plug. Specifications: 15½" Long x 10" Wide x 2½" Deep. Approx. shipping weight: 12 lbs. #19106.

Circular



Screws into any existing outlet. Made of aluminum, filtration is from a GBX type filter. Specifications 5¼" L x 4¼" D. Approx. shipping weight.5 lbs. #19104.

Fluorescent



Available with a red or green filter 24" L x 3" W x 6" H. 8 lbs.
Red 19109 R, Filter Only 19109 R-B,
Green 19109G, Filter Only 19109 G-B.

Wolf Film Duplicators



Wolf manufactures three sizes of x-ray duplicators, with screen sizes of 14" x 17", 10" x 12" and 6x12" respectively. Constructed of steel finished in durable baked enamel, our duplicators can be used with any available duplicating film.

To operate the duplicators, simply place the film to be duplicated on the plexiglass screen with a sheet of duplicating film on top, close the cover and set the exposure time and push the button. The red light will signal you when exposure is complete. The 14" x 17" surface of the **21524** accommodates all x-rays up to 14" x 17". The 10" x 12" surface of the **21522** is perfect for duplicating mammographs and podiatric radiographs, and the 6" x 12" **21520** duplicates any size up to 6" x 12".

Timers



Our large faced analog timer features a giant 8" diameter dial that gives you minute and second timing from 1 second to 60 minutes. It automatically shuts-off a duplicator or other electrical appliances, and has a clear buzzer signal. # 18900



Wolf's digital timer is a pocket sized 2 5/8" wide x 2 1/8" high x 1/2" thick. It attaches to anything, has a spring loaded clip that converts into a stand, a magnet for metal surfaces, and a hanging hole. Times from 100 minutes down to 1 second and emits a tone when timing period is completed. # 18906



Darkroom Flashlights

A handy item for any darkroom, GBX type filtered flashlights spot light items in a darkroom or for service technicians working in a darkroom environment. 6 inches long.

Rayovac # 21616 Mini Maglite # 21615

ID Printers



The **Auto Printer** has an adjustable exposure timer. Set the timer, insert ID card and press the exposure plate. # 17105. The **Manual ID Printer** insert the film and ID card, press exposure plate. # 17104

ID Printer Cards

1,000 per box. 3 x 5" Cards
Standard Printed with blanks for patient name and address. # 15100. **Personalized Standard** Same as above with customer imprint (3 line max.) # 15102. **Double Printed Standard** ID blocks printed top and bottom. Standard # 15150 Custom imprinted #15151

Thermometers



Dial Calibrated from 0° - 140°F, made of stainless steel, with a cup for floating in the tank. # 18903
Glass floats in the developing tank. It is calibrated from 0° - 100°F. # 18904
Stainless Steel has a hook for hanging on the tank, calibrated for accuracy between 20° and 150° fahrenheit. # 18905
Digital Runs on a standard watch battery (included) and comes with a lifetime warranty. Registers temperatures from - 40° to 300°F. Place it in the solution and press the button. The temperature is digitally displayed on a LCD. # 19107

"X-Ray In Use" Sign



Wolf's "X-Ray In Use" sign is also available with "MRI In Use", "Darkroom In Use" and "CT In Use" messages. Made of steel and lit by two 15 watt tubular bulbs, the 12" L x 7" H x 2 1/2" D box has your message on both sides.
 # 19110-XR # 19110-DR
 # 19110-MR # 19110-CT



Film Developing Hangers

Drop type construction prevents spilling and streaked negatives. 1-piece stainless steel construction with an interlock clip for perfect pin and hole alignment.

Industrial X-ray Developing Hangers, Two Film

4 1/2" X 10"	16317
4 1/2" X 17"	16319

Stainless Steel Medical X-ray Hangers

8" X 10"	16302
10" X 12"	16303
11" X 14"	16304
14" X 17" (non-drop)	16305
13 X 18 cm	16310

Developing Hanger Drying Rack



This wall mounted rack holds up to 12 developing hangers with wet films (any size hangers), allowing them to air dry. The Wolf drying rack is made of durable stainless steel and is 20" long and comes with 6 screws for mounting. # 16101