

TotalLook®

Easing the transition from Film to Digital



- **Uncompromised image quality with reduced file size**

- > Flawlessly converts film to digital image
 - Industry leading 42.3 micron spot size
 - Fine structure and spatial resolution
- > Intelligent Image Compression
 - Reduces the original file size by nearly 75%
 - Saves valuable storage/archiving space
 - Improves accessibility and distribution of images

- **Speeds workflow**

- > Easy comparative review of priors on single digital workstation
 - Accepts all views
 - Automatically loads priors into hanging protocols
 - Auto labels all views
- > Fastest scan time
 - 23 seconds for 18 x 24 cm films
- > Provides simple “Load & Go” functionality
 - Unaided batch processing of up to 60 films
- > Up to 12 views
- > Reliable with no daily maintenance



- **Seamless integration with PACS, Mammography Review Workstations, and Radiology Information Systems (RIS)**

- > Receive DICOM patient work lists
- > Send DICOM images to multiple destinations simultaneously
- > Filter patient work list to see only the cases you need