

Urine Catheter Kits

ALL YOU NEED TO COLLECT STERILE URINE SAMPLES SAFELY.

Smooth tapered catheter tip and eyes for low-trauma insertions

Sheath protects catheter from bending during shipping

Unique vented lid promotes urine flow without the need to loosen lid

Closed system reduces the risk of specimen contamination and false readings

Printed graduations on tube makes it easier for technician to see levels of urine

15cc centrifuge ready tubes

Catalog No. 410 **Female Catheter Kit**

- rigid catheter, 8 fr, 6" long
- vinyl exam gloves
- 3 iodophor PVP swabsticks
- 3gm lubricant packet
- towel
- patient ID label
- 1 blue plastic pouch
- packed 100/case (4 inner cartons of 25 ea.)

Catalog No. 411 **Adult Catheter Kit (Female or Male Use)**

- soft catheter, 8 fr, 9" long
- vinyl exam gloves
- 3 iodophor PVP swabsticks
- 3gm lubricant packet
- towel
- patient ID label
- 1 blue plastic pouch
- packed 50/case

Catalog No. 412 **Pediatric Catheter Kit**

- soft catheter, 8 fr, 9" long
- vinyl exam gloves
- 3 iodophor PVP swabsticks
- 3gm lubricant packet
- towel
- patient ID label
- instruction sheet
- 1 blue plastic pouch
- packed 50/case

Catalog No. 413 **Infant Catheter Kit**

- soft catheter, 5 fr, 9" long
- vinyl exam gloves
- 3 iodophor PVP swabsticks
- 3gm lubricant packet
- towel
- patient ID label
- instruction sheet
- 1 blue plastic pouch
- packed 50/case

Kits are packaged in a pouch organizer – procedure-ready.



busse
Hospital Disposables

LATEX
FREE



P.O. Box 11067, Hauppauge, NY 11788-0920
800.645.6526 Fax: 631.435.4721
www.busseinc.com mail@busseinc.com

© 2004 Busse Hospital Disposables E.O.E. - M/F/D/V 410