


ErgoFit® Hammock sling

The ErgoFit® Hammock sling is designed to be used with BHM Medical ceiling and floor lifts. The sling provides full head and neck support and double padding in the thigh region for comfort. The quick drying mesh fabric, makes it ideal for bathing as well as general transfers.

A six straps option is available.



| # Product | Description | Dimensions | |
|-----------|---|------------|--------|
| | | in | cm |
| THA-S | Small Hammock Sling | A | 42 106 |
| | Capacity | B | 30 76 |
| | 45-100 lb / 20-45 kg / 3.15-7.14 st | C | 26 66 |
| | | D | 9 23 |
| THA-M | Medium Hammock Sling | A | 49 124 |
| | Capacity | B | 37 94 |
| | 100-200 lb / 45-90 kg / 7.14-14.29 st | C | 34 86 |
| | | D | 11 28 |
| THA-L | Large Hammock Sling | A | 57 145 |
| | Capacity | B | 42 106 |
| | 200-400 lb / 90-180 kg / 14.29-28.57 st | C | 38 96 |
| | | D | 11 28 |



General Characteristics

- Guarantee against manufacturing defects: 1 year
- Maximum loading capacity: 600 lb / 272 kg / / 42.83 st
- Positioning handles at the back to help with transfers, and adequately position the patient above a chair, to simply and quickly transfers from a sitting to an inclined position and to help pivoting the user.
- Machine washable; easy to care for and clean.
- Double padding in the thigh section with ultra flexible nylon/polyester, which prevents skin irritations. Easily installed on a person already sitting or lying down.
- Strong nylon & polyester straps allowing the person to be in a sitting, semi-inclined or inclined position.
- User positioning system based on the color of the straps.
- Straps especially designed for the legs, prevents falls and allows the adjustment and positioning of the legs.
- Additional head support for better comfort.

Technical Specifications

Conform to ISO 10535 & CE marking

*Conditionally to a normal use of an appropriate BHM Medical patient lift and maintenance in accordance with the labels.



| INSTALLATION OPTIONS | | | Body position according to loop strap color combinations | | |
|---|--|----------------------|--|-----------|------|
| Bridge, leg closed with straps interlaced | Bridge, legs open with straps not interlaced | Cradled, legs closed | | Shoulders | Legs |
| | | | | | |
| | | | | Black | Blue |
| | | | | Gray | Blue |
| | | | | Blue | Blue |
| | | | | Blue | Gray |

