



BACKBOARD PLASTIC BLUE 16IN

Plastic backboard 72 inches long by 16 inches wide, load limit 350 lbs. Can be used for water rescue, colour blue.

Product Specifications

Dangerous Goods:	N/A
Depth (in.):	1 3/4
Length (in.):	72
Load Limit (lbs.):	350 lbs. (159 kg)
Model No.:	16
UNSPSC:	42171613
Weight (net):	11 lbs. (5 kg)
Width (in.):	16



BACKBOARD PLASTIC ORANGE 18IN

Plastic backboard 72 inches long by 18 inches wide, load limit 350 lbs. Can be used for water rescue, orange

Product Specifications

Dangerous Goods:	N/A
Depth (in.):	1 3/4
Length (in.):	72
Load Limit (lbs.):	350 lbs. (159 kg)
Model No.:	18
UNSPSC:	42171613
Weight (net):	12 lbs. (6 kg)
Width (in.):	18