



## SCRUB SINKS

**MATERIALS** - Stainless Steel, Type 300 Series, nickel-bearing.

**GAUGES** - Bowls – Windsor Models – 14 gauge. 8000 Series-14 gauge. Windsor-Brackets and hangers- 11 gauge. Aprons- 16 gauge. Base-8000 Series-16 gauge. Mounting brackets- 8000 Series- 12 gauge, “Z” bracket- 14 gauge.

### CONSTRUCTION

**Windsor Models** - Sinks are built with seamless construction, plus 3/4" (1.905) horizontal and vertical corners covered with spherical intersections. A continuous rolled rim of 1-1/2" (3.81) diameter is provided at the front and sides, with sides welded to the backsplash, which extends 8" (20.32) above the rim and returns, 2-1/2" (6.35) at the top and ends. Bowls are 8" (20.32) deep and pitched to a 1-1/2" (3.81) stainless steel waste outlet and tailpiece in the center. Flag brackets at the sides, under the rolled rims, and a “Z” shaped hanger, under the top rear of the backsplash, are furnished for installation.

Sinks provide scrub areas approximately 22" (55.88) wide, each with its own gooseneck outlet.

Options include foot operated mixing valves, wrist blade handled faucets or infrared sensor operated, tempered waterspouts complete with thermostatically controlled mixing valves.

**8000 Series** - Sinks are built with seamless construction with 1" (2.54) horizontal and vertical corners covered with spherical intersections. Bowls are 9-1/4"(23.495) deep at rim and 18-1/2"(46.99) from spout to drain allowing for a deep scrub area. The spout is positioned toward the rear of the sink to divert water splash away from the operator. Undersides of the bowl are sound deadened and enclosed by double walled end panels. Each sink is equipped with a stainless steel base that will house either an infrared sensor or knee panels for manual operation. Base extends 21"(53.34) from the basin and hangs 14-1/2"(36.83) from finished floor. Flag brackets under the sink basin and “Z” shaped hangers under the top rear of the backsplash are furnished for installation. Sinks provide scrub areas approximately 30" (76.20) wide x 21-1/2" (54.61) deep, each with its own swivel gooseneck outlet, sterilizable spray rose aerator, adjustable thermostatic mixing valve, check stop strainer, static sensing probe, flow control valve, 1-1/2" (3.81) stainless steel waste outlet and tailpiece.

**Welding** - All welding is accomplished by heli-arc method. All exposed welds ground smooth.

**Finish** - All exposed stainless steel surfaces are free from pit marks, weld seams and scale, and polished to a No. 4 satin finish.

**BLICKMAN (WINDSOR) SCRUB SINKS** include an 8" deep bowl and 22" scrub area with gooseneck faucet mounted to the back sink panel. All welded, seamless construction and covered corners means easy cleaning; no crevices for dirt to collect. Choose knee-action or infrared sensor controls; knee-action will be provided unless specified. Available in single, double, triple, or four-place sink units. Optional foot or elbow operated controls are available.



Model No.	Controls	Part No.	Width	Depth	No. places
7878SS	Knee-action	1317878000	25" (63.50)	18.5" (46.99)	1
7878SS-IF	Infrared	1317878001	25" (63.50)	18.5" (46.99)	1
7879SS	Knee-action	1317879000	47" (119.38)	18.5" (46.99)	2
7879SS-IF	Infrared	1317879001	47" (119.38)	18.5" (46.99)	2
7880SS	Knee-action	1317880000	69" (175.26)	18.5" (46.99)	3
7880SS-IF	Infrared	1317880001	69" (175.26)	18.5" (46.99)	3
7881SS	Knee-action	1317881000	91" (231.14)	18.5" (46.99)	4
7881SS-IF	Infrared	1317881001	91" (231.14)	18.5" (46.99)	4

**NOTE:** Metric dimensions are shown in centimeters.

Specifications may be improved and dimensions changed without prior notice.



**THE BLICKMAN 8000 SERIES SURGICAL SCRUB SINKS** offer all stainless steel construction with a clean design and low maintenance. OR staff can maintain standards of care for the surgical hand scrub by selecting water activation through an electronic infrared system with proximity sensor or knee-activated water. A swivel gooseneck spout with sterilizable spray rose aerator allows the user to control water flow while a deep sloping basin minimizes splash. Security features include automatic shutoff when temperature exceeds 115°F (46°C). An easy to remove panel provides access to plumbing and electrical components. A removable Plexiglas divider separates stations with multiple bays.

**OPERATIONS:**

**Infrared Units:** Water flow is activated through a continuous beam emitted from a sensor in the front of the sink. When the operator is in the range of the sink (approx. 18"), there is a 2 second delay before water flow is activated. It automatically shuts off after the operator moves out of the range of the sensor.

**Knee-Activated Units:** Water flow is activated through a one-touch on/off water valve. A removable, watertight, sloped control panel is located behind the faucet for accommodating mixing valve and soap spout. A thermostatically controlled mixing valve allows selection of water temperature up to 115°F.

**INSTALLATION:** Wall-mounted units are supported by a "Z" bracket and two main brackets supplied with the sink. It is recommended that a chair carrier be used for support. Connections required are 1/2" N.P.T. cold water, 1/2" N.P.T. hot water and 1-1/2" waste drain per station.



Model No.	Mounting	Controls	Part No.	Width	Depth	Height	No. places
8881SS	wall mount	knee-activated	1388812000	32" (81.28)	27-1/2" (69.85)	37-1/8" (94.3)	1
8881SS-IF	wall mount	infrared sensor	1388812001	32" (81.28)	27-1/2" (69.85)	37-1/8" (94.3)	1
8881SS-FL	floor mount	knee-activated	1388811000	32" (81.28)	27-1/2" (69.85)	51-1/2" (130.81)	1
8881SS-FL-IF	floor mount	infrared sensor	1388811001	32" (81.28)	27-1/2" (69.85)	51-1/2" (130.81)	1
8882SS	wall mount	knee-activated	1388822000	64" (162.56)	27-1/2" (69.85)	37-1/8" (94.3)	2
8882SS-IF	wall mount	infrared sensor	1388822001	64" (162.56)	27-1/2" (69.85)	37-1/8" (94.3)	2
8882SS-FL	floor mount	knee-activated	1388821000	64" (162.56)	27-1/2" (69.85)	51-1/2" (130.81)	2
8882SS-FL-IF	floor mount	infrared sensor	1388821001	64" (162.56)	27-1/2" (69.85)	51-1/2" (130.81)	2
8883SS	wall mount	knee-activated	1388832000	96" (243.84)	27-1/2" (69.85)	37-1/8" (94.3)	3
8883SS-IF	wall mount	infrared sensor	1388832001	96" (243.84)	27-1/2" (69.85)	37-1/8" (94.3)	3
8883SS-FL	floor mount	knee-activated	1388831000	96" (243.84)	27-1/2" (69.85)	51-1/2" (130.81)	3
8883SS-FL-IF	floor mount	infrared sensor	1388831001	96" (243.84)	27-1/2" (69.85)	51-1/2" (130.81)	3

**Options:** Knee-activated soap, scrub timer with display, top storage shelf, pedestal base, in-wall chair carrier.

**NOTE:** Metric dimensions are shown in centimeters.  
Specifications may be improved and dimensions changed without prior notice.



39 Robinson Road, Lodi, NJ 07644 • 800-247-5070 • Fax: 201-909-0832  
e-mail: [info@blickman.com](mailto:info@blickman.com) • [www.blickman.com](http://www.blickman.com)