

# OR Head Covering



## Spunbonded Bouffant Kool-Caps™

Made of lightweight but durable spunbonded fabric that permits increased airflow, keeping the wearer cool and comfortable during long surgical procedures.

### 21" Large Bouffant

**CATALOG NO. 327** White

**CATALOG NO. 328** Blue 

**CATALOG NO. 329** Green

*Packed 500/carton (5 dispensers of 100 each)*


### 24" Extra Large Bouffant for fuller hair styles

**CATALOG NO. 330** White

**CATALOG NO. 331** Blue 

**CATALOG NO. 332** Green

*Packed 500/carton (5 dispensers of 100 each)*



## Surgeon's Cap

Soft and absorbent Kaycel™ material, with wide panel sides, ventilated crown, and adjustable ties.

**CATALOG NO. 335** 

*500 per/carton (5 dispensers of 100 each)*



## Surgeon's Hood and Beard Cover

Designed to cover beards and longer hair, with elastic bottom for easy donning. Rayon.

**CATALOG NO. 336** 

*500/carton (5 dispensers of 100 each)*

## Rayon Bouffant Caps

Made of lightweight, non-woven rayon material that is comfortable and cool to wear and meets FTC requirements.

### 21" Large Bouffant

**CATALOG NO. 268** Blue 

**CATALOG NO. 270** White

*Packed 500/carton (5 dispensers of 100 each)*

# Shoe Covers



## Sur-Step™ Anti-Skid Style

Designed with exclusive anti-skid strips applied to the soles that cut down on accidental slips and falls.

- Made of durable spunbonded polypropylene
- Fit all sizes – small to large
- Elastic top ensures snug fit
- Liquid repellent and abrasion resistant
- Cool and comfortable to wear

Bulk

**CATALOG NO. 346** 

*Bulk packed, 150 pairs/carton*

Dispenser Packed

**CATALOG NO. 348** 

*Packed 150 pairs/carton (3 dispensers of 50)*

## Extra Large Style fits sizes 13 and above

Extra large, Bulk

**CATALOG NO. 347** Sur-Step™

*Bulk packed, 100 pairs/carton* 

Extra large, Bulk

**CATALOG NO. 349** Standard

*Bulk packed, 100 pairs/carton* 



## Standard Style

Features quality and durability at an economical price.

Bulk

**CATALOG NO. 338** 

*Bulk packed, 150 pairs/carton*

Dispenser Packed

**CATALOG NO. 340** 

*Packed 150 pairs/carton (3 dispensers of 50)*

# High Performance Isolation and Procedural Masks




## High Filtration Fluid Resistant Blue Mask with Face Shield

Exceeds 99% filtration efficiency at .1 microns. Combines the performance of the High Filtration Fluid Resistant Mask with an optically clear plastic wraparound face shield, provides peripheral protection for the face and eyes against body fluids. Lightly tinted to help eliminate glare. Anti-fog and anti-static treatment. Earloop design.

**CATALOG NO. 373** 

*Packed 100/carton (4 dispensers of 25 each)*



## High Performance Yellow Isolation Mask

Exceeds 99% filtration efficiency at .1 microns. Ultra soft earloop design. Enclosed, padded nosepiece.

Each layer is ultrasonically bonded to eliminate the sewing pinholes that permit bacteria to migrate from one side to the other.

**CATALOG NO. 370** 

*Packed 500/carton (10 dispensers of 50 each)*



## High Performance Blue Procedural Mask

Exceeds 99% filtration efficiency at .1 microns.

Ultra soft earloop design – comfortable to wear, easy to use! Enclosed, padded nosepiece for extra softness, extra comfort!

**CATALOG NO. 371** 

*Packed 500/carton (10 dispensers of 50 each)*



## High Filtration Fluid Resistant Blue Mask

Fluid Resistant Mask affords the highest level of fluid and respiratory protection.

Exceeds 99% filtration efficiency at .1 microns.

Earloop design. Enclosed, padded nosepiece.

**CATALOG NO. 372** 

*Packed 400/carton (10 dispensers of 40 each)*



## Fluid Resistant Pink Mask

99% filtration efficiency at .1 microns. Enclosed nosepiece. Earloop design. Spunbonded polypropylene.

**CATALOG NO. 374** 

*Packed 500/carton (10 dispensers of 50 each)*



## Fluid Resistant Blue Mask

Same as catalog no. 374 except blue. Spunbonded polypropylene.

**CATALOG NO. 375** 

*Packed 500/carton (10 dispensers of 50 each)*

*Test results available upon request.*

# Surgical Masks and Staff Protection Kits



## Surgical Mask, Blue

Has a soft, non-irritating inner facing and an adjustable aluminum nose piece, to assure maximum comfort during extended use. Composed of highly efficient filtering materials, the Busse Surgical Mask has a bacterial filtration efficiency (BFE) that exceeds the AORN standard of 95%.

**CATALOG NO. 383** 

*Packed 300/carton (6 dispensers of 50 each)*



## High Filtration Multi-Ply Mask, White

When an economical mask with a surgical filter is required. Earloop design for comfort and ease of donning.

**CATALOG NO. 381** 

*Packed 1000/carton (10 dispensers of 100 each)*

## Multi-Ply Mask, White

Our most economical, all purpose mask with earloop design for comfort and ease of donning.

**CATALOG NO. 380** 

*Packed 1000/carton (10 dispensers of 100 each)*



## Deluxe Staff Protection Kit

Contains:

- White, fluid repellent polyethylene coated, spunbonded polypropylene gown, with **full back**, elastic cuffs, and ties at neck and waist
- Fluid resistant earloop mask with clear plastic wraparound face shield
- Pair of examination gloves
- Pair of Sur-Step™ non-skid shoe covers
- 24" blue spunbonded cap
- Red waste bag with tie

**CATALOG NO. 198** 

*Packed 20/carton*

# Protective Isolation and Procedural Gowns





## Staff Protection Gowns

Busse Staff Protection Gowns are made of polyethylene coated, spunbonded polypropylene that is fluid repellent to liquids. All gowns have long sleeves, with either elastic or knit cuffs that fit snugly under gloves. They are available in either full back or open back styles.


Gowns are unisex and unisize to fit both men and women, with neck and waist ties for individual size adjustment.

Fluid resistant poly-coated gown, with **open back and elastic cuffs.**

**CATALOG NO. 199** White   
**CATALOG NO. 201** Yellow   
**CATALOG NO. 210** Blue

*Packed 50/carton*

Fluid resistant poly-coated gown, with **full back and elastic cuffs.**


**CATALOG NO. 200** White   
**CATALOG NO. 233** Yellow

*Packed 50/carton*

Fluid resistant poly-coated gown, with **open back and knit cuffs.**

**CATALOG NO. 230** Blue   
*Packed 50/carton*

Fluid resistant, poly-coated gown, with **full back and knit cuffs.**

**CATALOG NO. 203** Yellow   
**CATALOG NO. 205** White

*Packed 50/carton*



## Isolation Gowns

Busse Isolation Gowns (yellow) are made of fluid resistant, spunbonded polypropylene. Gowns have elastic cuffs and attached neck and waist ties to assure snug fit. Full back. Non-sterile.

Yellow, elastic cuffs.

**CATALOG NO. 202**   
*Packed 50/carton*

Full back, yellow, fluid resistant. Same as Cat. No. 202, but *sterile*. Excellent for use in burn areas or where a patient needs an extra measure of protection.

**CATALOG NO. 204**   
*Packed 20/carton*

# Protective Coveralls and Lab Coats





## Coveralls

Busse coveralls are made of white spunbonded polypropylene. Select catalog numbers available with polyethylene coating to increase fluid repellency.



Generously cut and styled for maximum comfort and mobility.

Nylon zipper front. Available with open or elastic cuffs at both wrist and ankle.


Fluid resistant, open cuffs.

**CATALOG NO. 212** Large   
**CATALOG NO. 213** X Large 

Fluid resistant, elastic cuffs.

**CATALOG NO. 214** Large   
**CATALOG NO. 215** X Large 

Fluid repellent poly-coated, elastic cuffs.

**CATALOG NO. 219** X Large 

*Packed 25/carton*

## Lab Coats

Busse disposable lab coats are made of fluid resistant spunbonded polypropylene and feature ribbed knit cuffs and collar. Lab coats have snap front closure and pockets. Available in blue or white.

Lab coat, Small, Medium.

**CATALOG NO. 224** Blue   
*Packed 30/carton*

Lab coat, Small, Medium.

**CATALOG NO. 226** White   
*Packed 30/carton*

Lab coat, Large, X Large.

**CATALOG NO. 227** Blue   
*Packed 30/carton*

Lab coat, Large, X Large.

**CATALOG NO. 228** White   
*Packed 30/carton*