

## Karapita Ivy Pole

### *ICU Model 0103 - Supplied with 1 bottle holder*

A heavy duty, precision engineered, one hand control unit for:

- ICU
- O.R.
- Coronary Care
- PARR
- Emergency
- Labour & Delivery
- Neonatology
- Pediatrics
- Urology
- Chemotherapy
- Day Care
- Burn Unit
- Dialysis



### Designed for use with IV track

External parts are constructed of type 302 stainless steel, satin finished to reduce glare. Tube wall thickness is .049. Bottle hooks are .250 dia. Rod. Lever is .3125 dia. Rod. Handle is constructed of 6063-T5 Aluminum anodized. Internal component parts are corrosion resistant stainless steel, aluminum and brass.