

PARALLEL BARS

FEATURES

- One piece stainless steel handrails with end bumpers
- Uprights are ultra-heavy, square, telescoping steel tubing
- All models, except wall mounted folding unit, feature width adjustment
- Platform models feature 1 1/2" thick, hardwood, plywood bases with ultra-durable, UV cured, satin finish
- Safety treads at both ends of platform models make it safer for patient to turn and exit wheelchair
- All models are for Physical Therapy use only

	Length	Width Range	Height Range
3-2007	7'	18"-28"	26"-39"
3-2010	10'	18"-28"	26"-39"
3-2012	12'	18"-28"	26"-39"

Platform Mounted Parallel Bars

- Stainless steel adjustable handrails
- 400 lbs. load capacity

1 1/2" Stainless steel handrails



Height is determined by recessed number strip

End bumpers


Chrome-plated indexing tube

Powder-coated, all-steel, base tube upright



Width adjustable handrails

Second knob removes side play



Spring-loaded plungers make adjustments easy

Tapered hardwood ends

PARALLEL BARS



	Length	Width	Height Range
3-3307	7'	28"	28"-42"

Wall Mounted Folding Parallel Bars

- Stainless steel handrails
- Requires only 14' of floor space when stored
- Folds up against the wall when not in use
- 400 lbs. load capacity



	Length	Width Range	Height Range
3-2106	10'	20"-34"	26"-39"

Bariatric Parallel Bars

- Height and width adjustment
- 3rd upright with through pin lock for extra support
- Extra width for bariatric patients
- 48" platform width
- 500 lbs. load capacity



	Length	Width Range	Height Range
3-4007	7'	adjusts 6"/side	26"-39"
3-4010	10'	adjusts 6"/side	26"-39"
3-4012	12'	adjusts 6"/side	26"-39"

Floor Mounted Parallel Bars

- Stainless steel adjustable handrails
- Same features as platform-mounted model, except it is floor-mounted
- 400 lbs. load capacity



	Length	Width Range	Height Range
3-3317	7'	16"-24"	22"-35"
3-3310	10'	16"-24"	22"-35"

Folding Parallel Bars

- Stainless steel handrails
- Unique double stabilizer bars stiffen the rails
- Spring-loaded plungers and sure lock positioning knobs
- 350 lbs. load capacity

PARALLEL BARS



	Length	Width	Height Range
3-5007	7'	20"	26"-39"
3-5010	10'	20"	26"-39"
3-5012	12'	20"	26"-39"

Hemiplegic Parallel Bars

- Floor mounted bars allow patients an uninterrupted bar to grasp with one hand
- 350 lbs. load capacity



-045

Optional Center Obduction Board

- Available with platform style only

SELECT SERIES STAIRCASES

FEATURES

- 1 1/2" diameter stainless steel handrails with end bumpers
- Fast easy height adjustments with self indexing spring action plungers and numbered height indicator strips (no holes to line up or count)
- Second knob removes any side-to-side play for a rigid upright
- Hardwood plywood base with durable, UV satin, topcoat finish
- Assembled with steel fasteners (No glued wood dowels)
- All edges are banded with tough, Slate Gray, PVC edge material
- 2 full width safety treads on each step
- Available in 30" and 36" widths
- 450 lbs. weight limit (both models)
- Reinforced bariatric load capacities



4-6500-30	Length	Width	Height
platform	24"	30"	12"
2 steps	10"	30"	6"
3 steps	10"	30"	4"
overall size	54"	34"	

Adjustable Small Staircase

- Adjusts to all ages from children to adult in 1 1/2" increments
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling



Adjust uprights by number strip on the uprights



Spring-loaded plungers make adjustments easy

4-6500-36	Length	Width	Height
platform	24"	36"	12"
2 steps	10"	36"	6"
3 steps	10"	36"	4"
overall size	54"	40"	

Adjustable Small Staircase

- Same as above with 36" width

4-6501-30	Length	Width	Height
platform	24"	30"	24"
each step	10"	30"	6"
overall size	54"	34"	

Adjustable One-Sided Staircase

- Adjusts in 1 1/2" increments to accommodate all ages, from children to adults
- Requires a minimum 9' ceiling

4-6501-36	Length	Width	Height
platform	24"	36"	24"
each step	10"	36"	6"
overall size	54"	40"	

Adjustable One-Sided Staircase


- Same as above with 36" width



CLASSIC SERIES STAIRCASES

FEATURES

- All hardwood plywood with ultra-durable, UV-cured, satin finish base
- All base edging finished with tough Slate Gray PVC edging (no exposed plywood core)
- Ultra-rigid, steel uprights with durable finish
- Handrails and uprights are pre-assembled
- Each model is available in 30" and 36" widths
- Completely assembled with steel fasteners (no glued dowel pins)
- Reinforced bariatric capacity
- All steps are a full 10" deep
- Built for years of use



* Available with optional Bus Step. See page 30.

CLASSIC SERIES STAIRCASES



4-5501-30	Length	Width	Height
each step	10"	30"	6"
platform	24"	30"	24"
overall size	56"	34"	

One-Sided Staircase

- Small corner staircase for limited space
- Requires a minimum 9' ceiling
- 500 lbs. load capacity

4-5501-36	Length	Width	Height
each step	10"	36"	6"
platform	24"	36"	24"
overall size	56"	40"	

One-Sided Staircase

- Same as above with 36" width



Straight Position



Corner Position

4-5500-30	Length	Width	Height
1 step	10"	30"	6"
2 steps	10"	30"	4"
platform	24"	30"	12"
overall size	54"	34"	

Small Up/Down Staircase

- Small straight staircase for limited space
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling
- 500 lbs. load capacity

4-5500-36	Length	Width	Height
1 step	10"	36"	6"
2 steps	10"	36"	4"
platform	24"	36"	12"
overall size	54"	40"	

Small Up/Down Staircase

- Same as above with 36" width

4-5080-30	Length	Width	Height
3 steps	10"	30"	6"
5 steps	10"	30"	4"
platform	30"	30"	24"

Convertible Staircase

- Converts from straight to corner arrangement
- 500 lbs. load capacity
- Requires a minimum 9' ceiling

4-5080-36	Length	Width	Height
3 steps	10"	36"	6"
5 steps	10"	36"	4"
platform	36"	36"	24"

Convertible Staircase

- Same as above with 36" width

VALUE SERIES STAIRCASES

FEATURES

- Comfort grip moulded hardwood handrails
- Hardwood, plywood base with durable, satin, topcoat finish
- All stair edges are banded with tough Slate Gray, PVC edge material
- 2 full-width safety treads on each step
- Rigid grip handrail
- Fully enclosed sides helps give patient a secure feeling
- Assembled with steel fasteners (No glued dowel pins)
- 350 lbs. weight limit



	Length	Width	Height
4-7501	56"	34"	58"
platform	24"	30"	24"
each step	10"	30"	6"

One-Sided Staircase with Full Sides

- Small corner staircase with three enclosed sides
- Requires a minimum 9' ceiling



	Length	Width	Height
4-7500	54"	34"	36"
platform	24"	30"	12"
1 step	10"	30"	6"
2 steps	10"	30"	4"

Small Straight Staircase with Full Sides

- Small straight staircase for limited space
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling



101

Bus Step Option

- Creates a 12" high step when in place
- Available option for any staircase with a 6" rise step
- MUST BE ORDERED WITH STAIRCASE

STAIRCASE STORAGE ISLANDS

Door Closed



Front



Door Open



Back



4-8010	Length	Width	Height
each step	10"	30"	6"
platform	24"	30"	24"
overall size	56"	34"	

Staircase Storage Island

- All hardwood plywood construction
- Small corner staircase for limited space
- Requires a minimum 9' ceiling
- Safety tread on all steps and platform
- 2-door storage, 1 on each side with an adjustable shelf for each
- 500 lbs. load capacity

4-8020	Length	Width	Height
each step	10"	30"	6"
platform	24"	30"	24"
overall size	63"	39"	

Staircase Storage Island

- All hardwood plywood construction staircase
- Small corner staircase for limited space
- Requires a minimum 9' ceiling
- Safety tread on all steps and platform
- All laminate bandrac and dumbbell and cuff strips
- 500 lbs. load capacity
- Holds up to 48 cuff weights, 22 dumbbells and 5 rolls of exercise band and/or tubing (not included)