

*Quality Care is at Your Fingertips*



---

## Surgical Sponges & Towels

---



## DEROYAL SURGICAL SPONGES & TOWELS

At DeRoyal, we know that clinician product needs can be exacting and variable, from gross dissection of tissue to breath-taking micro-surgical precision. DeRoyal Sponges and Towels span virtually every aspect of surgery – a diverse line using materials from a wide range of sources. Every time you pick up a DeRoyal Sponge or Towel, you hold quality in your hands and have excellence at your fingertips.

### CONTENT

#### Product Listing

	Page
Atraumatic Lint-Free Sponges.....	3 – 5
Ophthalmic Sponges.....	6
Dissecting Sponges.....	7 – 8
Tonsil and Adenoid Sponges.....	9
Surgical Gauze Sponges.....	10
Sterile Puffs.....	11
Miscellaneous Sponges.....	12 – 13
Neurology Sponges.....	14
Laparotomy Gauze Sponges.....	15
Gauze Sponges.....	16
Operating Room Towels.....	16 – 17
X-Ray Detectable Safety Towels.....	18
Sponge Management.....	19

# Quality Care is at Your Fingertips

## ATRAUMATIC LINT-FREE SPONGES

### Sinus Sponges

- Helps prevent lateralization of the middle turbinate
- Significantly reduces renewed ostiomeatal obstruction following FESS (Functional Endoscopic Sinus Surgery)
- Compressed configuration for easy insertion
- Expands anatomically so antibiotics stay in contact with surgical site
- Micro-pore size minimizes tissue growth
- Clears clots and debris upon removal, leaving a clean, open cavity
- Sterile
- Latex safe

Product No.	Description	Quantity
34-301	Sinus Dressing 3.5 cm long x 0.9 cm thick x 1.2 cm high	2/Pk, 10 Pks/Cs
34-302	Slim Sinus Dressing 3.5 cm long x 0.6 cm thick x 1.2 cm high	2/Pk, 10 Pks/Cs
34-303	Sinus Dressing 2.5 cm long x 0.9 cm thick x 2 cm high	2/Pk, 10 Pks/Cs
34-304	Slim Sinus Dressing 2.5 cm long x 0.6 cm thick x 2 cm high	2/Pk, 10 Pks/Cs
34-305	Sinus Dressing Kit 4.5 cm long x 1.3 cm thick x 2 cm high w/Airway Tube & Nasal Pack 7.5 cm long x 1.5 cm thick x 3.5 cm high	1 Kit/Pk, 10 Pks/Cs
34-306	Sinus Dressing 4.5 cm long x 0.9 cm thick x 2 cm high	2/Pk, 10 Pks/Cs
34-307	Slim Sinus Dressing 4.5 cm long x 0.6 cm thick x 2 cm high	2/Pk, 10 Pks/Cs



# Surgical Sponges & Towels

## ATRAUMATIC LINT-FREE SPONGES

### Nasal Sponges

- Provides gentle, even support
- Absorbs post-operative drainage following nasal surgery, ethmoidectomy procedures, or other endoscopic sinus surgeries
- Fiber free to reduce risk of adhesions or active rebleeding upon removal
- Absorbent to control post-nasal discharge – holds three times more fluid than gauze
- Optional airway tube aids in patient comfort for better breathing and pressure reduction
- Sterile
- Latex safe

Product No.	Description	Quantity
34-101	Small Standard Nasal Dressing w/String 4.5 cm long x 1.5 cm thick x 2 cm high	1/Pk, 10 Pks/Cs
34-102	Large Standard Nasal Dressing w/String 8 cm long x 1.5 cm thick x 2 cm high	1/Pk, 10 Pks/Cs
34-103	Small Standard Nasal w/Airway 4.5 cm long x 1.5 cm thick x 2 cm high	1/Pk, 10 Pks/Cs
34-104	Large Standard Nasal w/Airway 8 cm long x 1.5 cm thick x 3 cm high	1/Pk, 10 Pks/Cs
34-105	Anatomical Nasal Dressing w/String 8 cm long x 1.5 cm thick x 3 cm high	1/Pk, 10 Pks/Cs
34-106	Anatomical Nasal Dressing w/Airway & String 8 cm long x 1.5 cm thick x 3 cm high	1/Pk, 10 Pks/Cs
34-107	Slim Anatomical Nasal Dressing w/String 8 cm long x 0.75 cm thick x 3 cm high	1/Pk, 10 Pks/Cs



**ATRAUMATIC LINT-FREE SPONGES**

**Epistaxis Sponges**

- Provides quick and effective control of epistaxis
- Provides unsurpassed patient comfort
- Atraumatic removal prevents active rebleeding
- Compressed configuration slides easily into nasal cavity
- Absorbs three times more fluid than gauze
- Sterile
- Latex safe

Product No.	Description	Quantity
34-201	Large Epistaxis Dressing w/String 10 cm long x 1.5 cm thick x 2.5 cm high	1/Pk, 10 Pks/Cs
34-202	Small Epistaxis Dressing w/String 5.5 cm long x 1.5 cm thick x 2.5 cm high	1/Pk, 10 Pks/Cs



**Otology Sponges**

- Helps prevent stenosis of the external ear canal following otic surgery
- Compressed configuration allows easy insertion into the ear canal
- Keeps antibiotics in contact with the wound site through anatomical expansion
- Fiber free to minimize tissue ingrowth
- Sterile
- Latex safe

Product No.	Description	Quantity
34-401	Otitis External Wick, 9 mm x 15 mm	1 Pk, 10 Pks/Cs
34-402	Otitis External Wick, 9 mm x 24 mm	1 Pk, 10 Pks/Cs
34-403	Ear Packing, 12 mm x 15 mm	1 Pk, 10 Pks/Cs




# Surgical Sponges & Towels

## OPHTHALMIC SPONGES

### Lint-Free Ophthalmic Sponges

- Allows the ophthalmologist to perform completely lint-free fluid control
- A variety of premium absorbent lint-free ophthalmic sponges
- Extremely soft and fast-wicking
- Sterile
- Latex safe

Product No.	Description	Quantity
34-501	Corneal Shield, 7 mm diameter	1/Pk, 10 Pks/Cs
34-502	Cataract Drain, 140 cc	1/Pk, 10 Pks/Cs
34-503	Phaeco Drain, 400 cc	1/Pk, 10 Pks/Cs
34-504	Eye Wick, 8"	1/Pk, 10 Pks/Cs
34-505	Diamond Knife Cleaning Block	1/Pk, 10 Pks/Cs
34-506	Wick & Wipe Instrument Wipe w/Wick	1/Pk, 10 Pks/Cs
34-507	Eye Kit Wick & Wipe w/Corneal Light Shield	1/Pk, 10 Pks/Cs
25-201	Standard Instrument Wipe (1/Pk)	20 Pks/Bx, 6 Bxs/Cs



### Eye Spears

- Designed for microscopic sponging in the ocular area
- Fine-pore, cellulose fiber is virtually non-linting and extremely absorbent
- Fiber is compressed and die-cut into a "spear" shape and attached to a small plastic handle
- Sterile
- Latex safe

Product No.	Description	Quantity
30-049	Eye Spears	6 Pk, 10 Pks/Cs
30-049-6	Eye Spears	6 Pk, 50 Pks/Cs
30-049-10	Eye Spears	10 Pk, 50 Pks/Cs





**DISSECTING SPONGES**

**Rosebud Dissector**

- Excellent for blunt dissection and delicate surgery
- Smaller and firmer than a Peanut or Round Dissector
- Complete with the exclusive x-ray detectable radiopaque element
- Count holder, per pack, enables easy sponge count
- Sterile
- Latex safe

Product No.	Description	Quantity
30-202	1/4" in Foam Holder	5/Pk, 40 Pks/Cs



**Peanut Sponge**

- Excellent for delicate sponging and soft tissue dissection
- Soft, round construction may be used as is, or folded
- Complete with the exclusive x-ray detectable radiopaque element
- Count holder, per pack, enables easy sponge count
- Sterile
- Latex safe

Product No.	Description	Quantity
30-106	3/8" in Foam Holder	5/Pk, 50 Pks/Cs
30-010	Double Strung	5 Pk, 100 Pks/Cs



**Cherry Sponge**

- Largest of the Dissectors
- Very absorbent and ideal for general dissecting
- Complete with the exclusive x-ray detectable radiopaque element
- Count holder, per pack, enables easy sponge count
- Sterile
- Latex safe

Product No.	Description	Quantity
30-112	1/2" in Foam Holder	5/Pk, 50 Pks/Cs



# Surgical Sponges & Towels

## DISSECTING SPONGES

### Round (Hard) Dissector

- Also known as the "Beehive"
- Roundness and firmness make our hard dissectors beneficial for blunt dissection and ideal for surgical procedures
- Small hole in center of dissector makes it easier to grasp with forceps
- Complete with the exclusive x-ray detectable radiopaque element
- Count holder, per pack, enables easy sponge count
- Sterile
- Latex safe

Product No.	Description	Quantity
30-115	$\frac{3}{8}$ " in Foam Holder	5/Pk, 50 Pks/Cs



### Kittner (Blunt) Dissector

- Can be used as a longitudinal or folded sponge for dissection
- Constructed of high quality cotton twill
- Complete with the exclusive x-ray detectable radiopaque element
- Count holder, per pack, enables easy sponge count
- Sterile
- Latex safe

Product No.	Description	Quantity
30-101	$\frac{5}{16}$ " x $\frac{1}{4}$ " in Foam Holder	5/Pk, 50 Pks/Cs





# Quality Care is at Your Fingertips

## TONSIL & ADENOID SPONGES

### Single-Strung Tonsil Sponges

- Constructed with a 15", white, memory-free cotton string
- Covered with gauze and cotton filled for maximum absorbency
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-036	Small ¾"	5/Pk, 100 Pks/Cs
30-037	Medium 7/8"	5/Pk, 100 Pks/Cs
30-038	Large 1"	5/Pk, 100 Pks/Cs
30-039	X-Large 1¼"	5/Pk, 100 Pks/Cs



### Double-Strung Tonsil Sponges

- Constructed with a No. 3 blue thread for ease of retrieval
- Covered with gauze and cotton filled for maximum absorbency
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-040	Small ¾"	5/Pk, 100 Pks/Cs
30-041	Medium 7/8"	5/Pk, 100 Pks/Cs
30-042	Large 1"	5/Pk, 100 Pks/Cs
30-043	X-Large 1¼"	5/Pk, 100 Pks/Cs



### Round Adenoid Sponge

- Designed to control hemorrhaging following adenoid surgery
- Gauze-cotton filled with a No. 3 blue thread
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-044	Standard 1½"	5/Pk, 50 Pks/Cs



# Surgical Sponges & Towels

## SURGICAL GAUZE SPONGES

### Flat Stick Sponges

- Pre-folded with the edges tucked in
- May be used with a Forester clamp for deep stick sponging
- Excellent for absorption
- Ideal for prepping prior to surgery
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-022	1½"	10/Pk, 50 Pks/Cs
30-222	2"	10/Pk, 50 Pks/Cs
30-022NS	1½", Non-sterile	10/Pk, 50 Pks/Cs
30-222NS	2", Non-sterile	10/Pk, 50 Pks/Cs



### Round Stick Sponges

- Excellent absorbency rate
- Can be used as a prep sponge for surgery
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-026	Small ¾"	10/Pk, 50 Pks/Cs
30-025	Medium 1"	10/Pk, 50 Pks/Cs
31-567	Large 1¼"	10/Pk, 50 Pks/Cs
31-807	Large 1¼"	3/Pk, 50 Pks/Cs
30-122	X-Large 1½"	5/Pk, 50 Pks/Cs
30-024	X-Large 1½"	10/Pk, 50 Pks/Cs



### Triangle Sponge

- Ideally shaped for sponging
- May be easily rolled or folded for dissection or used as a prepping sponge
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-018	Medium	10/Pk, 50 Pks/Cs
30-019	Large	10/Pk, 50 Pks/Cs
30-020	Medium, Non-sterile	Bulk, 500/Cs
30-021	Large, Non-sterile	Bulk, 500/Cs



# Quality Care is at Your Fingertips

## STERILE PUFFS

### Sterile Cotton Balls

- Ideal for prepping, wound cleaning, and perineal care
- Sterile
- Latex safe

Product No.	Description	Quantity
31-090	Medium 1"	10/Pk, 100 Pks/Cs
30-030	Large 1 1/4"	5/Pk, 100 Pks/Cs



### Double-Strung Cotton Balls

- Constructed with a 15" blue silk thread for ease of retrieval
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-031	Tiny 1/4"	5/Pk, 50 Pks/Cs
30-176	X-Small 3/8"	5/Pk, 50 Pks/Cs
30-032	Small 1/2"	5/Pk, 50 Pks/Cs
30-033	Medium 3/4"	5/Pk, 50 Pks/Cs
30-034	Large 1"	5/Pk, 50 Pks/Cs
30-035	X-Large 1 1/4"	5/Pk, 50 Pks/Cs



### Sterile Rayon Balls

- 100% pure medical grade
- Provides the highest absorbency
- Virtually lint-free, uniformly shaped, and soft
- Will not cling to each other, for ease of use
- Sterile
- Latex safe

Product No.	Description	Quantity
31-937	3/4" Medium	5/Pk, 250 Pks/Cs
30-224	1 1/4" X-Large	5/Pk, 100 Pks/Cs



# Surgical Sponges & Towels

## MISCELLANEOUS SPONGES

### Double-Strung Cylindrical Sponges

- Highly absorbent
- Designed for tonsil, adenoid, and nasal surgery
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-045	1½" x ½"	5/Pk, 100 Pks/Cs



### Dental Rolls

- Highly absorbent
- Designed to keep the area dry, absorbing saliva and moisture
- Available in regular and quilted
- Sterile
- Latex safe


Product No.	Description	Quantity
30-096	Regular 1½" x ¾"	5/Pk, 100 Pks/Cs
30-419	Regular 3" x ¾"	4/Pk, 100 Pks/Cs
31-207	Quilted 1½" x ¾"	5/Pk, 25 Pks/Cs
32-1071	Regular 6" x ¾"	2/Pk, 100 Pks/Cs



### Ear Wicks (Cellulose)

- Atraumatic and virtually lint-free
- Designed to meet the surgeon's need in otolaryngology surgery
- Compressed and die-cut of highly absorbent fine-pore cellulose fiber
- Sterile
- Latex safe

Product No.	Description	Quantity
30-048	Ear Wicks	5/Pk, 20 Pks/Cs




## MISCELLANEOUS SPONGES

### Tracheotomy Sponge

- Made of highly absorbent 4 x 4 12-ply gauze sponge
- Ideal around any drainage area to protect from chafing
- Sterile
- Latex safe

Product No.	Description	Quantity
30-051	Tracheotomy Sponge	25/Pk, 12 Pks/Cs
30-052	Tracheotomy Sponge	25/Pk, 20 Pks/Cs



### Bronchial Sponge

- Highly absorbent cotton
- Complete with the exclusive x-ray detectable radiopaque element
- Sterile
- Latex safe

Product No.	Description	Quantity
30-205	Bronchial Sponge	5/Pk, 50 Pks/Cs



### Nasal Plugs

- Ideally shaped for easy insertion and comfortable use following surgery
- Sterile
- Latex safe

Product No.	Description	Quantity
32-267	Nasal Plugs	5/Pk, 250 Pks/Cs
30-047	Nasal Plugs	2/Pk, 500 Pks/Cs



# Surgical Sponges & Towels

## NEUROLOGY SPONGES

### Neuro Sponges in Patties & Strips


- High-quality rayon fibers prevent disintegration and linting, reducing the risk of granulomas
- Very pliable when wet, allowing patties and strips to conform to irregular contours
- X-Ray detectable element is secured on one side of pad
- Locator string attached to the patties increases visibility and assists in identification
- Patties come with counting card in each pack
- Sterile
- Latex safe

#### PATTIES

Product No.	Description	Quantity
30-054	1/2" x 1/2" Patties	10/Pk, 20 Pks/Cs
30-055	1/2" x 1" Patties	10/Pk, 20 Pks/Cs
30-056	1/2" x 1 1/2" Patties	10/Pk, 20 Pks/Cs
30-057	1/2" x 3" Patties	10/Pk, 20 Pks/Cs
30-057T	1/2" x 3" Patties w/Telfa®	10/Pk, 20 Pks/Cs
30-058	3/4" x 3/4" Patties	10/Pk, 20 Pks/Cs
30-059	1" x 1" Patties	10/Pk, 20 Pks/Cs
30-060	1" x 3" Patties	10/Pk, 20 Pks/Cs
30-061	1 1/2" x 1 1/2" Patties	10/Pk, 20 Pks/Cs
30-063	3" x 3" Patties	10/Pk, 20 Pks/Cs
30-067	1/4" x 1 1/2" Patties	10/Pk, 20 Pks/Cs
30-068	1/2" x 2" Patties	10/Pk, 20 Pks/Cs
30-069	1/4" x 3" Patties	10/Pk, 20 Pks/Cs
30-069T	1/4" x 3" Patties w/Telfa®	10/Pk, 20 Pks/Cs
30-072	3/4" x 3" Patties	10/Pk, 20 Pks/Cs
30-154	1/4" x 1/4" Patties	10/Pk, 20 Pks/Cs

#### STRIPS

Product No.	Description	Quantity
30-300	1/4" x 6" Strips	10/Pk, 20 Pks/Cs
30-301	1/2" x 6" Strips	10/Pk, 20 Pks/Cs
30-302	1" x 6" Strips	10/Pk, 20 Pks/Cs
30-303	3/4" x 6" Strips	10/Pk, 20 Pks/Cs
30-304	1 1/2" x 6" Strips	10/Pk, 20 Pks/Cs
30-305	3" x 6" Strips	10/Pk, 20 Pks/Cs
30-306	2" x 6" Strips	10/Pk, 20 Pks/Cs
30-307	3 1/2" x 6" Strips	10/Pk, 20 Pks/Cs



\* Telfa is a registered trademark of Tyco HealthCare Group, L.P.

# Quality Care is at Your Fingertips

## LAPAROTOMY GAUZE SPONGES

### Cotton Lap Sponges

- Consists of 100% cotton
- X-ray detectable
- Pre-washed consisting of 24/20 gauze
- Wrapped in CSR and placed in peel pouch for easy sterile presentation
- Manufactured in Mexico, a NAFTA partner country
- Sterile
- Latex safe

Product No.	Description	Quantity
1-418	4" x 18" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-836	8" x 36" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-8108	8" x 108" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-1212	12" x 12" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-1818	18" x 18" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-1818R	18" x 18" Pre-washed, w/Rings	5/Pk, 20 Pks/Cs
1-1836-5	18" x 36" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
1-1836-3	18" x 36" Pre-washed, w/o Rings	3/Pk, 25 Pks/Cs
1-1836	18" x 36" Pre-washed, w/o Rings	2/Pk, 25 Pks/Cs



### Lap Sponges

- Pre-washed consisting of 24/20 gauze
- Economically priced for a variety of uses
- Sterile
- Latex safe

Product No.	Description	Quantity
9-1818	18" x 18" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
200-1818	18" x 18" Pre-washed, w/o Rings	5/Pk, 40 Pks/Cs
200-1818R	18" x 18" Pre-washed, w/ Rings	5/Pk, 40 Pks/Cs
9-836	8" x 36" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
9-1212	12" x 12" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
200-1212	12" x 12" Pre-washed, w/o Rings	5/Pk, 40 Pks/Cs
9-418	4" x 18" Pre-washed, w/o Rings	5/Pk, 20 Pks/Cs
200-418	4" x 18" Pre-washed, w/o Rings	5/Pk, 40 Pks/Cs
200-418R	4" x 18" Pre-washed, w/Rings	5/Pk, 40 Pks/Cs





# Surgical Sponges & Towels

## GAUZE SPONGES

### Vaginal & Perineal Packing

- Pre-washed, 4-ply gauze
- Vaginal packing complete with the exclusive x-ray detectable radiopaque element
- Folded and packaged in peel pouch
- Sterile
- Latex safe

Product No.	Description	Quantity
1-136	Vaginal Packing, w/X-Ray Element, 1" x 36"	1/Pk, 100 Pks/Cs
1-237	Vaginal Packing, w/X-Ray Element, 2" x 36"	1/Pk, 50 Pks/Cs
1-272	Vaginal Packing, w/X-Ray Element, 2" x 72"	1/Pk, 50 Pks/Cs
1-424	Vaginal Packing, w/X-Ray Element, 4" x 96"	1/Pk, 20 Pks/Cs
1-472	Vaginal Packing, w/X-Ray Element, 4" x 72"	1/Pk, 100 Pks/Cs
1-412	Vaginal Packing, w/X-Ray Element, 4" x 48"	1/Pk, 20 Pks/Cs
1-436	Vaginal Packing, w/X-Ray Element, 4" x 36"	1/Pk, 100 Pks/Cs
1-250	Perineal Packing, w/o X-Ray Element, 2" x 5 Yds	1/Pk, 50 Pks/Cs
1-230	Perineal Packing, w/o X-Ray Element, 2" x 3 Yds	1/Pk, 50 Pks/Cs
1-150	Perineal Packing, w/o X-Ray Element, 1" x 5 Yds	1/Pk, 25 Pks/Cs
1-554	Perineal Packing, w/o X-Ray Element, 1/2" x 54"	1/Pk, 50 Pks/Cs



Vaginal Packing



Vaginal Packing w/Loop



Perineal Packing




Perineal Packing w/Loop

## OPERATING ROOM TOWELS

### Fenestrated Towels

- High-quality "huck" towels
- Pre-washed and de-linted
- 26" x 17"
- 58/25 mesh
- 65 grams by weight
- Colors will not bleed
- 100% cotton
- Sterile
- Latex safe

Product No.	Description	Quantity
70-150	Towel, White, 1 1/2" Fenestration	1/Pk, 20 Pks/Cs
32-082	Towel, Blue, 2 1/2" Fenestration	1/Pk, 50 Pks/Cs
70-103	Towel, White, 3" Fenestration	1/Pk, 20 Pks/Cs
70-102	Towel, White, 3 1/2" Fenestration	1/Pk, 80 Pks/Cs
32-998	Towel, Blue, 30" x 34" w/5" x 10" Fenestration	1/Pk, 40 Pks/Cs



# Quality Care is at Your Fingertips

## OPERATING ROOM TOWELS

### O.R. Towels

- Pre-washed
- Low linting
- Manufactured in Mexico, a NAFTA partner country
- Sterile
- Latex safe

	Blue	
Qty/Pack	Product No.	Qty/Case
1/Pk	60-102	80 Pks/Cs
2/Pk	60-106	40 Pks/Cs
4/Pk	60-101	20 Pks/Cs
6/Pk	60-103	12 Pks/Cs
8/Pk	60-100	10 Pks/Cs
10/Pk	60-107	8 Pks/Cs
12/Pk	60-109	6 Pks/Cs



### Standard O.R. Towels

- High-quality "huck" towels
- Pre-washed, de-linted
- 26" x 16"
- 58/25 mesh
- 65 grams by weight
- Colors will not bleed
- 100% cotton
- Sterile
- Latex safe

	Blue	Green	White	
Qty/Pack	Product No.	Product No.	Product No.	Qty/Case
1/Pk	63-102	63-202	61-402	80 Pks/Cs
2/Pk	63-106	63-206	61-406	40 Pks/Cs
4/Pk	63-101	63-201	61-401	20 Pks/Cs
5/Pk	63-108	63-208	61-408	16 Pks/Cs
6/Pk	63-103	63-203	61-403	12 Pks/Cs
8/Pk	63-100	63-200	61-400	10 Pks/Cs
10/Pk	63-107	63-207	61-407	8 Pks/Cs
12/Pk	63-109	N/A	N/A	6 Pks/Cs
30/Pk	63-110	N/A	N/A	2 Pks/Cs
Bulk Non-sterile	63-001	63-002	63-003	100/Cs
Bulk Sterile	63-103B	N/A	N/A	6 Pks/Cs



# Surgical Sponges & Towels

## X-RAY DETECTABLE SAFETY TOWELS

### O.R. Safety Towels

- Preferred for wound packing by nurses
- X-ray detectable radiopaque element
- Pre-washed
- Low-linting
- Manufactured in Mexico, a NAFTA partner country
- Sterile
- Latex safe

	Blue	White	
Qty/Pack	Product No.	Product No.	Qty/Case
1/Pk	N/A	68-402	80 Pks/Cs
2/Pk	68-106	68-406	40 Pks/Cs
4/Pk	68-101	68-401	20 Pks/Cs
5/Pk	68-108	68-408	16 Pks/Cs
6/Pk	68-103	68-403	12 Pks/Cs



### Standard O.R. Safety Towels

- X-ray detectable
- Pre-washed
- Sterile
- Latex safe

	Blue	
Qty/Pack	Product No.	Qty/Case
2/Pk	67-106	40 Pks/Cs
4/Pk	67-101	20 Pks/Cs
5/Pk	67-108	16 Pks/Cs
6/Pk	67-103	12 Pks/Cs



## SPONGE MANAGEMENT

### Sponge Counter Bag

- Allows continuous sponge counting and containment to minimize exposure to bloodborne pathogens
- Easy open pockets for the simple disposal of used sponges
- Accommodates five Laparotomy sponges or ten 4 x 4s
- Clear polyethylene pockets allow sponges to be viewed for visual estimate of patient blood loss
- Additional bags can be placed over filled ones as necessary



Product No.	Description	Quantity
25-500	Sponge Counter Bags, Clear Back, Non-Sterile	50 Bgs/Crtn, 5 Crtns/Cs
25-551	Sponge Counter Bags, Blue Back w/Adhesive Seal, Non-Sterile	50 Bgs/Crtn, 5 Crtns/Cs
25-501	Mounting Bracket	1 each