

### Patient Protection Set

- Set or order the individual aprons
- Four sizes, from large adult to infant
- .50mm Pb equivalent.
- Infab's exclusive belt with buckle closure
- Space efficient rack
- Solid aluminum, powder coated rack
- Rack dimensions - 13" wide and 8" deep



Model	Description
683400	Storage Rack with 4 Aprons (Toddler, Child, Average & Large Adult)
683401	8" W X 10" Toddler Apron Only
683402	12" W X 12" Child's Apron Only
683403	14" W X 16" Average Adult Apron Only
683404	16" W X 20" Large Adult Apron Only
683405	Storage Rack with 1 Average Adult Apron and 1 Large Adult Apron
683406	Storage Rack Only

### Deluxe Mobile Apron Rack

- Heavy gauge steel
- Powder coated
- Durable and long lasting finish
- Hospital grade casters
- Smooth rolling, even under heavy loads
- Holds up to ten protection aprons
- Optional vest and skirt adapter (Part Number 683000)
- Optional glove rack - three pair of gloves (shown)
- Shipped assembled

Model	Description
683407	Deluxe Mobile Apron Rack
683408	Deluxe Mobile Apron Rack With Glove Rack



### Mini Mobile Apron Rack

- Heavy gauge cold rolled steel
- Powder coated
- Attractive, durable, long lasting finish
- Hospital grade casters
- Holds up to six protection aprons, a pair of lead gloves and two pairs of radiation protection glasses
- Compact size (78" H x 23" W x 24" D)
- Aprons and other accessories available separately

Model	Description
683418	Mini Mobile Apron Rack
683419	Mini Mobile Apron Rack with six model 683450 Hangers included



## Wall Mounted Apron Rack

- Heavy gauge steel
- Powder coated
- Scratch resistant
- Easy to install
- Requires only minimal space
- 36" wide x 8" tall
- Extends from the wall only 26"
- Stores up to 10 aprons
- Optional vest and skirt adapter (Part Number 683000)



Model	Description
683429	10 Arm Swivel Apron Rack Wall Mounted

## Budget Saver Mobile Apron Rack

- Heavy duty
- Economically priced
- Heavy gauge steel
- Powder coated
- Durable, hospital grade casters
- Stores up to ten protection aprons
- Optional vest and skirt adapter (Part Number 683000)


Model	Description
683427	Budget Mobile Apron Rack
683428	Budget Mobile Apron Rack With Glove Rack



**Deluxe 5 Arm Apron Rack**

- Store five aprons
- Immediate access to individual aprons
- Heavy duty
- Coated steel mounting bracket
- Arms with nylon bushings
- Racks fold flush to wall
- Swing out for easy access
- Racks can be mounted side by side
- Optional vest and skirt adapter (Part Number 683000)

Model	Description
683100	5 Arm Folding Right Hand Folding
683101	5 Arm Folding Left Hand Folding
683105	Weak Wall Support/8-3/4" X 20"
683000	Vest and Skirt Rack Adapter



**Wall Tower Rack for 683100 and 683101 racks**

- Heavy duty
- Handles toughest hospital environments
- Heavy gauge steel
- Powder coated
- Scratch resistant
- Retain its good looks for years
- Easy to install
- Hardware provided
- Requires only minimal space
- 18" wide x 63" tall



Model	Description
683115	Weak Wall Support
683415	Rack/Tower W/ R.H. & L.H.Racks

### Apron & Glove Combo Rack

- Quick, convenient storage
- Sturdy
- Compact
- Inexpensive unit
- Mounts quickly and easily
- Heavy gauge steel mounting plate
- Half inch outside diameter tubing
- Electro static powder coated



Model	Description
683420	Lead Apron and Glove Storage Rack

### Apron and Glove Combo Rack

- Rounded corners to prevent accidental apron tears or skin punctures
- Quick, convenient storage
- Single apron and a pair of lead gloves
- Sturdy
- Compact
- Inexpensive unit
- Mounts quickly and easily
- Heavy gauge steel mounting plate
- Half inch outside diameter tubing
- Electro static powder coated



Model	Description
683430	Apron and Glove Combination Rack

### Protective Accessories Rack

- Storage unit for two or more aprons, a pair of glasses and a pair of lead gloves
- Mounts quickly and easily
- 18" wide
- Extends a mere 4.5" from the wall
- Contains no sharp edges
- Perfectly suited to small areas
- Powder coated heavy steel construction



Model	Description
683102	Protective Accessory Rack Only

**Six or Eight Peg Apron Rack**



- Available in 2 peg, 6 peg or 8 peg configurations
- Inexpensive but very effective
- Heavy gauge steel plate
- Electro static powder coated
- Tilted storage bars, aprons can't slip off
- Requires minimal space
- Mounts to existing walls



- Add Locking System to your Apron Rack by adding -LS to part number



Model	Description
683412	Two Peg Wall Mounted Rack For Aprons
683410	Six Peg Wall Mounted Rack For Aprons
683411	Eight Peg Wall Mounted Rack For Aprons

**Deluxe Apron Hanger**

- Securely hold a radiation protection apron
- Constructed from heavy gauge steel tube
- Powder coated
- "No slip" tubing shield



Model	Description
683450	Deluxe Custom Apron Hanger