

Radiation Protection Glasses

Where function meets fashion!



ISOLATE YOUR EYES
COCOONS



nicole miller

bebe
eyes

DAKOTA SMITH
EYEWEAR

An essential part of a complete personal radiation protection plan, Infab Corporation glasses are manufactured in the USA and are available in a wide range of styles and models, from Industrial type safety frames, to top designers fashion frames, and all with your choice of .50mm or .75mm Pb Equivalent lenses.

Most frames can be fitted with your prescription from single vision to Bi-Focal (extra charges apply) or choose from a wide variety of "Fitover" style glasses that can be comfortably worn over your existing prescription glasses.

All Infab brand radiation protection glasses come with their own high quality carrying case and are lightweight, comfortable, durable and offer essential protection all at a truly affordable price and in stock for immediate delivery.

Note:

Although often thought to be a personal choice, OSHA regulation 1910.132(a) actually states:

"Protective equipment, including personal protective equipment for eyes, face, head, and extremities, protective clothing, respiratory devices, and protective shields and barriers, shall be provided, used, and maintained in a sanitary and reliable condition wherever it is necessary by reason of hazards of processes or environment, chemical hazards, radiological hazards, or mechanical irritants encountered in a manner capable of causing injury or impairment in the function of any part of the body through absorption, inhalation or physical contact". (Emphasis added)

Team



Team is a stylish wrap around safety frame made of TR-90 high impact resistant plastic. ANSI Z87.1+ Approved. Lightweight and comfortable, can be fit with most prescriptions and includes a high quality carrying case. "Team" glasses are .75 Pb Lead Equivalency, and are available from stock in non prescription in six great colors.

Model	Color	Weight
860001	Black	48.5g/1.7 oz
860002	Blue/Yellow	48.5g/1.7 oz
860003	Silver/Purple	48.5g/1.7 oz
860004	Red	48.5g/1.7 oz
860005	Silver	48.5g/1.7 oz
860006	Taupe	48.5g/1.7 oz

Achilles



The Achilles Goggle is a vented radiation protection goggle with neoprene foam on the inside specifically designed to protect against fluid splashing and perspiration on your face. The elastic adjustable head band makes this frame easy to adjust and extremely comfortable to wear even for extended periods. NOTE: This is a larger frame design is not recommended for petite face shapes and size.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Four great colors
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
800075	Red	57.4g / 2.02oz.
800076	Blue	57.4g / 2.02oz.
800077	Plated Silver	57.4g / 2.02oz.
800079	Shiny Black	57.4g / 2.02oz.

Zone



The Zone is a comfortable and lightweight yet sturdy "uni fit" (male/female) nylon frame. Don't let the sharp-looking streamlined design fool you, this frame is all function. Rubber nose pads and rubberized temple bars provide comfortable grip while the wrap around design ensures that your eyes are protected. All lenses are standard with .75 Pb Lead Equivalency, High-Lite Transmission Lenses with 20% More Visible Light, 1.80 High Index Lenses.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
870001	Silver	48.5g/1.7oz
870002	Red	48.5g/1.7oz
870003	Black	48.5g/1.7oz

Cyberflex Lead Glasses



Cyberflex lead glasses are a high quality safety frame offering excellent protection, large viewing area, adjustable temple bar and built-in side shields for added splash protection. Can also be ordered with radiation protective side shields.

- Adjustable temple bar
- Oversized protective lens area
- Splash protection
- Affordable price
- Optional .50mm Pb side shields (add -S)
- RX-able depending on suitability of prescription

Model	Color	Weight
683206	Clear/Blue	66.9g / 2.36oz.
683206-S	Clear/Blue - Side Shields	66.9g / 2.36oz.

Mirage



The all new Mirage frame fits well, offers a large protective area and is lightweight, yet extremely tough and durable. Along with the added protection of an ANSI Z87 safety rating, this gloss black frame comes with perforated rubber nosepieces for increased breathability. Mirage radiation protection glasses come standard with .75 Pb Lead Equivalency lenses and hard shell carrying case

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
875000		Black	48.5g/1.7oz

Fitovers, Lotus



- Wrap around Radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Velcro carrying case with sport clip
- Detachable neck cord
- Micro-fiber cleaning cloth
- **Designed to be worn over medium, oval or round prescription eyewear frame shapes**
- Non Rx

Model	Brand	Color	Weight
800304	Fitovers	Metallic Black	67.5g / 2.38oz
800305	Fitovers	Brushed Horn	67.5g / 2.38oz
800306	Fitovers	Smoke Marble	67.5g / 2.38oz
800307	Fitovers	Claret Stripe	67.5g / 2.38oz

Fitovers, Razor



- Wrap around Radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Velcro carrying case with sport clip
- Detachable neck cord
- Micro-fiber cleaning cloth
- **Designed to be worn over Small, Square or rectangle prescription eyewear frame shapes**
- Non Rx

Model	Brand	Color	Weight
800230	Fitovers	Black	65.6g / 2.31oz.
800235	Fitovers	Olive/Charcoal	65.6g / 2.31oz.
800240	Fitovers	Gun Metal	65.6g / 2.31oz.

Fitovers, Rhombus



- Wrap around Radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Velcro carrying case with sport clip
- Detachable neck cord
- Micro-fiber cleaning cloth
- **Designed to be worn over medium prescription eyewear frame shapes**
- Non Rx

Model	Brand	Color	Wight
800245	Fitovers	Matte Black	68.7g / 2.42oz
800250	Fitovers	Leopard/Black	68.7g / 2.42oz
800255	Fitovers	Blue Marble	68.7g / 2.42oz

Fitovers, Navigator



- Wrap around Radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Velcro carrying case with sport clip
- Detachable neck cord
- Micro-fiber cleaning cloth
- **Designed to be worn over medium-large prescription eyewear frame shapes**
- Non Rx

Model	Brand	Color	Weight
800200	Fitovers	Matte Black	69.8g / 2.46oz
800205	Fitovers	Tortoiseshell	69.8g / 2.46oz
800210	Fitovers	Olive Crystal	69.8g / 2.46oz

Fitovers, EcoFit



- Wrap around Radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Velcro carrying case with sport clip
- Detachable neck cord
- Micro-fiber cleaning cloth
- **Designed to be worn over prescription eyewear**
- Non Rx

Model	Brand	Color	Weight
832000	Fitovers	Black	92g
832001	Fitovers	Tortoise	92g

Flex-Fit Lead Glasses



- Lighter weight lead lens and frame
- Splash protection
- Optional leaded side shields (add-S)
- .50mm lead equivalent
- Rx-able depending on suitability of prescription

Model	Color	Weight
683208	Smokey Grey	68.5g / 2.42oz.
683208-S	Smokey Grey	68.5g / 2.42oz.

Incredibles



- Lightweight Acrylic Frames
- Full .50mm of front lens lead protection
- Comfortable to wear
- Padded Ear Temple Bars
- Built-In Retaining Strap
- Available in six colors
- Lens sizes: 61mm
- Rx-able depending on suitability of prescription

Model	Color	Weight
683260	Silver	48.5g / 1.7oz.
800064	Soft Matte Black	48.5g / 1.7oz.
800065	Sand Metallic Brown	48.5g / 1.7oz.
800066	Red	48.5g / 1.7oz.
800067	Black Graphic	48.5g / 1.7oz.
800068	Grey Carbon	48.5g / 1.7oz.

Maxx 30



The Maxx 30 frame is the first protective eyewear to receive the American Ophthalmology Associations seal of acceptance. The Maxx 30 Meets ANSI Z87.1+ standards along with ASTM F803 impact resistance standards. It has a unisex shape and temple bars and nose bridges that have soft inside rubberized padding for extra comfort. Every pair comes complete with their own semi-rigid case and adjustable nylon strap. All lenses are standard with .75 Lead Equivalency

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
830250	Black	48.5g/1.7oz
830251	Red	48.5g/1.7oz
830252	Silver	48.5g/1.7oz
830253	Brown	48.5g/1.7oz

Narcissus



This frame has a flexible padded temple design and a discrete inner cushion lining to ensure ultimate comfort.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
800084	Matte Black	56.7g / 2oz.
800085	Matte Grey	56.7g / 2oz.
800086	Shiny Gunmetal	56.7g / 2oz.

Premium Wraps



The only radiation protection frames on the market with lead glass lenses on the side shields for superior protection and excellent front and peripheral vision. Spring hinges, formable side temple bars and a fully adjustable nose piece make these glasses very versatile and comfortable. Large front lens area is ideal for regular and Bi-Focal prescriptions.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
683213	Tortoise Shell	71.3g / 2.52oz.

Xerxes



This frame offers a wrap-around shape that delivers more protection from radiation with nose pads for added comfort.

- High quality carrying case
- Always in stock with non prescription lenses
- Three great colors
- Lightweight and comfortable
- Large lens area for maximum protection
- NON PRESCRIPTION ONLY!!!!

Model	Color	Weight
800081	Shiny Black	60g / 2,11oz.
800082	Metallic Dove	60g / 2,11oz.
800083	Light Bronze	60g / 2,11oz.

Zephyr



These frames are sporty and comfortable with the padded nose bridge and sleek design. The Zephyr can hold higher prescriptions or can be worn as planos!

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Color	Weight
800069	Red	53.6g / 1.89oz.
800070	Blue	53.6g / 1.89oz.
800071	Plated Silver	53.6g / 1.89oz.
800073	Shiny Black	53.6g / 1.89oz.

Harley Davidson, 1000



These sporty safety frames from Harley Davidson offer fantastic wrap-around protection with large lenses and soft, adjustable temple bars for ultimate comfort.

- Available in three colors: Black, Gunmetal and Silver
- Lens sizes: 60mm
- Non Rx

Model	Brand	Color	Weight
800093	Harley Davidson	Black	69.6g / 2.45oz.
800094	Harley Davidson	Gunmetal	69.6g / 2.45oz.
800095	Harley Davidson	Silver	69.6g / 2.45oz.

Nicole Miller, Éclair



Beautifully shaped, the extra curves will set off high cheek bones and the contrasting colors add for extra fun, this frame is purely sinful!

- Available in two colors: Black Palm and Ambre
- Lens sizes: 51mm
- Rx-able depending on suitability of prescription

Model	Brand	Color	Weight
800308	Nicole Miller	Black Palm	48.9g / 1.7oz.
800309	Nicole Miller	Ambre	48.9g / 1.7oz.

Nicole Miller, Glamazon



Be inspired by the dramatic patterns of this elegant Nicole Miller frame. Global-fit bridge, spring hinges and triple zyl lamination.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800021	Nicole Miller	Leopard	49.5g / 1.92oz.
800022	Nicole Miller	Zebra	49.5g / 1.92oz.

Cocoons, Mini Line



Cocoons feature patented optical technology unique. All frames are engineered to comfortably float over prescription eyewear, but function equally as well when worn as conventional radiation protection glasses.

- High quality carrying case
- Always in stock
- Three great colors
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800900	Cocoons	Black	69.8g/2.46oz
800901	Cocoons	Tortoise	69.8g/2.46oz
800902	Cocoons	Sand	69.8g/2.46oz

Cocoons, Stream Line



Cocoons feature patented optical technology. All frames are engineered to comfortably float over prescription eyewear, but function equally as well when worn as conventional radiation protection glasses.

- High quality carrying case
- Always in stock
- Three great colors
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Wight
800905	Cocoons	Black	69.8g/2.46oz
800906	Cocoons	Tortoise	69.8g/2.46oz
800907	Cocoons	Sand	69.8g/2.46oz

Cocoons, Wide Line



Cocoons feature patented optical technology unique to sunglasses. All frames are engineered to comfortably float over prescription eyewear, but function equally as well when worn as conventional radiation protection glasses.

- High quality carrying case
- Always in stock
- Three great colors
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800910	Cocoons	Black	69.8g/2.46oz
800911	Cocoons	Tortoise	69.8g/2.46oz N/A
800912	Cocoons	Sand	69.8g/2.46oz

Cocoons, Slim Line



Cocoons feature patented optical technology unique to sunglasses. All frames are engineered to comfortably float over prescription eyewear, but function equally as well when worn as conventional radiation protection glasses.

- High quality carrying case
- Always in stock
- Three great colors
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800920	Cocoons	Black	69.8g/2.46oz
800921	Cocoons	Sand	69.8g/2.46oz
800922	Cocoons	Tortoise	69.8g/2.46oz

Bebe, Alter Ego



A fashion frame with a great work ethic, the Alter Ego from Bebe is a very high quality frame that works very well with radiation protection lenses. If you are looking for something comfortable yet functional, but aren't a fan of the industrial type of frames so prevalent in radiation protection eyewear, these are the ones for you.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800023	Bebe	Black List	42.3g / 1.5oz.
800024	Bebe	Brown	42.3g / 1.5oz.

Bebe, Goddess



Another high fashion frame with a great work ethic, the Goddess from Bebe is a very high quality frame that works very well with radiation protection lenses. Ideally suited to smaller, narrower faces. Like their counterpart the Alter Ego, if you are looking for something comfortable yet functional, but aren't a fan of the industrial type of frames so prevalent in radiation protection eyewear, these are the ones for you.

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800033	Bebe	Jet	48.7g / 1.71oz.
800034	Bebe	Scarlet	48.7g / 1.71oz.

Dakota Smith, Wild West



This Dakota Smith frame has an athletic look and feel with its wrap style and wide temples

- Prescription compatible (for most RX)
- High quality carrying case
- Always in stock with non prescription lenses
- Most prescriptions filled in two weeks or less
- Lightweight and comfortable
- Large lens area for maximum protection

Model	Brand	Color	Weight
800263	Dakota Smith	Amber	66.55g / 2.35oz.
800264	Dakota Smith	Black	66.55g / 2.35oz.

Standard Glasses Retainer



Soft cloth glasses retainer slips easily, yet securely on almost all protective eyewear. Slip not locking ball adjuster allows just the right amount of tension and holds firm until you release it. Easy on/easy off convenience, washable and long lasting. Comfortable, will not pull hair rub or chafe skin and is non allergenic.

Model	Color
683525	Blue

NEW Flip Top Radiation Protection Glasses



“Flip Top” radiation goggles are an all new product, designed to reduce scatter radiation by as much as 70% (actual reduction is dependent upon imaging technique used)

0.7mm Pb Equivalent acrylic shield is lightweight (45grams) offers a large protective area (lens area 2.5” (6cm) tall x 6.25 (16cm) long) and flips open when not needed. These new goggles also have an adjustable nose bridge and can be used as a stand alone radiation reduction goggle or fits easily over your existing glasses.

Available from stock in stylish, heavy duty amber acrylic with black and gray accents. Non RX



Model	Color	Weight
800275	Black/Grey	45g

NEW Laser Safety Glasses

All Laser eyewear is available in your choice of Glasses or Goggles. Please specify with order.



Model	Type of laser	Range OD Values in ()	Lens Color	VLT%*
900000	Homium.YAG,Co2	900-1000nm(5),1000-1500nm(7),1500-2400nm(7),2900-10600nm(7)	CLEAR	90
900001	Alexandrite Diode	730-855nm(5),755-830nm(7)	PINK	32,6
900002	Alx.Did. Hi Light Xmsn	632nm(2),650-690(3)nm,690-1000nm(6)	BLUE	55,5
900003	Argon Alignment	480-532nm(2,3),470-555nm(1.5)	RED	34
900004	Argon/KTP	190-532nm(5),5000-11000nm(6)	ORANGE	46,4
900005	Argon/KTP	190-532nm(5),5000-11000nm(6)	ORANGE	46,4
900006	Co2/Eximer	190-370nm(5),10600nm(6)	CLEAR	82
900007	Co2/Erbium	190-340nm(7),2760-3470nm(7),5000-11000nm(7)	CLEAR	85
900008	DIODE	600-900nm(6)	GREEN	27,9
900009	DIODE 810nm	800-830nm (5)	YEL/GRN	68
900010	Dye (CE Approved)	570-595nm(5),583-592nm(6)	PURPLE	34,4
900011	Intense Pulse Light	550-1100nm(1.5),620-1050nm(3)	GREEN	9
900012	HEliumNEon Alignmt	190-380nm(5),633nm(2)	BLUE	36,6
900013	YAG	190-380nm(5),1064nm(7),900-1100nm(5)	GREEN	58,7
900014	RUBY	680-710nm(4),690-700nm(7)	TEAL	33,8
900015	YAG Double Harmncs	190-534nm(5),910-1070nm(6)	BROWN	26,2
900016	Alx.Did. Hi Light Xms	700-1200nm(6)	BLUE	52,5
900017	M.W.YAG/Alx Did	190-532nm(7),825-890nm(4),890-965nm(5),965-1064nm(7),10600NM(5+)	AMBER	10,6
900018	Dye,Diode,HeNe,Ruby	513-680nm (6), 576-668nm(7)	BLUE 90	14

Prescription Order Form

Company Name: _____

Telephone number: _____

Company Purchase Order #: _____

Patient's Name: _____ Date: _____

Ship To Address: _____

City: _____ State: _____ Zip: _____

Glasses Model #: _____

Lead Equivalency:

.50mm .75mm

Distant:

	Sph.	Cyl.	Axis	Dec.	Prism	Base
R						
L						

Add for read:

	Segment Height	Segment Width	Segment Dec.	Total Dec.
R				
L				
	PD -	Distant	Near	

IMPORTANT: For proper prescription centering PD MUST BE INCLUDED.
 In order to process prescription orders, we require the dispensing opticians signature!

Dispensing Opticians Signature _____

Please photo copy this form and mail or fax to:

Infab Corporation
 3651 Via Pescador
 Camarillo, CA 93012
 Tel: (805) 987-5255 Fax: (805) 482-8424

