

Crutches

Invacare’s line of crutches provides support, strength and comfort for consumers. Crutches range in size from junior to tall adult, allowing a custom fit for a wide variety of consumers.



Invacare® Bariatric Crutches

Model nos. 8130-A, 8130-T

Strong and comfortable, Invacare’s bariatric crutches range in size from junior to tall adult. Large, double push buttons allow for height adjustments in one-inch increments. Crutch frame offers durable steel construction.

Features

- 1000 lb. weight capacity
- Push-button height adjustment
- Patient height guide on crutch for easy adjustment
- Contoured underarm design
- Arm pads (model no. 7131), handgrips (model no. 7126) and two-inch crutch tips (model no. 6135) also available as accessories

Bariatric Crutch Specifications

Model no.	Size	*Approximate client height	Adjustment range	Weight per crutch	Weight capacity	Units/carton	Limited warranty
8130-A	adult	5'2" - 5'10"	44" - 52"	4 lb. 6 oz.	1000 lb.	1 pair	3 years
8130-T	tall adult	5'10" - 6'6"	52" - 60"	5 lb.	1000 lb.	1 pair	3 years

*Client height is approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

Invacare® Bariatric Crutch Accessories For Model No. 8130 Crutches



Model no. 6135



Model no. 7131



Model no. 7126

Rubber Tips

- Extra large, 2" outside diameter (O.D.) base

Arm Pads

- Contoured underarm design

Handgrips

- Solid handgrips



Model no. 8115-T Model no. 8115-A Model no. 8115-J



Push handle lock for quick handgrip adjustment.



Model no. 8120-T Model no. 8120-A Model no. 8120-J



Model no. 6153

Invacare® Quick-Change Crutches

Model no. 8115-J – Junior
 Model no. 8115-A – Adult
 Model no. 8115-T – Tall Adult

Invacare quick-change crutches offer a solid one-piece crutch frame in lightweight anodized aluminum with push-button height adjustment and contoured handgrips.

Invacare® Quick-Adjust Crutches

Model no. 8120-J – Junior
 Model no. 8120-A – Adult
 Model no. 8120-T – Tall Adult

The one-piece handgrip assembly eliminates the need for wing nuts and other assembly pieces. Plus, the push-button height adjustment and patient height guide allow the crutch to be fitted to the proper user height. This patented model comes completely assembled for the consumer's convenience.

Invacare® Forearm Crutches

Model no. 6153 – Adult

Forearm crutches feature vinyl-coated one-inch solid aluminum cuffs, cushioned vinyl handgrips and extra large crutch tips.

Quick-Change/Quick-Adjust Crutch Specifications

Model no.	Size	*Approximate client height	Adjustment range	Weight per pair	Weight capacity	Units/carton	Limited warranty
8115-J	junior	4' 6" - 5' 2"	36" - 44"	3 lb. 4 oz.	350 lb.	8 pairs	3 years
8115-A	adult	5' 2" - 5' 10"	44" - 52"	3 lb. 8 oz.	350 lb.	8 pairs	3 years
8115-T	tall adult	5' 10" - 6' 6"	52" - 60"	3 lb. 10 oz.	350 lb.	8 pairs	3 years
8120-J	junior	4' 6" - 5' 2"	36" - 44"	3 lb. 7 oz.	250 lb.	8 pairs	3 years
8120-A	adult	5' 2" - 5' 10"	44" - 52"	3 lb. 11 oz.	250 lb.	8 pairs	3 years
8120-T	tall adult	5' 10" - 6' 6"	52" - 60"	4 lb.	250 lb.	8 pairs	3 years

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

Forearm Crutch Specifications

Model no.	Size	*Approximate client height	Height adj. to center of handle	Top of cuff to center of handle	Weight per pair	Weight capacity	Units/carton	Limited warranty
6153	adult	5' 0" - 6' 2"	28" - 37"	9 1/2" - 11"	4 lb. 8 oz.	300 lb.	6 pairs	3 years

*All client height ranges are approximations. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

Invacare® Crutch Accessories for Model No. 8115 and 8120 Crutches



Model nos. 6134, 6135

Rubber Tips

Model no. 6134
 Rubber Tips – Large 1 3/4" O.D. Base
 Model no. 6135
 Rubber Tips – Extra Large 2" O.D. Base

- Help grip the surface and absorb shock
- Reinforced with metal inserts for longer-lasting wear



Model no. 6131

Arm Pads

Model no. 6131 – Arm Pads

- Add comfort to armrests



Model nos. 6125*, 6126** (not shown)

Handgrips

Model nos. 6125, 6126 – Solid Handgrips

- Padded to keep hands and wrists comfortable
- Contoured grips that limit rotation

Model no.	Only available for Model nos.
* 6125	8120-J, 8120-A and 8120-T crutches
** 6126	8115-J, 8115-A and 8115-T crutches