

KIMBERLY-CLARK\* Headwear



Kimberly-Clark offers a variety of protective headwear, with styles and protection levels for use throughout the hospital – from the cafeteria to the operating room. Bouffant caps, surgical caps and surgical hoods keep workers comfortable and well-protected, and help keep you in compliance with OSHA standards.

**Bouffant Caps.** Lightweight and easy to wear, Kimberly-Clark Bouffant Caps are available in both a 1-layer breathable Spunbond fabric and a 3-layer SMS fabric. The caps are designed to cover the hair in no fluid contamination or in light fluid contamination environments.

**Surgical Caps.** Two surgical cap styles are available to protect workers who face varying degrees of risk, depending on the procedure.

- KAYCEL\* fabric caps are best suited for environments in which no fluid contamination is expected.
- Protective 3-layer fabric caps are the choice when light fluid contamination from splashing or spraying is anticipated.

**Surgical Hoods.** Made from a protective 3-layer fabric, Kimberly-Clark Surgical Hoods come in two models – a tie-at-the-neck style for light fluid contact, and a full-coverage hood for moderate fluid contact.

From protective apparel to surgical products – everything we do, everything we make, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

**VALUEselect® Bouffant Cap**

Coverage: Hair  
 Materials: 1-layer Spunbond fabric  
 Exposure Anticipated\*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69064	VALUEselect® Bouffant Cap, Blue	Regular	500	1 (Bulk)
69067	VALUEselect® Bouffant Cap, Blue	Regular	300	3
69081	VALUEselect® Bouffant Cap, White	Regular	500	1 (Bulk)
69084	VALUEselect® Bouffant Cap, White	Regular	300	3
69071	VALUEselect® Bouffant Cap, Blue	Large	500	1 (Bulk)
69074	VALUEselect® Bouffant Cap, Blue	Large	300	3
69090	VALUEselect® Bouffant Cap, White	Large	500	1 (Bulk)
69094	VALUEselect® Bouffant Cap, White	Large	300	3

**Bouffant Cap**

Coverage: Hair  
 Materials: 1-layer Spunbond fabric  
 Exposure Anticipated\*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69054	Bouffant Cap, Blue	Regular	500	1 (Bulk)
69221	Bouffant Cap, Blue	Regular	300	3
70535	Bouffant Cap, Cushion Band, Blue	Regular	125	4
69220	Bouffant Cap, White	Regular	300	3
69222	Bouffant Cap, Green	Regular	300	3
69053	Bouffant Cap, Blue	Large	500	1 (Bulk)
69231	Bouffant Cap, Blue	Large	300	3
69230	Bouffant Cap, White	Large	300	3
69232	Bouffant Cap, Green	Large	300	3
69004	Bouffant Cap, Blue	X-Large	300	3
60212	Bouffant Cap, Cushion Band, Blue	X-Large	250	3

**Protective Bouffant Cap**

Coverage: Hair  
 Materials: Protective 3-layer SMS fabric  
 Exposure Anticipated\*: Potential for light fluid contamination from splashing or spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69051	Protective Bouffant Cap, Floral	Large	240	3
69055	Protective Bouffant Cap, White	X-Large	210	3
69057	Protective Bouffant Cap, Contemporary	Large	240	3
69058	Protective Bouffant Cap, Contemporary	Regular	240	3
60400	Protective Bouffant Cap, Cushion Band, Rayon, Blue	Adjustable	75	4

**Non-Elastic Tie Cap**

Coverage: Hair  
 Materials: Protective 3-layer SMS fabric  
 Exposure Anticipated\*: Potential for light fluid contamination from splashing or spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69060	Non-Elastic Tie Cap, Blue	Universal	300	3

**Surgical Cap**

Coverage: Hair  
 Materials: KAYCEL® Fabric  
 Exposure Anticipated\*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69240	Surgical Cap, Blue	Universal	300	3

**Protective Surgical Cap**

Coverage: Hair  
 Materials: Protective 3-layer SMS fabric on sides, 1-layer fabric on top  
 Exposure Anticipated\*: Potential for light fluid contamination from splashing and spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69520	Protective Surgical Cap, Blue	Universal	300	3

**Protective Surgical Hood – Tie Neck**

Coverage: Hair  
 Materials: Protective 3-layer SMS fabric, 1-layer fabric on top  
 Exposure Anticipated\*: Potential for light fluid contamination from splashing and spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69110	Protective Surgical Hood, Blue	Universal	300	3

- OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate protective clothing that prevents contamination of an employee's skin or clothing by blood or other potentially infectious materials.



**Commitment to Excellence**

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at [www.kchealthcare.com](http://www.kchealthcare.com).

\* Registered Trademark or Trademark of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates. © 2001 KCC. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Kimberly-Clark's Neenah Paper.

KL-3638-B H5121

## KIMBERLY-CLARK\* Footwear



The Kimberly-Clark footwear product line features innovative fabrics and styling to deliver the right amount of protection for the task, from basic shoe covers to full-coverage boots.

**Basic Shoe Covers.** BASIC\*, Heavy-Duty or X-TRA TRACTION\* shoe covers come in universal and extra large sizes with or without traction.

**ANKLE GUARD\* Shoe Covers.** Designed to cover the shoe and go over the ankle, ANKLE GUARD\* does the job when moderate fluid contact and/or splashing and spraying above the shoe are anticipated.

- Made of Kimberly-Clark's own protective 3-layer fabric.
- Feature elastic at the ankle for better fit and protection.

**HI GUARD\* Full Coverage Boot.** Step up coverage with KIMBERLY-CLARK\* Knee High Shoe Covers, which cover the shoe and lower leg.

- Depending on the risk, choose from Regular or Ultra HI GUARD\* Full-Coverage Boots.
- Made from a protective 3-layer fabric, and come in a range of sizes.
- Feature a plastic film coating, traction slips, elastic at the ankle and ultra-sonically bonded seams.

From protective apparel to surgical products – everything we do, everything we make, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

**VALUEselect® Shoe Cover**

Coverage: Shoe only  
 Materials: 1-layer Spunbond fabric  
 Exposure Anticipated\*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69097	VALUEselect® Shoe Cover	Regular	500	1 (Bulk)
69101	VALUEselect® Shoe Cover	Regular	100	3
69104	VALUEselect® Shoe Cover	X-Large	500	1 (Bulk)
69108	VALUEselect® Shoe Cover	X-Large	100	3
69111	VALUEselect® Shoe Cover with Traction	Regular	500	1 (Bulk)
69114	VALUEselect® Shoe Cover with Traction	Regular	100	3
69117	VALUEselect® Shoe Cover with Traction	X-Large	500	1 (Bulk)
69121	VALUEselect® Shoe Cover with Traction	X-Large	100	3

**BASIC® Shoe Cover**

Coverage: Shoe only  
 Materials: Protective, 3-layer SMS fabric  
 Exposure Anticipated\*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69561	BASIC® Shoe Cover	Universal	300	3
69562	Value Pack	Universal	400	1 (Bulk)
69005	Value Pack	X-Large	380	1 (Bulk)

**X-TRA TRACTION® Shoe Cover**

Coverage: Shoe only  
 Materials: Protective 3-layer SMS fabric, adhesive traction strips  
 Exposure Anticipated\*: Light fluid contamination

Stock #	Description	Size	Items/Case	Boxes/Case
69252	X-TRA TRACTION® Shoe Cover	Universal	300	3
69254	X-TRA TRACTION® Shoe Cover	X-Large	240	3
69352	BASIC X-TRA TRACTION®	Universal	300	2
69015	BASIC X-TRA TRACTION®	X-Large	300	2

**KIMBERLY-CLARK® Heavy-Duty Shoe Cover**

Coverage: Shoe only  
 Materials: Heavy-duty protective 3-layer SMS fabric  
 Exposure Anticipated\*: Light fluid contamination

Stock #	Description	Size	Items/Case	Boxes/Case
69551	Heavy-Duty Shoe Cover	Universal	300	3
69324	Heavy-Duty, X-TRA TRACTION®	Universal	300	2
69345	Heavy-Duty, X-TRA TRACTION®	X-Large	240	2

**ANKLE GUARD® Shoe Cover**

Coverage: Over the ankle  
 Materials: Protective 3-layer SMS fabric, adhesive traction strips  
 Exposure Anticipated\*: Moderate fluid contamination with splashing or spraying above the shoe

Stock #	Description	Size	Items/Case	Boxes/Case
69253	ANKLE GUARD® Shoe Cover	Universal	270	3
69353	ANKLE GUARD® Shoe Cover	X-Large	270	3

**HI GUARD® Regular Full-Coverage Boot**

Coverage: Shoe and lower leg  
 Materials: Protective, 3-layer SMS fabric, coated fabric on foot area, adhesive traction strips  
 Exposure Anticipated\*: Potential for moderate fluid contamination with splashing

Stock #	Description	Size	Items/Case	Boxes/Case
69571	HI GUARD® Regular Full-Coverage Boot	Universal	150	3
69671	HI GUARD® Regular Full-Coverage Boot	X-Large	150	3

**HI GUARD® Ultra Full-Coverage Boot**

Coverage: Shoe and lower leg  
 Materials: Plastic film coating over 3-layer SMS fabric, two 1" urethane foam traction strips  
 Exposure Anticipated\*: Potential for heavy fluid contamination with splashing

Stock #	Description	Size	Items/Case	Boxes/Case
69572	HI GUARD® Ultra Full-Coverage Boot	Universal	120	4
69672	HI GUARD® Ultra Full-Coverage Boot	X-Large	120	4

\* OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate protective clothing that prevents contamination of an employee's skin or clothing by blood or other potentially infectious materials.

**Caution:** Footwear users should be aware that a risk of slipping exists. Hospital floors may be waxed or wet. Good technique and OSHA require that footwear be removed when contaminated or upon leaving the work area. Do not use nonconductive shoe covers in locations where flammable anesthetics are administered.

**Note:** Shoe dimensions can vary from street shoes to athletic shoes. When selecting a shoe cover size, please consider the outside dimensions of the shoe to ensure a proper fit.



**Commitment to Excellence**

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at [www.kchealthcare.com](http://www.kchealthcare.com).

\* Registered Trademark or Trademark of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates. © 2001 KCC. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Kimberly-Clark's Neenah Paper.

KL-3638-C H5121