

Aquaflor™ Hydrogel Wound Dressings



Aquaflor™ Hydrogel Wound Dressing

Order Code	Description	Ship Case
8884476139	Aquaflor™ 3" Disk Hydrogel Wound Dressing 5/Box	50
8884476154	Aquaflor™ 4.75" Disk Hydrogel Wound Dressing 5/Box	50

Features and Benefits

Features	Benefits
Advanced formulation is extremely absorbent and will not dry out	Longer wear time and less dressing changes
Unique disk shape maximizes wound coverage and helps fill shallow cavities	Protects wounds as it provides a moist environment
Dressing may be cut to size	Allows customization
Translucent gel allows wound visualization	No need to remove dressing to observe wound
Cooling sensation helps reduce pain upon dressing change	Will not cause trauma upon dressing removal

Frequently Asked Questions

Question	Answer
How should the dressing be applied?	Dressing should overlap the wound by 1cm (1/4 inch).
How often should the dressing be changed?	Change every 1-7 days, depending on drainage.
What kind of secondary dressings can be used?	Transparent film, island dressings, roll gauze, tapes, gauze, ABD pads

Curafil™* Hydrogel Impregnated Gauze



Curafil™* Hydrogel Impregnated Gauze

Order Code	Description	Ship Case
9255	Curafil™* Hydrogel Impregnated Gauze 2" X 2", Pad (.25 oz gel) 25/Box	100
9256	Curafil™* Hydrogel Impregnated Gauze 4" X 4", Pad (1.0 oz gel) 25/Box	100
9257	Curafil™* Hydrogel Impregnated Gauze 4" X 8", Pad (2.0 oz gel) 25/Box	100
9259	Curafil™* Hydrogel Impregnated Gauze 1" X 36" Packing Strip 12/Box	72

Features and Benefits

Features	Benefits
Clear Hydrogel impregnated in a gauze sponge	Provides optimal environment for moist wound healing
Viscous formula	Less mess and waste
Unique glycerin formulation	Will not dry out
Soft gel impregnated into soft sponge	Conforms to wound bed

Frequently Asked Questions

Question	Answer
What is the wear time?	One to seven days depending on level of drainage

Curafil™* Gel Wound Dressing



Curafil™* Amorphous Hydrogel Wound Dressing

Order Code	Description	Ship Case
9250	Curafil™* Gel Wound Dressing 1/2 Ounce Tube	24 Tubes
9251	Curafil™* Gel Wound Dressing 1 Ounce Tube	12 Tubes
9252	Curafil™* Gel Wound Dressing 3 Ounce Tube	12 Tubes

Features and Benefits

Features	Benefits
High moisture content	Promotes moist wound healing environment
Contains Germaben	Extends the life of the gel in the tube
Clear and viscous	Allows wound visualization, conforms to the wound

Frequently Asked Questions

Question	Answer
How is product applied?	Gel may be applied directly to the wound bed, or on a sterile dressing or Q-Tip™* swabs
What is the wear time?	Gel should be changed daily depending on wound drainage