

Under Cast Padding – Covidien/ Kendall



WEBRIL™* 100% Cotton Undercast Padding

WEBRIL™* Undercast Padding is constructed of 100% cotton fibers for soft cushioning to be used with either synthetic or plaster casting material.

Order Code	Description	Ship Case
1418	WEBRIL™* 100% Cotton Undercast Padding, 2" X 4 Yds, 24/Bag	72 Rolls
2059	WEBRIL™* 100% Cotton Undercast Padding, 3" X 4 Yds, 12/Bag	72 Rolls
3175	WEBRIL™* 100% Cotton Undercast Padding, 4" X 4 Yds, 12/Bag	72 Rolls
3489	WEBRIL™* 100% Cotton Undercast Padding, 6" X 4 Yds, 6/Bag	36 Rolls
2283	WEBRIL™* 100% Cotton Undercast Padding, 2" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2394	WEBRIL™* 100% Cotton Undercast Padding, 3" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2502	WEBRIL™* 100% Cotton Undercast Padding, 4" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2554	WEBRIL™* 100% Cotton Undercast Padding, 6" X 4 Yds, 1 roll/pouch, Sterile	24 Rolls

Features and Benefits

Features	Benefits
Constructed of 100% purified cotton fibers for soft and natural cushioning	Ideal for use with synthetic and plaster casting material
Processed to minimize bulk	Minimize dimensional changes that lead to movement of padding under cast
Retains consistency wet and dry	Product will not collapse as padding absorbs moisture
Shortest fibers, hook & loop effect	Prevents pressure points
Shortest fibers, hook & loop effect	Combination of mild stretch and cohesiveness, holds padding in place without shifting or bunching
Available in a sterile blister pouch for postoperative casting applications	Patient safety



WEBRIL™* II 100% Cotton Undercast Padding

WEBRIL II™* Undercast Padding is constructed of 100% cotton fibers with crimped finish for extra loft and conformability.

Order Code	Description	Ship Case
4095	WEBRIL™* II 100% Cotton Undercast Padding, 2" X 4 Yds, 24/Bag	72 Rolls
4152	WEBRIL™* II 100% Cotton Undercast Padding, 3" X 4 Yds, 12/Bag	72 Rolls
4221	WEBRIL™* II 100% Cotton Undercast Padding, 4" X 4 Yds, 12/Bag	72 Rolls
4519	WEBRIL™* II 100% Cotton Undercast Padding, 6" X 4 Yds, 6/Bag	36 Rolls
2666	WEBRIL™* II 100% Cotton Undercast Padding, 2" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2754	WEBRIL™* II 100% Cotton Undercast Padding, 3" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2847	WEBRIL™* II 100% Cotton Undercast Padding, 4" X 4 Yds, 1 roll/pouch, Sterile	50 Rolls
2944	WEBRIL™* II 100% Cotton Undercast Padding, 6" X 4 Yds, 1 roll/pouch, Sterile	24 Rolls

Features and Benefits

Features	Benefits
Constructed of 100% purified cotton fibers for soft and natural cushioning	Ideal for use with synthetic and plaster casting material
Processed to minimize bulk	Minimize dimensional changes that lead to movement of padding under cast
Retains consistency wet and dry	Product will not collapse as padding absorbs moisture
Shortest fibers, hook & loop effect	Prevents pressure points
Shortest fibers, hook & loop effect	Combination of mild stretch and cohesiveness, holds padding in place without shifting or bunching
Available in a sterile blister pouch for postoperative casting applications	Patient safety



TENDEROL™ Synthetic Undercast Padding

TENDEROL™ Synthetic Undercast Padding is constructed of 100% polyester fibers, available in multiple sizes.

Order Code	Description	Ship Case
6242	TENDEROL™ Synthetic Undercast Padding, 2" X 4 yds, 24 Per Bag, 3 Bags/Case	72 Rolls
6243	TENDEROL™ Synthetic Undercast Padding, 3" X 4 yds, 12 Per Bag, 6 Bags/Case	72 Rolls
6244	TENDEROL™ Synthetic Undercast Padding, 4" X 4 yds, 12 Per Bag, 6 Bags/Case	72 Rolls
6246	TENDEROL™ Synthetic Undercast Padding, 6" X 4 yds, 6 Per Bag, 6 Bags/Case	36 Rolls

Features and Benefits

Features	Benefits
Constructed of 100% polyester fibers	Absorbs excess moisture, then dries quickly, retaining its bulk for a consistent, comfortable cast fit