

# Rocker & Rotator

## ROCKER

The LW Blood Rocker gently rocks blood sample tubes and small vials. The gentle motion assures mixing without damage to fragile cell suspensions.

### Features

- 16-tube capacity
- Space-saving design
- Washable tube pad
- 24 rocking cycles per minute

### Specifications

- Speed: 24 cycles per minute
- Motor: 120V AC 3W
- Height: 4" (10.2cm)
- Length: 13.5" (34.3cm)
- Width: 5.25" (13.3cm)
- Weight: 4.8lbs (2.2kg)



## Lab Equipment

Model #	Description
MXL-BLF7-16T1	Blood Rocker
RTL-BLV7-2501	Rotator

## ROTATOR

The Serology Rotator is designed to give years of dependable service, even in the busiest laboratory. It's the ideal space-saving instrument for procedures that require constant or timed rotation (mixing) such as RPR, VDRL, latex tests and for procedures that require a nonstandard speed.

### Features

- Variable speed (45-250 rpm)
- Ideal for all routine testing
- 30-minute timer
- Built-in analog tachometer
- Timed or continuous operations
- Space-saving design
- Suction-cupped feet
- Hold-down spring for small flasks
- Includes moisture cover for RPR
- Orbital diameter: 22mm / 0.87"

### Specifications

- Speed: 45-250 rpm
- Motor: 120V AC 10W
- Height: 6" (15.2cm)
- Length: 9" (22.9cm)
- Width: 9" (22.9cm)
- Weight: 9.7lbs (4.4kg)
- Platter: 9" x 9" (22.9cm x 22.9cm)



P 770.270.1394

F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com