

CxR

Horizontal centrifugation offers advantages over traditional angled-rotor centrifuges, such as tighter straight-line gel separations, fewer re-spins and re-mixes, and fewer high-potassium results. The new CxR centrifuge delivers 1800 G-Force for quick spin down of serum-separator blood tubes. Low-maintenance features: self contained rotor, whisper-quiet operation, brushless no-maintenance motor, automatic braking, simple operation...set the timer and go!

ADVANTAGES

- Low profile Saves Counter Space
- Tight gel line separation
 - easier reads, and reduced re-mixes
- 6 Tube Capacity (accommodates 3-15ml tubes)
- Large 17.5mm x 133mm tube capacity

SPECIFICATIONS

- Fixed speed - 3400rpms
- 1800 G-Force
- Brushless maintenance-free motor
- 6 Place (3-15 ml tube) Rotor
- Horizontal Spin Down
- Front Access Timer Control Knob

CONTROLS

- 15 minute mechanical timer
- On/off switch on back

WARRANTY

One year on materials and workmanship

DIMENSIONS AND WEIGHT

- Height: 8.75"
- Length: 15"
- Width: 13.5"
- Weight: 23.6 lbs
- Height w/ lid open: only 19.75"



Centrifuge

CxR

Model #	Description
CXC-06SF-1501	6 place horizontal swing out centrifuge

•Economically priced for reference labs!



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30045

LWScientific.com